Para-Pak® TRICHROME BLUE

REF 400501

IVD

Rx Only

INTENDED USE

Trichrome Blue is a stain which provides excellent differentiation of microsporidium spores in fresh and preserved stool, urine sediment, nasopharyngeal aspirates and other body fluids.

SUMMARY AND EXPLANATION OF THE TEST

Microsporidia are obligate intracellular spore-forming protozoa belonging to the phylum Microspora. Species from five genera: *Nosema*, *Encephalitozoon*, *Pleistophora*, *Septata*, and *Enterocytozoon* have been implicated as disease-causing parasites in humans, especially in immunocompromised individuals such as AIDS patients.¹⁻³

Definitive diagnosis has been impeded by the difficulty of detecting these protozoa. In 1992, Weber et.al.⁴ described a procedure using a modified Trichrome stain for the staining of Microsporidia from formalin-fixed stool and duodenal aspirates. While the Weber stain constituted a significant improvement, the general lack of background color makes the screening of large number of slides very difficult and time-consuming.

Ryan et.al. (1993) further modified the Trichrome stain, substituting the fast green in the Weber stain with aniline blue.⁵ The last modification stained the background material a deep blue color, which nicely contrasted with the pinkish-red stained microsporidia. Detailed internal structure can also be seen using this stain. Therefore, this method opens up a practical means for routine screening of microsporidia in fresh and preserved stools, urine sediments, induced sputum and nasopharyngeal aspirates.⁵ It has also been found that the increase of staining temperature can improve intensity and shorten time of staining.⁶

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this kit is listed on the outer box.

1. Chromotrope 2R
2. Aniline Blue
3. Phosphotungstic Acid
4. Glacial Acetic Acid
5. Deionized H₂O

MATERIALS NOT PROVIDED

1. Positive and Negative Controls
2. 10% Formalin, collection containers
3. Glass slides, coverslips
4. Staining containers, stain rack
5. Absolute Methanol
6. Acid Alcohol (0.45% acetic acid in 90% Ethanol)
7. 95% Ethanol
8. HEMO-DE (Xylene substitute)
9. Mounting medium
10. Microscope (100X oil immersion objective)

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. The reagents may cause irritation to skin, eyes, and respiratory tract. Avoid breathing vapor and eye/skin contact.
3. Every specimen represents a potential source of infectious material and should be handled accordingly.
4. Never return used stain to the original bottle. Discard appropriately after use.

HAZARD and PRECAUTIONARY STATEMENTS

 Trichrome Blue	Signal word Danger Hazard statements H314 - Causes severe skin burns and eye damage H332 - Harmful if inhaled H412 - Harmful to aquatic life with long lasting effects Precautionary Statements - EU (S28, 1272/2008) P280 - Wear protective gloves/ protective clothing/ eye protection/ face protection P260 - Do not breathe dust/fume/gas/mist/vapors/spray P303 + P361 + P353 - IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/shower P280 - Wear eye protection/ face protection P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing P310 - Immediately call a POISON CENTER or doctor/ physician
---	---

SHELF LIFE AND STORAGE

The expiration date of the Trichrome-blue stain is indicated on the kit label. Store at 15-30 C. Excessive heat and cold should be avoided.

TEST PROCEDURE

A. Specimen Preparation

Stool specimens preserved in 10% formalin, 5% formalin, or sodium acetate-acetic acid formalin (SAF) may be used. The preservative-stool mixture must be allowed to stand for a minimum of 30 minutes at room temperature to ensure adequate fixation. Urine specimens can be sedimented by centrifugation at 15,000 xg for 30 seconds or 400 xg for 10 minutes.³ Then resuspend in 10% buffered formalin at one-tenth (1/10) of the original volume. For best results, 2-3 mL of preserved stool specimens may be concentrated prior to processing using the formalin-ethyl acetate sedimentation procedure (centrifuge concentrate for 10 minutes at 500 xg).⁷ Specimens should be handled according to established guidelines.^{8,9}

B. Slide Preparation

1. With a calibrated loop or pipet, place 10 µL of the fixed stool concentrate onto a glass slide.
2. Spread thinly over an area approximately 45 by 25 mm.
3. Allow slides to dry completely (usually room temperature overnight or 37 C for two hours).
4. Fix smears in methanol for five minutes and allow to air dry.

C. Staining

- 1a. Heat Trichrome-blue stain* by placing in a 37 C incubator for at least one hour. Then add slides to stain for 30 minutes at 37 C. (Recommended procedure)
- OR
- 1b. Trichrome-blue stain can be used at room temperature. Place slides in stain for 90 minutes.

Note: All subsequent dehydration steps (2-6) should be carried out at room temperature.

2. Rinse with acid alcohol for 10 seconds.
3. Rinse with 95% ethanol for 10 seconds.
4. Dehydrate twice in 95% ethanol or five minutes each.
5. Dehydrate in HEMO-DE cleaning agent for five minutes.
6. Mount with Permout, then coverslip.**

*Hand-shake the Trichrome-blue stain well before dispensing. If precipitates are noticeable, filter before use.

**If permanent record is desired, avoid exposure to light.

INTERPRETATION OF RESULTS

Microsporidia: Spore size approximately 1.5 X 0.9 µm (as large as 3.5 X 2.0 µm for Septata). Spores display a characteristic refractile ovoid shape with a clear vacuole-like zone. The spore wall stains a bright pinkish-red, while the spore contents remain either transparent or demonstrate a pinkish-red belt-like stripe.

Yeasts: Larger and more variable in size. Stain reddish or light blue and more uniformly; sometimes show budding.

Other parasites may stain gray to dark blue. Background material and artifacts will stain blue.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

A. **Positive and Negative Control**

Microsporidium positive and negative control samples or QC slides should accompany each run. Criteria of acceptance is the same as in "Interpretation of Results".

B. **Quality of the reagents should be evaluated each time a patient's specimen is tested.**

1. All staining dishes should be covered to prevent evaporation of reagents.
2. Depending on the volume of slides stained, reagents will need to be changed on an as-needed basis (for example, after 30 slides).
3. The purpose of the first alcohol rinse after decolorizing is to halt decolorization. Since carryover of the acid alcohol will occur, change this solution frequently.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

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REFERENCES

1. Shadduck JA and Greeley E. Microsporidia and human infections. Clin Micro Rev. 2(2):158-165. 1989
2. Orenstein JM. Microsporidiosis in the acquired immunodeficiency syndrome. J Parasitol. 1991;77(6):843-864.
3. Aldras AM, Orenstein JM, Kotler DP, Shadduck JA and Didier ES. Detection of microsporidia by indirect immunofluorescence antibody test using polyclonal and monoclonal antibodies. J Clin Micro 1994;32(3):608-612.
4. Weber R, Bryan RT, Owen RL, Wilcox CM, Gorelkin L, Visvesvara GS and The enteric opportunistic infections working group improved light-microscopical detection of microsporidia spores in stool and duodenal aspirates. N Engl J Med 1992;32(63):161-166.
5. Ryan NJ, Sutherland G, Coughlan K, Globan M, Doultree J, Marshall J, Baird RW, Pedersen J, and Dwyer B. A new trichrome-blue stain for detection of microsporidial species in urine, stool and nasopharyngeal specimens. J Clin Micro 1993;31(12):3264-3269.
6. Kokoskin E, Gyorkos TW, Campus A, Cedilotte L, Purtil T and Ward B. Modified technique for efficient detection of microsporidia. Clin Micro 1994;32(4):1074-1075.
7. Yang J and Scholten T. A fixative for intestinal parasites permitting the use of concentration and permanent staining procedures. Am J Clin Pathol. 1977;67:300-304.
8. Hausler Jr WJ, Herrman KL, Isenberg HD and Shadomy HJ. (Eds). Manual of clinical microbiology, 5th Ed. 1991. ASM.
9. Garcia LS and Bruckner DA. Diagnostic medical parasitology, 2nd ed. Elsevier, New York. 1993.

SN10902

REV. 10/22

 Manufactured By	Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com
--	--

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/package of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.