**DNA Amplification Assay for the Detection of
Herpes Simplex Virus Type 1 and Herpes Simplex Virus Type 2**

REF 480650

IVD

R_x Only

INTENDED USE

The Alethia HSV 1&2 DNA amplification assay, performed on the Alethia reader is a qualitative in vitro diagnostic test for the direct detection and differentiation of herpes simplex virus type 1 (HSV-1) and herpes simplex virus type 2 (HSV-2) DNA in cutaneous and mucocutaneous lesion specimens from male and female patients suspected of Herpetic infections.

Alethia HSV 1&2 utilizes loop-mediated isothermal DNA amplification (LAMP) technology to detect HSV-1 and HSV-2 by targeting segments of the herpes simplex virus 1 and herpes simplex virus 2 genomes. Results from Alethia HSV 1&2 are used as an aid in the diagnosis of HSV infection in symptomatic patients.

The assay is intended for use in hospital, reference or state laboratory settings. This device is not intended for nonlaboratory point-of-care use.

WARNING: Alethia HSV 1&2 is not FDA cleared for use with cerebrospinal fluid (CSF) or to aid in the diagnosis of HSV infections of the central nervous system (CNS). The device is not intended for prenatal screening.

SUMMARY AND EXPLANATION OF THE TEST

Alethia HSV 1&2 is based on loop-mediated amplification (LAMP)^{1,2} technology. The assay individually targets regions of the HSV-1 and HSV-2 genomes in two separate Test Devices. The HSV-1 target is a 208 base pair (bp) sequence of the HSV-1 glycoprotein G (US4) gene. The HSV-2 target is a 189 bp sequence of the HSV-2 glycoprotein G (US4) gene.

Loop-mediated amplification uses specially designed primers to provide for specific and continuous isothermal DNA amplification. A by-product of amplification is magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. Reaction solution absorbance characteristics are monitored by the Meridian Alethia Incubator/Reader. Changes in reaction solution absorbance characteristics created by precipitation of magnesium pyrophosphate indicate the presence of target DNA. The absence of target DNA results in no significant change in sample absorbance.

The Alethia HSV 1&2 kit includes Alethia Sample Preparation Apparatus III (SMP PREP III), Alethia HSV 1 Test Devices, Alethia HSV 2 Test Devices, Mineral Oil, Heat Treatment Tubes, and 50 µL transfer pipettes. The Alethia SMP PREP III, used for specimen dilution and preparation, includes a Tris-buffered solution containing formalin-treated *E. coli* harboring *Staphylococcus aureus* DNA. The Alethia HSV 1 Test Device contains one lyophilized amplification reagent bead in each of two chambers: a TEST chamber with HSV-1-specific primers and a CONTROL chamber with *S. aureus*-specific primers. The Alethia HSV 2 Test Device contains one lyophilized amplification reagent bead in each of two chambers: a TEST chamber with HSV-2-specific primers and a CONTROL chamber with *S. aureus*-specific primers. The *S. aureus* DNA in SMP PREP III and the *S. aureus*-specific primers in the Test Device CONTROL chambers function as the Internal Control for the assay. During specimen preparation, each patient specimen is added to the SMP PREP III containing Assay Control and combined with the *S. aureus* DNA prior to amplification. Addition of *S. aureus* DNA to the patient sample allows for parallel processing of target DNA and control DNA through amplification and detection. The Internal Control monitors amplification inhibition, assay reagent performance and sample processing effectiveness. The control *S. aureus* target must be amplified and detected in the final reaction or the test is considered invalid and results are not reported.

The Alethia reader monitors changes in absorbance characteristics by measuring transmission of light through the TEST and CONTROL reaction solutions. Light transmission is determined at the assay Run Start (Signal_{initial}, S_i) and at the assay Run End (Signal_{final}, S_r). The Alethia reader calculates the change in light transmission between Run End and Run Start (S_r:S_i) and compares the ratio to a fixed cut-off value.

Fixed cut-off values for the TEST chamber are used to report sample results. TEST chamber S_r:S_i ratios less than 82% are reported as 'POSITIVE'; TEST chamber S_r:S_i ratios greater than or equal to 82% are reported as 'NEGATIVE'. *Numerical values are not reported.*

Fixed cut-off values for the CONTROL chamber are used to determine validity. CONTROL chamber S_r:S_i ratios less than 90% are considered valid and allow for reporting of TEST chamber results (POSITIVE, NEGATIVE). CONTROL chamber S_r:S_i ratios greater than or equal to 90% are considered invalid and prevent reporting of TEST chamber results. Invalid CONTROL chamber reactions are reported as 'INVALID'. *Numerical values are not reported.*

More stringent cut-off criteria are applied to the CONTROL chamber reaction to ensure amplification is not inhibited, reagents are performing as intended and that sample processing was performed appropriately.

BIOLOGICAL PRINCIPLES

Herpes simplex virus causes chronic, recurrent infections and viral meningitis in adults and children. Infection with herpes simplex virus type 1 (HSV-1) is relatively common, with an overall seroprevalence approaching 54% of the United States population.³ Herpes simplex virus Type 2 (HSV-2) is the primary cause of genital herpes. It is estimated that 50 million people in the US and >400 million worldwide are infected with HSV-2 genital herpes, or about 1 out of every 5 people.^{4,5} While either HSV-1 or HSV-2 is capable of oral or genital infections, HSV-1 is typically associated with oral lesions and HSV-2 is typically associated with genital lesions. However, the frequency of HSV-1 genital infection is increasing.^{3,6} HSV transmission to neonates can occur during birth, causing symptoms ranging from disseminated disease (CNS and visceral organs) to infection of the skin, eye, and mouth.⁷

HSV infection can be diagnosed using microscopy, cell culture, molecular, and serological methods;⁸ cell culture and PCR are the preferred diagnostic methods recommended by the Centers for Disease Control.⁴ The sensitivity of viral culture is low, particularly for recurrent or healing lesions, due to limited quantity of viable virus.⁴ The HSV type cannot be directly determined from viral culture, requiring additional typing methods (such as ELVIS).⁸ NAAT methods are more sensitive than culture, have the ability to differentiate between HSV types, and provide more rapid results.⁸

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- Alethia Sample Preparation Apparatus III (SMP PREP III) (item #2792-02):** Tris-buffered solution containing formalin-treated *E. coli* harboring *S. aureus* DNA and sodium azide (0.09%) as a preservative.
- Alethia HSV 1 Test Device (item #2865):** Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) and either HSV-1-specific primers (TEST chamber) or *S. aureus*-specific primers (CONTROL chamber).
- Alethia HSV 2 Test Device (item #2870):** Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) and either HSV-2-specific primers (TEST chamber) or *S. aureus*-specific primers (CONTROL chamber). **The Alethia HSV 2 Test Device is visually identified by an orange band on the Test Device closure tab.**
- Mineral Oil (item #14298-01)** (bottle with dropper tip)
- Alethia Heat Treatment (HT) Tubes (item #2335)**
- Transfer Pipettes (item #11523A)**

MATERIALS PROVIDED SEPARATELY

Alethia HSV 1&2 External Control Kit, Catalog Number: 479960

MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
2. DNase/RNase-free, aerosol resistant pipette tips
3. Viral transport medium (2-3 mL): MicroTest™ M4®, MicroTest™ M4RT®, MicroTest™ M6™, UniTranz-RT™ Transport System, Copan UTM™ Virus Collection, HealthLink® UTM or identical transport medium*, BD Universal Viral Transport Medium (UVT), ViraTrans™ VTM, Hardy Diagnostics VTM.
**HealthLink UTM is equivalent to Quest V-C-M Medium*

EQUIPMENT NOT PROVIDED

1. Dry-bath with 12 mm heat block capable of 95 C
2. Digital thermometer with max/min temperature memory (eg, Traceable® Lollipop™ Waterproof/Shockproof Thermometer)
3. Vortex mixer
4. Interval timer
5. Micropipette capable of dispensing 50 µL
6. Alethia reader, Meridian Bioscience, Inc. Catalog Number: 610189

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Viral transport media containing protein stabilizers, such as MicroTest™ M5® and calcium alginate swabs, are not acceptable for Alethia HSV 1&2 specimen collection and transport.
3. Do not interchange the Sample Preparation Apparatus or Test Devices between lots. Mineral Oil and Heat Treatment Tubes are interchangeable provided they are within the assigned expiration date when used.
4. Follow Biosafety Level 2 and Good Laboratory practices during testing.⁹ Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
5. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
6. Quality Control Programs for Molecular Testing Laboratories, including proper use and care of equipment, should be employed.¹⁰
7. The Alethia HSV 1&2 Test Devices contain lyophilized reagents. The protective pouches should not be opened until ready to perform the assay.
8. The Alethia HSV 1&2 Test Devices include a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
9. Dispose of used Alethia HSV 1&2 Test Devices immediately after processing, leaving the device latch securely in place. Do NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.
10. Any serious incident that has occurred in relation to the device should be reported to Meridian Bioscience, Inc., 3471 River Hills Drive, Cincinnati, Ohio 45244 USA or Technical Support Center 800.343.3858 and competent authority of the EU Member State in which the clinician and/or patient is established.
11. IMPORTANT: See SDS for additional safety and hazard information.

HAZARD AND PRECAUTIONARY STATEMENTS


There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-30 C.

PROCEDURAL NOTES

The transfer pipette included in the Alethia HSV 1&2 kit is individually wrapped and capable of dispensing 50 µL. The pipette is diagrammed below.

REAGENT PREPARATION

Ensure kit reagents are at room temperature (19-30 C) before use. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

SPECIMEN COLLECTION AND PREPARATION

Sample type: Human cutaneous and mucocutaneous lesion swab samples in viral transport media. The following viral transport media are acceptable for use: MicroTest™ M4®, MicroTest™ M4RT®, MicroTest™ M6™, UniTranz-RT™ Transport System, Copan UTM™ Virus Collection, HealthLink® UTM or identical transport medium*, BD Universal Viral Transport Medium (UVT), ViraTrans™ VTM, Hardy Diagnostics VTM.

**HealthLink UTM is equivalent to Quest V-C-M Medium*

NOTE: Viral transport media with protein stabilizers, ie, MicroTest™ M5® and calcium alginate swabs, are not acceptable sample types for use.

Sample Collection: Cutaneous and mucocutaneous lesion sample collection should be performed in accordance with institutional guidelines for collection of clinical specimens for herpes simplex virus infection.

Place swab(s) in a viral transport medium. Samples should be stored refrigerated (2-8 C) after collection and during transportation to the laboratory. Samples should be tested as soon as possible, but may be stored refrigerated (2-8 C) up to 7 days prior to testing. Do not freeze samples. Do not store at room temperature.

SPECIMEN PREPARATION

NOTE: Ensure that the Alethia reader is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the Alethia reader Operator's Manual for further information regarding instrument set-up and operation.

NOTE: Ensure specimens are at room temperature (19-30 C) before specimen preparation.

1. Vortex the collected swab in transport medium for approximately 45-60 seconds.
2. Using a transfer pipette, transfer 50 µL of the swab medium into one Alethia SMP PREP III. Discard the transfer pipette. Label the SMP PREP with the specimen identification.
3. Replace SMP PREP III cap and vortex for approximately 10 seconds.
4. Remove the tip cap from the SMP PREP III and squeeze the entire contents of the SMP PREP III into a Heat Treatment Tube. Label the tube with the specimen identification.
5. Repeat SPECIMEN PREPARATION steps for all samples to be processed. Processed samples may be held up to 2 hours at room temperature (19-30 C) prior to heat treatment.
6. Heat each Heat Treatment Tube containing sample in a dry-bath/heat block at 95 ± 5 C for 10 ± 2 minutes. Monitor heat-treatment step with digital thermometer and interval timer.
7. Remove each Heat Treatment Tube from the dry-bath/heat block. Heat-treated samples may be held at room temperature (19-30 C) for up to 15 minutes prior to testing.

NOTE: The amount of sample prepared at the end of SPECIMEN PREPARATION is enough to complete both Alethia HSV 1 and HSV 2 testing.

TEST PROCEDURE

NOTE: A maximum of 5 samples can be processed for both HSV 1 and HSV 2 in a single Alethia reader run. **HSV 2 Test Devices are identified by an orange band on the Test Device closure tab.**

1. Vortex heat-treated samples for approximately 10 seconds.
2. Use 1 Alethia HSV 1 Test Device and 1 Alethia HSV 2 Test Device per sample. Remove the devices from their protective pouches. Carefully open the device by lifting the latch, holding the chambers such that the lyophilized reagents will not fall out upon opening. Place devices on a flat surface or in a rack that can accommodate the devices.
3. Using a micropipette, transfer 50 µL of the sample to both the TEST (Left/White Bead) and CONTROL (Right/Yellow Bead) chambers of the Alethia HSV 1 Test Device and the HSV 2 Test Device. Take care to not introduce air to the reaction mixtures. Do not mix reactions with pipette.
4. Add 1 drop of Mineral Oil to both the TEST and CONTROL chambers of the Alethia HSV 1 and HSV 2 Test Devices. Close the Alethia Test Devices and fasten the latches securely.
5. Tap each device on the bench top or mix to remove air bubbles. Carefully examine each Test Device for rehydration of the Control/Test Beads, for air bubbles left in the chamber and liquid in the top of the device. If undissolved beads, air bubbles or liquid in the top of the device are noted, tap the device on the bench top and repeat visual inspection. Amplification and detection should be initiated within 15 minutes.
6. Repeat TEST PROCEDURE steps for all samples to be tested.
7. Insert the Alethia Test Devices into the Alethia reader and initiate run. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Patient Specimen	POSITIVE	HSV 1 TEST DEVICE: Sample contains HSV-1 target DNA. HSV 2 TEST DEVICE: Sample contains HSV-2 target DNA.
	NEGATIVE	HSV 1 TEST DEVICE: No HSV-1 target DNA detected. HSV 2 TEST DEVICE: No HSV-2 target DNA detected.
	INVALID	No reportable result. Repeat the test using the original sample. Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure.
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia reader performing correctly.
	NEGATIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia reader performing correctly.
	INVALID	No reportable result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia reader Well. OR The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. Each device contains an internal control that controls for amplification inhibition, assay reagents and sample processing effectiveness. Internal control DNA is present in the SMP PREP III and is processed through all steps of the procedure. Primers for amplification of internal control DNA are present in the Control Chamber of the Alethia Test Device.
2. The heat treatment step is monitored with an external thermometer and interval timer. Use the max/min temperature memory of the thermometer to ensure that a temperature of 95 ± 5 C is maintained. Use the interval timer to ensure that heat-treatment duration is 10 ± 2 minutes.
3. Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
4. Alethia HSV 1&2 External Control Reagents are supplied separately (Catalog 479960). It is recommended that the reactivity of each new lot and each new shipment of Alethia HSV 1&2 be verified on receipt or before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The Alethia HSV 1&2 test kit should not be used in patient testing if the external controls do not produce the correct results.
5. A separate Test Device must be used for each external control reagent.

EXPECTED VALUES

The prevalence of HSV-1 and HSV-2 was calculated during 2014-2015 clinical studies based on the age of patient and by lesion anatomical location. The overall incidence during the study was 20.5% (237/1155) for HSV-1 and 19.4% (224/1156) for HSV-2. Specific prevalence data is given below.

Combined Study - Cutaneous Prevalence by Sample Anatomical Location (N=306)						
Source	HSV 1			HSV 2		
	Total #	Total Positive	Prevalence	Total #	Total Positive	Prevalence
Genital – Penis	92	7	7.6%	92	28	30.4%
Skin Lesion	214	47	22.0%	214	27	12.6%
Combined Study - Mucocutaneous Prevalence by Sample Anatomical Location						
Source	HSV-1 (N=849)			HSV-2 (N=850)		
	Total #	Total Positive	Prevalence	Total #	Total Positive	Prevalence
Anorectal	47 (1*)	7	14.9%	46 (2*)	9	19.6%
Genital -- Vaginal/Cervical	624 (2*)	112	17.9%	626	158	25.2%
Nasal	18	9	50.0%	18	0	0.0%
Ocular	20	0	0.0%	20	0	0.0%
Oral Lesion	135	54	40.0%	135	2	1.5%
Urethral	5	1	20.0%	5	0	0.0%

* Number of samples producing invalid Alethia result, which could not be resolved and therefore, were excluded from the analysis.

Combined Study - Cutaneous Prevalence by Age						
Age	HSV-1 (N=306)			HSV-2 (N=306)		
	Total #	Total Positive	Prevalence	Total #	Total Positive	Prevalence
≤ 5 years	38	12	31.6%	38	1	2.6%
6 to 11 years	14	7	50.0%	14	1	7.1%
12 to 21 years	51	14	27.5%	51	4	7.8%
22 to 59 years	166	18	10.8%	166	36	21.7%
≥ 60 years	37	3	8.1%	37	13	35.1%
Not Provided	0	0	0.0%	0	0	0.0%
Combined Study - Mucocutaneous Prevalence by Age						
Age	HSV-1 (N=849)			HSV-2 (N=850)		
	Total #	Total Positive	Prevalence	Total #	Total Positive	Prevalence
≤ 5 years	47	8	17.0%	47	0	0.0%
6 to 11 years	12	0	0.0%	12	0	0.0%
12 to 21 years	174 (1*)	46	26.4%	175	42	24.0%
22 to 59 years	550 (1*)	111	20.2%	551	116	21.1%
≥ 60 years	63 (1*)	18	28.6%	63 (1*)	11	17.5%
Not Provided	3	0	0.0%	2 (1*)	0	0.0%

* Number of samples producing invalid Alethia result, which could not be resolved and therefore, were excluded from the analysis.

LIMITATIONS OF THE PROCEDURE

1. Performance of Alethia HSV 1&2 has been established with prospective, male and female cutaneous and mucocutaneous lesion specimens only. Performance with other specimen types (frozen specimens or CSF) has not been established.
2. The assay is not intended to be used for prenatal screening.
3. Medications containing the active ingredient Zincum Gluconicum produced invalid results when tested at 7% v/v with Alethia HSV 1&2. Lower concentrations were not tested. The effect of Zincum Gluconicum on Alethia HSV 1&2 results have not been determined using clinical specimens.
4. Casein at concentrations greater than 5 mg/mL was found to interfere with the Alethia HSV 1&2 assay.
5. Due to confirmed unacceptable results, refrigerated storage (2-8 C) will not be considered an acceptable storage condition for SMP PREP III processed specimens to be tested with the Alethia HSV 1&2 DNA Amplification Assay.
6. Negative results do not preclude infection and should not be used as the sole basis for treatment or other patient management decisions. This product can only be used with the Alethia instrument.
7. Alethia HSV 1&2 is a qualitative assay and does not provide quantitative values or information about organism load.
8. The detection of nucleic acids is dependent upon proper specimen collection, handling, transportation, storage and preparation. Failure to observe proper procedure in any one of these steps can lead to incorrect results.
9. Organism nucleic acid may persist *in vivo*, independent of organism viability. The Alethia HSV 1&2 does not distinguish between viable and nonviable organisms.
10. As with all molecular based diagnostic tests, (A) False-negative results may occur from the presence of inhibitors, mutations or polymorphisms in the HSV target region, technical error, sample mix-up or low numbers of organisms in the clinical specimen; (B) False-positive results may occur from the presence of crosscontamination by target organisms, their nucleic acids or amplified product, and from nonspecific signals.

SPECIFIC PERFORMANCE CHARACTERISTICS

The performance characteristics of the Alethia HSV 1&2 DNA Amplification Assay were established in clinical studies conducted between October 2014 and March 2015 at seven independent clinical test sites representing four geographically distinct regions throughout the United States. Performance characteristics of the assay were compared to the reference method ELVIS® HSV ID and D³ Typing Test System.

Samples included fresh cutaneous and mucocutaneous lesion swab specimens from patients suspected of herpes simplex virus type 1 (HSV-1) or type 2 (HSV-2) infection. Samples classified as cutaneous included skin and genital – penis lesions. Samples categorized as mucocutaneous included anorectal, genital-vaginal/cervical, nasal, ocular, oral, and urethral sites.

A total of 1158 eligible leftover, deidentified lesion swab specimens from symptomatic male and female patients were evaluated. Samples were collected from patients ranging in age from 1 day to 89 years. From the 1158 eligible samples, 1156 samples gave valid ELVIS results. The ELVIS HSV ID and D³ Typing Test System is unable to detect HSV-1 if HSV-2 has been identified first in coinfected specimens. Therefore, if a specimen was HSV-2 positive by ELVIS it was removed from the calculation of the HSV-1 performance analysis. One hundred and eighty-one (181) specimens were identified as HSV-2 positive by ELVIS and were removed from the HSV-1 performance analysis. In addition, one sample tested by the Alethia HSV 1&2 DNA Amplification Assay for HSV-1 gave an initially invalid result and could not be re-tested; this sample was excluded from the HSV-1 performance analysis. Thus, a total of 974 specimen results were included in the performance analysis for HSV-1. From the 1156 samples that gave valid ELVIS results, one sample tested by the Alethia HSV 1&2 DNA Amplification Assay for HSV-2 gave an initially invalid result and could not be re-tested; this sample was excluded from the HSV-2 performance analysis. Thus, a total of 1155 specimen results were included in the performance analysis for HSV-2. The tables summarize the performance of the Alethia HSV 1&2 DNA Amplification Assay for all sites combined.

A total of 723 (74.3%) female and 247 (25.3%) male specimens were tested for HSV-1 and a total of 873 (75.6%) female and 277 (24.0%) male specimens were tested for HSV-2. There were four (4) specimens tested for HSV-1 (0.4%) and five (5) tested for HSV-2 samples with unknown gender (0.4%). There is no expectation that assay performance is influenced by gender.

Combined Sites: HSV-1 Cutaneous (N=264)

		Reference Method			Alethia	Performance			
		Pos	Neg	Total					95% CI
Alethia HSV 1&2	Pos	48	6 ^a	54	0 (0)	Sensitivity	48/51	94.1%	84.1-98.0%
	Neg	3 ^b	207	210	0 (0)	Specificity	207/213	97.2%	94.0-98.7%
	Total	51	213	264	0 (0)				

^a 6/6 samples identified as HSV-1 positive by an alternative, FDA-cleared molecular assay.

^b 1/3 samples identified as HSV-1 negative by an alternative, FDA-cleared molecular assay.

^c Initial invalid results are reported within the parentheses. The final number of invalid samples remaining after repeat testing is shown before the parenthesis.

Combined Sites: HSV-1 Mucocutaneous (N=710)

		Reference Method			Alethia	Performance			
		Pos	Neg	Total					95% CI
Alethia HSV 1&2	Pos	152	28 ^b	180	0 (1)	Sensitivity	152/160	95.0%	90.5-97.5%
	Neg	8 ^c	522	530	0 (5)	Specificity	522/550	94.9%	92.7-96.5%
	Total	160	550	710	0 (6) ^a				

^a There were six initial INV samples by Alethia. Five repeated as Alethia negative (ELVIS negative); one repeated as Alethia HSV 1 positive (ELVIS HSV 1 positive).

^b 19/28 samples identified as HSV-1 positive by an alternative, FDA-cleared molecular assay; 3 samples could not be tested.

^c 7/8 samples were identified as HSV-1 negative by an alternative, FDA-cleared molecular assay; 1 sample could not be tested.

^d Initial invalid results are reported within the parentheses. The final number of invalid samples remaining after repeat testing is shown before the parenthesis.

Combined Sites: HSV-2 Cutaneous (N=306)

		Reference Method			Alethia	Performance			
		Pos	Neg	Total					95% CI
Alethia HSV 1&2	Pos	42	13 ^b	55	0 (0)	Sensitivity	42/42	100%	91.6-100.0%
	Neg	0	251	251	0 (1)	Specificity	251/264	95.1%	91.8-97.1%
	Total	42	264	306	0 (1) ^a				

^a There was one initial INV sample by Alethia. The sample repeated as Alethia negative (ELVIS negative).

^b 8/13 samples were identified as HSV-2 positive by an alternative, FDA-cleared molecular assay; 1 sample could not be tested.

^c Initial invalid results are reported within the parentheses. The final number of invalid samples remaining after repeat testing is shown before the parenthesis.

Combined Sites: HSV-2 Mucocutaneous (N=849)

		Reference Method			Alethia	Performance			
		Pos	Neg	Total					95% CI
Alethia HSV 1&2	Pos	137	31 ^b	168	0 (0)	Sensitivity	137/139	98.6%	94.9-99.6%
	Neg	2 ^c	679	681	0 (1)	Specificity	679/710	95.6%	93.9-96.9%
	Total	139	710	849	0 (1) ^a				

^a There was one initial INV sample by Alethia. The sample repeated as Alethia negative (ELVIS negative).

^b 24/31 samples were identified as HSV-2 positive by an alternative, FDA-cleared molecular assay; 4 samples could not be tested.

^c 1/2 samples were identified as HSV-2 negative by an alternative, FDA-cleared molecular assay.

^d Initial invalid results are reported within the parentheses. The final number of invalid samples remaining after repeat testing is shown before the parenthesis.

ANALYTICAL SENSITIVITY

Analytical sensitivity or Limit of Detection (LoD) was determined using two strains each of HSV-1 and HSV-2. The LoD was confirmed using 60 replicates and a stated probability (eg, 95%, where 57/60 replicates are positive) of obtaining positive responses.

HSV Strain	LoD (TCID ₅₀ /mL)
HSV-1 HF (VR-260)	7.20×10^3
HSV-1 McIntyre (VR-539)	9.89×10^4
HSV-2 G (VR-734)	1.20×10^3
HSV-2 MS (VR-540)	1.60×10^3

The final assigned assay LoD is 9.89×10^4 TCID₅₀/mL for HSV-1 and 1.60×10^3 TCID₅₀/mL for HSV-2.

In addition, twenty (20) HSV-1 and twenty (20) HSV-2 clinical samples from a variety of anatomical locations (ie, oral, genital, anorectal, etc.) were evaluated at the assay LoD for HSV-1 (9.89×10^4 TCID₅₀/mL) and HSV-2 (1.60×10^3 TCID₅₀/mL). All clinical samples were detected by the Alethia HSV 1&2 assay at LoD, or lower, except one HSV-1 sample, which was detected at 7.20×10^7 TCID₅₀/mL.

REPRODUCIBILITY

Reproducibility studies were carried out by three of the seven participating Clinical Sites. Blind-coded panels of eighteen (18) samples were supplied to participating laboratories. Samples were randomly sorted within each panel to mask sample identities. The panels included contrived HSV-1 (strain HF) and HSV-2 (strain MS) samples manufactured as moderate positive samples (HSV-1: 2.16×10^4 TCID₅₀/mL or HSV-2: 4.80×10^3 TCID₅₀/mL), low positive samples (HSV-1: 1.08×10^4 TCID₅₀/mL or HSV-2: 2.40×10^3 TCID₅₀/mL), or two HSV-contrived negative samples (HSV-1 near cut-off: 308 TCID₅₀/mL, HSV-1 high negative: 29.7 TCID₅₀/mL; HSV-2 near cut-off: 24 TCID₅₀/mL, HSV-2 high negative: 2.2 TCID₅₀/mL). The panel also included two HSV-1 and HSV-2 negative samples and a Positive and Negative Control. Testing was performed by different operators at each site on the same day (intra-assay variability) for five days (inter-assay variability). Three lots of Alethia HSV 1&2 and five Alethia instruments were used in this study. Positive and Negative Controls were tested with each panel. The results are provided in the table below:

Reproducibility Study Summary for HSV-1								
Sample Type	Site 1		Site 2		Site 3		Total	
	Percent Agreement		Percent Agreement		Percent Agreement		Percent Agreement	
HSV-1 Moderate Positive	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
HSV-1 Low Positive	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
HSV-1 & HSV-2 Near Cut-Off	20/30	66.7%	26/30	86.7%	20/30	66.7%	66/90	73.3%
HSV-1 & HSV-2 High Negative	29/30	96.7%	30/30	100.0%	29/30	96.7%	88/90	97.8%
HSV-1 Negative	60/60	100.0%	60/60	100.0%	60/60	100.0%	180/180	100.0%
Negative Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Positive Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%

Reproducibility Study Summary for HSV-2								
Sample Type	Site 1		Site 2		Site 3		Total	
	Percent Agreement		Percent Agreement		Percent Agreement		Percent Agreement	
HSV-2 Moderate Positive	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
HSV-2 Low Positive	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
HSV-1 & HSV-2 Near Cut-Off	25/30	83.3%	29/30	96.7%	25/30	83.3%	79/90	87.8%
HSV-1 & HSV-2 High Negative	29/30	96.7%	30/30	100.0%	29/30	96.7%	88/90	97.8%
HSV-2 Negative	60/60	100.0%	60/60	100.0%	60/60	100.0%	180/180	100.0%
Negative Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Positive Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%

CROSSREACTIVITY STUDIES

Crossreactivity studies employed positive (HSV-1 strain HF and HSV-2 strain MS) and negative specimens inoculated with bacterial or fungal organisms to a minimum concentration of 1.0×10^6 CFU/mL (1.0×10^6 copies/mL where genomic DNA was used) or virus at a minimum of 1.0×10^5 TCID₅₀/mL (1.0×10^6 copies/mL where genomic DNA or RNA was used). None of the following organisms or their genetic material reacted with Alethia HSV 1&2:

Acinetobacter calcoaceticus, *Acinetobacter lwoffii*, *Bacteroides fragilis*, *Bordetella bronchiseptica*, *Bordetella pertussis*, *Candida albicans*, *Candida glabrata*, *Candida guilliermondii*, *Candida krusei*, *Candida lusitanae*, *Candida parapsilosis*, *Candida tropicalis*, *Chlamydia trachomatis*, *Chlamydomydia pneumoniae*, *Clostridium difficile*, *Clostridium perfringens*, *Corynebacterium diphtheriae*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Enterococcus faecium*, *Escherichia coli* (ESBL), *Fusobacterium nucleatum*, *Gardnerella vaginalis*, *Haemophilus ducreyi*, *Haemophilus influenzae* (Type A), *Klebsiella pneumoniae*, *Lactobacillus acidophilus*, *Legionella pneumophila*, *Mobiluncus curtisii*, *Mobiluncus mulieris*, *Moraxella catarrhalis*, *Mycoplasma hominis*, *Mycoplasma orale*, *Mycoplasma pneumoniae*, *Mycoplasma salivarium*, *Neisseria gonorrhoeae*, *Neisseria meningitidis*, *Prevotella melaninogenica*, *Proteus mirabilis*, *Pseudomonas aeruginosa*, *Salmonella enteritidis*, *Salmonella typhimurium*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Staphylococcus saprophyticus*, *Streptococcus agalactiae*, *Streptococcus mitis*, *Streptococcus mutans*, *Streptococcus pneumoniae*, *Streptococcus pyogenes*, *Streptococcus salivarius*, *Toxoplasma gondii*, *Treponema palladium*, *Trichomonas vaginalis*, *Ureaplasma urealyticum*, Adenovirus, Coronavirus, Cocksackievirus, Cytomegalovirus, Echovirus, Enterovirus, Epstein Barr virus, Influenza A virus, Influenza B virus, Hepatitis B virus, Hepatitis C virus, Human herpes 6 virus, Human herpes 7 virus, Human herpes 8 virus, Human immunodeficiency virus type 1, Human metapneumovirus, Human papilloma virus, Measles virus, Mumps virus, Parainfluenza virus, Respiratory syncytial virus, Rubella virus, Varicella zoster virus.

Human genomic DNA was nonreactive at 1.0×10^6 copies/mL. In addition, there was no competitive inhibition observed between the organisms listed above and with HSV-1 or HSV-2 in the Alethia HSV 1&2 assay.

TESTS FOR INTERFERING SUBSTANCES

Interference testing was performed in the presence of worst-case concentrations of chemical and biological substances introduced directly into contrived positive (HSV-1 strains: McIntyre and HF; HSV-2 strains: G and MS) and negative samples. No interference was observed with the following substances:

- (1) At concentrations of 7% w/v or v/v: Abreva® (10% Docosanol), Balneol® Hygienic Cleansing Lotion, Carmex® Original Lip Balm (1.7% camphor, 0.7% menthol), Desitin® (40% zinc oxide), Douche (CVS Pharmacy® Disposable), Fluoride toothpaste (Crest®, 0.243% sodium fluoride), K-Y® Brand Jelly, Lanacane® (0.2% benzethonium chloride, 20% benzocaine), Lip Clear® Lysine+® (1.2% zinc oxide), Miconazole 3 (CVS™ Yeast Infection Relief: 2% miconazole nitrate), Mouthwash (Listerine® Original: 0.092% eucalyptol, 0.042% menthol, 0.060% methyl salicylate, 0.064% thymol), Preparation H® Hemorrhoidal Ointment (14% mineral oil, 74.9% petrolatum, 0.25% phenylephrine HCl), Releev® (0.13% benzalkonium chloride), Tioconazole, Vagisil® Regular Strength (5% benzocaine, 2% resorcinol), Yeast Guard® Gel Treatment (Candida albicans, 27X HPUS; Candida parapsilosis, 27X HPUS; Pulsatilla, 27X HPUS), Feces, Seminal fluid, Urine, Whole blood, Buffy coat.
- (2) At concentrations of: 60 µg/mL – Mucus (Mucin, bovine submaxillary gland type I-S); 1.25 mg/mL – Cornstarch; 3.3 mg/mL – Albumin; 5 mg/mL – Acetaminophen, Casein, Chlorpheniramine maleate; 7 mg/mL – Acyclovir; 10 mg/mL – Acetylsalicylic acid, Dextromethorphan hydrobromide.

Casein at concentrations greater than 5 mg/mL was found to interfere with the assay.

Cold-EEZE® Cold Remedy plus Sore Throat (zincum gluconicum 2X) was tested at 7% v/v and produced invalid results in all replicates.

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REFERENCES

1. Nagamine K, Hase T, Notoni T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. *Mol Cell Probes* 2002;16:223-29.
2. Mori Y, Kitao M, Tomita N, Notoni T. Real-time turbidimetry of LAMP reaction for quantitating template DNA. *J Biochem Biophys* 2004;59:145-47.
3. Bradely H, Markowitz L, Gibson T, McQuillan M. Seroprevalence of herpes simplex virus Types 1 and 2 – United States, 1999-2010. *J Infect Dis* 2014; 209:325-333.
4. Centers for Disease Control and Prevention, *Sexually transmitted diseases treatment guidelines, 2010*. Morbidity and Mortality Weekly Report, 2010: 59(RR-1220-25).
5. Looker KJ, Margaret AS, Turner KM, Vickerman P, Gottlieb SL, Newman LM. Global estimates of prevalent and incident herpes simplex virus type 2 infections in 2012. *PLoS ONE* 2015; 10(1): e114989. doi:10.1371/journal.pone.0114989.
6. Bernstein DI, Bellamy AR, Hook EW 3rd, Levin MJ, Wald A, Ewell MG, et al. Epidemiology, clinical presentation, and antibody response to primary infection with herpes simplex virus type 1 and type 2 in young women. *Clin Infect Dis* 2013; 56(3):344-351.
7. Kimberlin DW, Baley J. Guidance on management of asymptomatic neonates born to women with active genital herpes lesions. *American Academy of Pediatrics, Pediatrics* 2013;131(2):e635-e646.
8. Anderson NW, Buchan BW, Ledebore NA. Light microscopy, culture, molecular, and serologic methods for detection of herpes simplex virus. *J Clin Microbiol* 2014;52(1):2-8.
9. US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.
10. CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute, 2006.

SN11028

REV. 02/26

 Manufactured By	Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com
EC REP	Meridian Bioscience Europe, SRL Via Dell'Industria 7, 20035 Villa Cortese (Milano) ITALY www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+39) 0331.433636 E-mail: info@meridianbioscience.eu Technical Support: MBE-TechService@meridianbioscience.eu Customer Service/Orders: • For Italian Customers: ordini@meridianbioscience.com • For Distributors / International Customers: Export.CustomerService@meridianbioscience.eu
UK Authorised Representative	Launch Diagnostics Limited Lakeview West Crossways Business Park Galleon Boulevard Dartford Kent DA2 6QE United Kingdom
CH REP	MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/package of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	Meridian products carrying the European Conformity (CE) mark fulfill the requirements of Directive 98/79/EC or the Regulation 2017/746 on in-vitro diagnostic medical devices / I prodotti Meridian recanti il marchio di Conformità Europea (CE) soddisfano i requisiti della Direttiva 98/79/CE o del Regolamento 2017/746 sui dispositivi medico-diagnostici in vitro / Les produits Meridian portant la marque de Conformité européenne (CE) sont conformes aux exigences de la Directive 98/79/CE ou du Règlement 2017/746 relatifs aux dispositifs médicaux de diagnostic in vitro / Los productos de Meridian que llevan la marca de conformidad europea (CE) cumplen los requisitos de la Directiva 98/79/CE o del Reglamento 2017/746 sobre los productos sanitarios para diagnóstico in vitro / Produkte von Meridian mit der CE-Kennzeichnung erfüllen die Anforderungen der EU-Richtlinie 98/79/EG bzw. der Verordnung (EU) 2017/746 über In-vitro-Diagnostika	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenzubereitung, in dem sich Probenverdünnungspuffer befindet
		STERILE R	Sterilization by gamma irradiation / Sterilizzazione con raggi gamma / Stérilisation par irradiation aux rayons gamma / Esterilizado por irradiación gamma / Sterilisation durch Gammastrahlen
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	RoHS	Restriction of Hazardous Substances / Restrizione all'uso di sostanze pericolose / Limitation de substances dangereuses / Restricción de Substancias Nocivas / Beschränkung der Verwendung bestimmter gefährlicher Stoffe
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		Cautions, consult accompanying documents / Attenzione, vedere le istruzioni per l'uso / Attention voir notice d'instructions / Atención, ver instrucciones de uso / Achtung, Begleitdokumente beachten
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	STERILE EO	Sterilization by ethylene oxide / Sterilizzazione con ossido di etilene / Stérilisation par oxyde d'éthylène / Esterilizado por óxido de etileno / Sterilisation durch Ethylenoxid
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> saggi / Contenu suffisant pour <n> tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung		ETL Registered Mark Certified / Marchio di certificazione registrato a livello nazionale / Certifié Conforme ETL / Marca de Certificación Registrada Nacional / ETL Konform beglaubigt
	Female / Femminile / De sexe féminin / Hembra / Frau		Recycle – do not dispose of as general waste / Riciclare – non eliminare come rifiuto generico / Recycler – ne pas jeter dans une poubelle / Recicle – no desecho como basura general / Recycling- dieses Produkt nicht über den Hausmüll entsorgen
	Single Use Only / Prodotto Monouso / A usage unique / Para Un Solo Uso / nur für die einmalige Anwendung	TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	HT TUBE	Heat Treatment Tube / Provetta per il Trattamento termico / Tube pour le traitement thermique / Tubo de tratamiento de calor / Röhrrchen zur Hitzebearbeitung
	LASER RADIATION: Avoid Exposure to Beam / RADIAZIONE LASER: Evitare l'esposizione al raggio / RAYONNEMENT LASER: Eviter toute exposition au faisceau / Radiación Laser: Evite Exposición a los Rayos / LASERSTRAHLUNG: Direkten Kontakt mit dem Strahl vermeiden		HOT SURFACE: Keep hands away from Hot Surfaces / Superficie calda: tenere le mani lontane dalle superfici calde / SURFACES CHAUDES: Ne pas toucher les surfaces chaudes / Superficie Caliente: Mantenga las manos alejadas de la superficie caliente / Heiße Oberfläche: Kontakt mit heißen Oberflächen vermeiden
	CAUTION: Laser Radiation / ATTENZIONE: Radiazione Laser / AVERTISSEMENT: Rayonnement laser / Precaución: Radiación Laser / WARNUNG: Laserstrahlung	IPX-0	CAUTION: Protect from water / ATTENZIONE: Proteggere dall'acqua / AVERTISSEMENT: Protéger de l'humidité / Precaución: Proteja del agua / WARNUNG: Vor Feuchtigkeit schützen
	CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Risikogefahr	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
MIN OIL	Mineral Oil / Olio Minerale / Huile Minérale / Aceite Mineral / Mineralöl	MEDIA	Media / Terreno di trasporto / Milieux / Medio / Medium
ST TUBE	Screw Top Tube / Provetta con tappo a vite / Tube à bouchon vissé / Tubo con tapa de rosca / Röhrrchen mit Schnappverschluss	COL	Sample Preparation Column / Colonna di preparazione del campione / Colonne pour la préparation de l'échantillon / Columna de preparación de muestra / Säule zur Probenaufarbeitung
BUF SMP	Sample Buffer / Soluzione tampone per il campione / Tampón de l'échantillon / Tampón de muestra / Probenpuffer	PRE REAG	Pretreatment Reagent / Reagente di Pretrattamento / Réactif de prétraitement / Reactivo de pretratamiento / Reagenz für die Vorbehandlung
Rx Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	SMP PREP	Sample Preparation / Preparazione del campione / Préparation de l'échantillon / Preparación de Muestra / Probenzubereitung
TUBE	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß	CH REP	Swiss Authorized Representative / Mandatario svizzero / Mandataire Suisse / Representante Autorizado Suizo / Schweizer Bevollmächtigter

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.