premier®

C. difficile GDH

Enzyme immunoassay for the detection of
Clostridium difficile GDH in stool specimens

REF 611096

IVD

Rx Only

INTENDED USE

Premier *C. difficile* GDH is a qualitative enzyme immunoassay screening test to detect *Clostridium difficile* antigen, glutamate dehydrogenase, in fecal specimens from symptomatic persons suspected of having *C. difficile* infection (CDI). This test does not distinguish between toxigenic and non-toxigenic strains of *C. difficile*. Samples from symptomatic patients that produce positive results with this test must be further tested with an assay designed to detect toxigenic *C. difficile* strains and assist with the diagnosis of CDI.

SUMMARY AND EXPLANATION OF THE TEST

Clostridium difficile has become the most important etiologic agent in nosocomial diarrhea worldwide. *C. difficile* infection (CDI) can range from mild diarrhea to life-threatening toxic megacolon. Classic risk factors for nosocomial CDI include antibiotic use and old age. Recently however CDI has been seen more frequently outside the hospital environment without the typical risk factors historically associated with *C. difficile* infection.¹ *C. difficile* produces two toxins (toxin A and toxin B) that have both been associated with disease caused by *C. difficile*. Strains of *C. difficile* that do not harbor the toxin A and toxin B genes are considered to be non-pathogenic. All strains of *C. difficile* produce glutamate dehydrogenase (GDH), also known as common antigen.^{2,3} Non-molecular laboratory approaches to *C. difficile* detection can be divided into two general groups: toxin-based and non-toxin based. Toxin-based assays include enzyme-immunoassay (EIA) for toxin A or toxins A and B and cytotoxin, which primarily detects toxin B. Non-toxin based assays, such as EIA for *C. difficile* GDH, while inherently more sensitive, must be followed up with testing for toxigenic *C. difficile*.

BIOLOGICAL PRINCIPLES

Premier *C. difficile* GDH is a qualitative enzyme immunoassay screening test to detect *Clostridium difficile* antigen, glutamate dehydrogenase, in fecal specimens from symptomatic persons suspected of having *C. difficile* infection (CDI). Breakaway microwells are coated with rabbit anti-GDH antibodies. Diluted patient specimen and a horseradish peroxidase conjugated anti-GDH antibody are added to the microwells and incubated. Upon completion of the incubation, a wash step is performed to remove unbound material. If *C. difficile* GDH is present, an antibody-enzyme complex is formed. A chromagenic substrate is added to the microwells and incubated. A blue color develops in the presence of bound enzyme. Stop solution is added, changing the initial blue reaction to yellow. Test results are interpreted spectrophotometrically.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Premier C. difficile GDH Microwells:** Polyclonal antibody-coated microwells. The antibodies are specific to *C. difficile* GDH.
2. **Premier C. difficile GDH Enzyme Conjugate:** HRP-conjugated polyclonal antibodies specific to *C. difficile* GDH in a buffered protein solution containing 0.1% ProClin® and 0.03% gentamicin as preservatives.
3. **Premier 20X Wash Buffer II:** Concentrated wash buffer containing 0.2% thimerosal as a preservative. Reagent is diluted before use.
4. **Premier Substrate I:** Buffered solution containing urea peroxide and tetramethylbenzidine at pH 5.0.
5. **Premier Stop Solution I:** 1M Phosphoric acid.
6. **Premier C. difficile GDH Sample Diluent/Negative Control:** Buffered protein solution containing 0.1% ProClin® and 0.03% gentamicin as preservatives.
7. **Premier C. difficile GDH Positive Control:** *C. difficile* GDH in a buffered protein solution containing 0.1% ProClin® and 0.03% gentamicin as preservatives.
8. Transfer pipettes
9. Microwell strip holder
10. Microwell plate sealers

MATERIALS NOT PROVIDED

1. Disposable latex gloves
2. Test tubes (eg, 10 x 75 mm or 12 x 75 mm) for dilution of sample
3. Distilled or deionized water
4. Squir bottle
5. Graduated cylinder for making 1X Wash Buffer
6. Absorbent paper
7. Waste container with disinfectant and/or autoclavable biohazard bags
8. Vortex mixer
9. Interval timer
10. EIA microplate reader capable of reading absorbance at 450 nm or 450/630 nm*
11. StatFax™-2200 Incubator/Shaker (StatFax™ is a trademark of Awareness Technology, Inc.) (optional)*
12. Semiautomated microplate washer (optional)*

* Note: It is the operator's responsibility to validate StatFax™, semiautomated plate washers and readers prior to their use with this product.

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Do not interchange the Microwells, Enzyme Conjugate, Premier Substrate I, or Positive Control reagents between lots. (The Sample Diluent/Negative Control, Premier 20X Wash Buffer II and Premier Stop Solution I are interchangeable provided the reagents are within their assigned expiration dates when used.)
3. Do not use kit components beyond labeled expiration date.
4. Do not use vials that lack a label, a lot number, or an expiration date.
5. Do not use any reagent if it is discolored or turbid. Discoloration or turbidity may be a sign of microbial contamination.
6. Allow reagents to warm to 21–27 °C before use.
7. All reagents should be mixed gently before use.
8. Hold reagent vials vertically at suitable distance above the well to insure proper drop size and delivery.
9. Replace colored caps on correct vials.
10. Do not reuse microwells.
11. Unused microwells must be placed back inside resealable pouch. It is important to protect strips from moisture.
12. The transfer pipettes provided with this kit must be used for specimen preparation and transfer. Use one per specimen.
13. Avoid splashing when dispensing diluted stool into microwells by placing the transfer pipette tip about halfway down the well and dispensing slowly down the side of well.
14. Microwell washing is to be performed precisely as directed in the assay procedure. **Inadequate washing may be the cause of elevated background in any EIA protocol.**
15. All reagents except the Premier 20X Wash Buffer II are provided already diluted to the proper concentration.
16. Any deviation below or above set incubation times may affect sensitivity and specificity and should be avoided.
17. Stool must be mixed thoroughly (regardless of consistency) to ensure a representative sample prior to pipetting.

WARNINGS

1. Patient specimens may contain infectious agents and should be handled and disposed of as potential biohazards.
2. Dispose of used wash buffer and all test materials in an appropriate container. Treat waste as a potential biohazard.
3. The Positive Control reagent contains *C. difficile* GDH. It should be handled as a potential biohazard.
4. Avoid skin contact with Premier Stop Solution I. Flush with water immediately if contact occurs.

HAZARD and PRECAUTIONARY STATEMENTS


 Premier Stop Solution I	Signal Word Danger Hazard Statements H412 - Harmful to aquatic life with long lasting effects H314 - Causes severe skin burns and eye damage Precautionary Statements - EU (§28, 1272/2008) P260 - Do not breathe dust/fume/gas/mist/vapors/spray P264 - Wash face, hands and any exposed skin thoroughly after handling P280 - Wear protective gloves/ protective clothing/ eye protection/ face protection P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing P303 + P361 + P353 - IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/shower P363 - Wash contaminated clothing before reuse P304 + P340 - IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing P310 - Immediately call a POISON CENTER or doctor/ physician P301 + P330 + P331 - IF SWALLOWED: Rinse mouth. Do NOT induce vomiting P405 - Store locked up P501 - Dispose of contents/ container to an approved waste disposal plant.
 Premier 20X Wash Buffer II	Signal Word Danger Hazard Statements H302 – Harmful if swallowed H311 – Toxic in contact with skin Precautionary Statements - EU (§28, 1272/2008) P280 – Wear protective gloves/protective clothing
 Premier C. difficile GDH Sample Diluent / Negative Control Premier C. difficile GDH Enzyme Conjugate Premier C. difficile GDH Positive Control	Signal Word Danger Hazard Statements H314 – Causes severe skin burns and eye damage H317 – May cause an allergic skin reaction H319 – Causes serious eye irritation Contains Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl-with 2-methyl-3(2H)-isothiazolone Precautionary Statements - EU (§28, 1272/2008) P280 – Wear protective gloves/protective clothing/ eye protection/ face protection P260 – Do not breathe dust/fume/gas/mist/vapors/spray P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P310 – Immediately call a POISON CENTER or doctor/physician P303 + P361 + P353 – IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/shower P280 – Wear eye protection/face protection

SHELF LIFE AND STORAGE

The kit expiration date is indicated on the kit label. Store the kit at 2–8 C and return the kit promptly to the refrigerator after each use.

PROCEDURAL NOTES

The Premier C. *difficile* GDH transfer pipette is diagrammed below.

REAGENT PREPARATION

- Bring the entire kit, including microwell pouch, to 21–27 C before use.
- Prepare 1X Wash Buffer as needed. For example: 10.0 mL of Premier 20X Wash Buffer II + 190.0 mL distilled or deionized water is sufficient to wash one strip. 1X Wash Buffer can be stored at 21–27 C for up to three months.

SPECIMEN COLLECTION AND PREPARATION

Human stool samples, unpreserved: Samples should be received in an airtight transport container and stored at 2–8 C prior to testing. Samples should be tested as soon as possible, but may be held up to five days at 2–8 C. (See SPECIMEN PREPARATION section for instructions on diluting samples.) Samples that will not be tested within five days should be frozen immediately upon receipt and stored at ≤ –20 C until tested. Specimens may be held frozen for up to two months and may be frozen and thawed twice.

SPECIMEN PREPARATION

Mix stool as thoroughly as possible prior to pipetting.

- Human stool samples, unpreserved:**
 - With the dropper assembly (or equivalent), add 200 µL of Sample Diluent/Negative Control to a small test tube.
 - Liquid/Semi-Solid stools:**
 - Using a transfer pipette provided with the kit, add 50 µL (first mark from pipette tip) of thoroughly mixed stool to the Sample Diluent/Negative Control test tube.
 - Vortex the mixture for a minimum of 15 seconds.
 - Save the transfer pipette in the test tube for future use.
 - Stool diluted in Sample Diluent:** If needed, stool diluted in Sample Diluent can be held at 2–27 C for up to 24 hours before testing provided the tube is sealed.
 - Stool specimens may be centrifuged after dilutions. Centrifuge at approximately 2750 x G for 3 minutes or until solid matter separates from liquid. Proceed with the assay after recovering the supernatant.

TEST PROCEDURE

1. After the pouch has reached temperature (21–27 C), break off the required number of microwells (1 well for each specimen, plus 1 positive and 1 negative control well per batch). Place the microwells in the microwell strip holder and record the location of all wells. Unused microwells must be resealed in the pouch immediately.
2. Using the specimen transfer pipette, add 100 µL of diluted stool (second mark from the tip of the pipette) to the appropriate well. (Place the pipette tip halfway into well and let the sample slowly run down side of well.)
3. Add 2 free falling drops of Positive Control. With a transfer pipette, add 100 µL (second mark from the tip of the pipette) of Sample Diluent/Negative Control to the appropriate wells.
4. Add 1 free falling drop (approximately 50 µL) of Enzyme Conjugate to each well. Firmly shake/swirl the plate for 30 seconds.
5. Cut the plate sealer to size and press it firmly onto the top of the microwells to seal. Incubate the plate for 50 minutes at 35–39 C. **Alternatively, laboratories equipped with a heated plate shaker (StatFax™-2200) can incubate and rotate the plate for 20 minutes at 37 C at 1000 rpm (setting #5).**
6. Carefully remove the plate sealer and wash the wells:
 - a. **Manual method:**
 - i. Dump the plate contents firmly into a biohazard receptacle.
 - ii. Invert and bang the plate on a clean stack of paper towels.
 - iii. Fill all wells with 1X Wash Buffer, directing a stream of buffer to the sides of the wells to avoid foaming.
 - iv. Repeat the wash cycle (dump, bang on fresh towels, fill) 4 to 6 more times for a total of 5-7 wash cycles. After the last fill, dump and bang plate on fresh towels hard enough to remove as much excess wash buffer as possible, but do not allow wells to completely dry at any time.
 - b. **Semi-automated wash method:**
 - i. Dump the plate contents firmly into a biohazard receptacle.
 - ii. Invert and bang the plate on a clean stack of paper towels, then place empty plate on washer device.
 - iii. Fill the wells to the top (approx. 300-350 µL/well) with 1X Wash Buffer and then aspirate. The washer manifold should be adjusted to ensure no foaming occurs during the filling of the wells and that the wells are thoroughly aspirated after each wash.
 - iv. Repeat step iii 4 more times. Following the last wash, test wells should be thoroughly aspirated to remove as much moisture as possible, but do not allow wells to completely dry at any time.
7. Add 2 free-falling drops (approximately 100 µL) of Premier Substrate I to each well. Firmly shake/swirl the plate for 30 seconds. Incubate the plate for 10 minutes at 21–27 C.
8. Add 2 free-falling drops (approximately 100 µL) of Premier Stop Solution I to each well. Firmly shake/swirl the plate for 30 seconds. **Note:** The initial color of a positive reaction is blue, which changes to yellow upon addition of Premier Stop Solution I.
9. Clean the underside of all wells with a lint-free tissue.
10. Inspect and record reactions. Test results should be read using a spectrophotometric reader. Zero EIA reader on air. Read absorbance at 450 nm or 450/630 nm within 30 minutes of adding Premier Stop Solution I.
11. Suspect failure of the washing method or device if the Negative and/or Positive Controls consistently produce out-of-specification results. Increasing the number of washes, washing more vigorously, decanting more thoroughly or recalibrating washing devices should correct the problem. If problems persist, contact Meridian's Technical Support Services at 800-343-3858 (USA) or your authorized Meridian distributor for assistance.

INTERPRETATION OF RESULTS

Spectrophotometric Single Wavelength (450 nm)

Negative:	< 0.200
Positive:	≥ 0.200
Negative Control:	< 0.150
Positive Control:	≥ 0.600

Spectrophotometric Dual Wavelength (450/630 nm)

Negative:	< 0.150
Positive:	≥ 0.150
Negative Control:	< 0.100
Positive Control:	≥ 0.600

If a Negative Control is < 0.000, reblank the plate reader to air and reread the plate.

A positive result indicates the presence of *C. difficile* GDH. A negative result indicates the absence of *C. difficile* GDH, or that the level of GDH is below what can be detected by the assay. The magnitude of the OD above the cutoff is neither indicative of the severity or extent of *C. difficile* infection, nor can it be correlated to an endpoint titer. Extremely strong positive samples may yield either an intense yellow color or a purple precipitate within a few minutes of stopping the reaction. In this case, the spectrophotometer may yield an "out" reading. This reading is considered a positive.

If the frequency of low positive results (OD between 0.200 and 0.250 for single wavelength and between 0.150 and 0.200 with dual wavelength) is greater than 5% of specimens tested, this may indicate insufficient washing. More vigorous washing or increasing the washes to seven washes in step 6 of the procedure is recommended.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing or leakage. Do not use contaminated or suspect reagents.
2. The performance of Premier *C. difficile* GDH should be verified using the Positive Control and Sample Diluent/Negative Control with each test batch. See the INTERPRETATION OF RESULTS section above for a description of the expected results for control reagents. Tests should be considered invalid when either control reagent does not produce the expected results. In such cases, repeat tests and controls. If, on repeat testing, the expected reactions are still not observed with in-date reagents, call Meridian's Technical Support Services at 800-343-3858 (USA) or your authorized Meridian distributor for assistance.
3. The controls are used to monitor reagent reactivity. Failure of the controls to produce the expected results can mean that one or more of the reagents are no longer reactive at the time of use, the test was not performed correctly, or that reagents or samples were not added. **If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.**
4. Specimen matrix interference has not been observed in this assay as samples are significantly diluted before testing in Sample Diluent. For this reason, the Positive and Negative Controls supplied as part of this assay are prepared in matrices similar to the Sample Diluent. If control materials that are identical in composition to test samples are preferred, the user can prepare these by diluting known positive and negative specimens in Sample Diluent according to the SPECIMEN PREPARATION section of this insert. Add 100 µL of user prepared controls to test wells.

EXPECTED VALUES

The frequency of antibiotic-associated diarrhea caused by *C. difficile* is dependent on several factors including: patient population, type of institution and epidemiology. The reported incidence of *C. difficile* infection in patients suspected of having antibiotic-associated diarrhea is 15-25%⁴ although different facilities may find positivity rates outside this range.

LIMITATIONS OF THE PROCEDURE

1. For in vitro diagnostic use only.
2. The test is qualitative and no quantitative interpretation should be made with respect to the values.
3. A positive result does not define the existence of *C. difficile* disease or the presence of toxigenic *C. difficile*, but only indicates the detection of the organism. All Premier *C. difficile* GDH positive stools should be tested to verify the presence of toxigenic *C. difficile*.
4. Overincubation of the tests may lead to an increase in false-positive test results. Conversely, incubation for periods less than those defined in this insert can result in an increase in false-negative tests. Follow incubation times defined in this insert.
5. *Staphylococcus aureus* (Cowan strain I) and *Clostridium sporogenes* were found to be crossreactive with Premier *C. difficile* GDH.
6. This test is for use with unformed unpreserved stool samples only. Performance characteristics of other clinical specimen types have not been established.

SPECIFIC PERFORMANCE CHARACTERISTICS

Premier *C. difficile* GDH was evaluated by five independent test sites located in the Southwestern, Southeastern, and Midwestern regions of the United States. A total of 733 qualified patient samples were evaluated; all specimens were prospectively collected samples. Samples were collected from males (46%) and females (53%). In the case of 1% of the patients, the gender was not known. The age groups of the patients ranged from 22 days to 99 years of age. No differences in test performance were observed based on patient age or gender. The following tables show the assay performance by clinical site and patient age.

Table 1. Performance characteristics by clinical site

Site	Positive Samples			Negative Samples		
	Premier GDH/ Culture	Sensitivity %	95% CI	Premier GDH/ Culture	Specificity %	95% CI
Total Sites	108/117	92.3%	86.0 – 95.9%	590/616	95.8%	93.9 – 97.1%
Site 1	13/16	81.3%	57.0 – 93.4%	84/87	96.6%	90.3 – 98.8%
Site 2	28/30	93.3%	78.7 – 98.2%	132/140	94.3%	89.1 – 97.1%
Site 3	44/46	95.7%	85.5 – 98.8%	147/153	96.1%	91.7 – 98.2%
Site 4	15/15	100.0%	79.6 – 100.0%	169/175	96.6%	92.7 – 98.4%
Site 5	8/10	80.0%	49.0 – 94.3%	58/61	95.1%	86.5 – 98.4%

- Discrepant specimens were evaluated using an FDA-cleared ELISA test for the detection of *C. difficile* GDH.
- Sixteen of the 26 false positive specimens were positive when tested with another FDA cleared GDH assay.
- Eight of the nine false negative specimens were negative when tested with another FDA cleared GDH assay.

Table 2. Performance characteristics by patient age

Patient Age	Positive Samples			Negative Samples		
	Premier GDH/ Culture	Sensitivity %	95% CI	Premier GDH/ Culture	Specificity %	95% CI
0-28 days	0/0	N/A	N/A	1/1	100.0%	20.7 – 100.0%
29 days to 2 years	22/24	91.7%	74.2 – 97.7%	56/60	93.3%	84.1 – 97.4%
> 2 years to < 12 years	24/25	96.0%	80.5 – 99.3%	106/109	97.2%	92.2 – 99.1%
12 years to < 18 years	9/11	81.8%	52.3 – 94.9%	59/61	96.7%	88.8 – 99.1%
18 years to 21 years	3/4	75.0%	30.1 – 95.4%	22/23	95.7%	79.0 – 99.2%
> 21 years	50/53	94.3%	84.6 – 98.1%	345/361	95.6%	92.9 – 97.3%
Not Defined	0/0	N/A	N/A	1/1	100.0%	20.7 – 100.0%

REPRODUCIBILITY

Assay precision, intra-assay variability and inter-assay variability were assessed with a reference panel prepared from pools of negative samples spiked with *C. difficile* GDH antigen. The reproducibility panel consisted moderately positive (n=3), low positive (n=3), high negative (n=3), and negative (n=1) specimens. The low positive and high negative specimens were prepared near the assay limit of detection. Each sample was evaluated twice per day for five days by three different laboratories. Reproducibility was 100% with no intra-assay and inter-assay variability for samples prepared above or below the limit of detection for the assay.

CROSSREACTIVITY STUDIES

Crossreactivity studies were performed with positive and negative stool specimens inoculated with bacterial or fungal organisms to a final concentration of 1.2×10^8 CFU/mL or virus concentration greater than 1×10^5 TCID₅₀/mL. None of the following organisms in stool reacted with Premier *C. difficile* GDH:

Aeromonas hydrophila, *Bacillus cereus*, *Bacillus subtilis*, *Bacteroides fragilis*, *Campylobacter coli*, *Campylobacter fetus*, *Campylobacter jejuni*, *Candida albicans*, *Citrobacter freundii*, *Clostridium butyricum*, *Clostridium bifementans*, *Clostridium histolyticum*, *Clostridium novyi*, *Clostridium perfringens*, *Clostridium septicum*, *Clostridium tetani*, *Enterobacter aerogenes*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Escherichia coli* O157:H7, *Escherichia hermannii*, *Escherichia fergusonii*, *Helicobacter pylori*, *Klebsiella pneumoniae*, *Lactococcus lactis*, *Listeria monocytogenes*, *Peptostreptococcus anaerobius*, *Plesiomonas shigelloides*, *Porphyromonas asaccharolytica*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Pseudomonas fluorescens*, *Salmonella* Group B, *Salmonella* Group C, *Salmonella* Group D, *Salmonella* Group E, *Serratia liquifaciens*, *Serratia marcescens*, *Shigella boydii*, *Shigella flexneri*, *Shigella sonnei*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Vibrio parahaemolyticus*, *Yersinia enterocolitica*, Adenovirus Type 40, Adenovirus Type 41, Coxsackievirus Strain B4, Echovirus Strain 11, Rotavirus Strain WA

Stools spiked with *Staphylococcus aureus* (Cowan strain I) and *Clostridium sporogenes*, were found to be crossreactive with Premier *C. difficile* GDH.

ANALYTICAL SENSITIVITY

The analytical sensitivity of the Premier *C. difficile* GDH test is 8 ng/mL when spiked in stool.

ASSAY REACTIVITY

The following *C. difficile* stock cultures from different sources were tested and produced positive reactions at 5.7×10^7 CFU/mL with Premier *C. difficile* GDH:

Toxigenic *C. difficile* strains: 8864, 10463, 43598, 2004052, 2004111, 2004118, 2004205, 2004206, 2005070, 2005257, 2005325, 2006240, 2007431, 2007435, 2007858, 2008016, 2008029, 2008162, 2008188, 2008341, 2008351, 2009018, 2009065, 2009066, 2009099, 2009132, 2009155, 2009277, B1, B17, B18, BK6, CF1, G1, J7, K12, Y1

Non-toxigenic *C. difficile* strains: 11186, 234, 586, 611, 620, 2C62, 2C165, C122, UNC19904, X15076

TESTS FOR INTERFERING SUBSTANCES

The following substances, at the specified saturated solvent/diluent concentrations, do not interfere with test results in the final concentrations listed:

Barium sulfate at 5 mg/mL, Stearic acid/Palmitic acid (fecal fat) at 2.65mg/mL/1.3 mg/mL, Hemoglobin at 3.2 mg/mL, Imodium AD® (Loperamide HCl) at 6.67×10^{-3} mg/mL, Kaopectate® (Bismuth subsalicylate) at 0.87 mg/mL, Metronidazole at 12.5 mg/mL, Mucin at 3.33 mg/mL, Mylanta® (Ammonium hydroxide w/ magnesium hydroxide) at 4.2 mg/mL, Pepto Bismol® (Bismuth subsalicylate) at 0.87 mg/mL, Polyethylene glycol 3350 at 79.05 mg/mL, Prilosec® (Omeprazole) at 0.5 mg/mL, Simethicone at 0.625 mg/mL, Tagamet® (Cimetidine) at 0.5 mg/mL, Tums® (Calcium carbonate) at 0.5 mg/mL, Vancomycin hydrochloride at 2.5 mg/mL, White blood cells at 5%, Whole blood at 25%

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REFERENCES

1. Dawson LF, Valiente E, Wren B. *Clostridium difficile* - A continually evolving and problematic pathogen. Infect Genet Evol 2009;9:1410-1417.
2. Lysterly DM, Barroso LA, Wilkins TD. Identification of the latex-reactive protein of *Clostridium difficile* as glutamate dehydrogenase. J Clin Micro 1991;29:2639-2642.
3. Willis, DH, Kraft JA. Confirmation that the latex-reactive protein of *Clostridium difficile* is a glutamate dehydrogenase. J Clin Micro 1992;30:1363-1364.
4. Bartlett JG, Gerding D. Clinical Recognition and Diagnosis of *Clostridium difficile* infection. Clin Inf Dise 2008;46:S12-18.

SN11082

REV. 01/24

 Manufactured By	Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com
	Meridian Bioscience Europe, SRL Via Dell'Industria 7, 20035 Villa Cortese (Milano) ITALY www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+39) 0331.433636 E-mail: info@meridianbioscience.eu Technical Support: MBE-TechService@meridianbioscience.eu Customer Service/Orders: - For Italian Customers: ordini@meridianbioscience.com - For Distributors / International Customers: Export.CustomerService@meridianbioscience.eu
AUSTRALIAN SPONSOR	Emergo Australia Level 20, Tower II Darling Park 201 Sussex Street Sydney, NSW 2000 Australia

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/package of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.