ImmunoCard[®] STAT_® Mono Test

1

1 DROP
1 GOCCIA
1 GOUTTE
1 GOTA
1 TROPFEN

For serum, plasma,
or whole blood
Samples in tubes:

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-

ImmunoCard[®] STAT! Mono Test

CLIA Complexity:
Whole Blood – Waived
Serum or Plasma – Non-Waived

REF 755725

Rx Only

IVD In Vitro Diagnostic Medical Device

FOR LABORATORY AND PROFESSIONAL IN VITRO DIAGNOSTIC USE ONLY.

INTENDED USE

The ImmunoCard STAT! Mono Test is intended for the qualitative detection of infectious mononucleosis heterophile antibodies in serum, plasma or whole blood as an aid in the diagnosis of infectious mononucleosis.

SUMMARY AND EXPLANATION OF TEST

The diagnosis of infectious mononucleosis (IM) is suggested on the basis of the clinical symptoms of fever, sore throat and swollen lymph glands. The highest incidence of symptomatic IM occurs during late adolescence (15-24 years of age). Infectious mononucleosis is caused by the Epstein-Barr Virus (EBV).^{1,2} The laboratory diagnosis of IM is based on the detection of IM heterophile antibodies. These heterophile antibodies are directed against antigens found in bovine, sheep and horse erythrocytes. The ImmunoCard STAT! Mono Test utilizes an extract of bovine erythrocytes to give the required sensitivity and specificity.

PRINCIPLES OF TEST

The ImmunoCard STAT! Mono Test uses color immunochromatographic dipstick technology with bovine erythrocyte extract coated on the membrane. In the test procedure, serum, plasma or whole blood is mixed with the Diluent. Then the Test Stick is placed in the mixture and the mixture migrates along the membrane. If the specific IM heterophile antibody is present in the sample, it will form a complex with the bovine erythrocyte extract conjugated color particles. The complex will then be bound by bovine erythrocyte extract immobilized on the membrane and a visible blue Test Line will appear to indicate a positive result.

KIT CONTENTS AND STORAGE

- 25 Test Sticks
- 25 Test Tubes
- 25 Transfer Pipettes
- 27 Capillary Pipettes
- 1 Diluent (contains buffer with 0.2% sodium azide)
- 1 Mono Positive Control (contains rabbit anti-beef stroma in tris buffer with 0.2% sodium azide and 0.05% gentamycin sulfate preservatives)
- 1 Mono Negative Control (contains goat albumin in tris buffer with 0.2% sodium azide)
- 1 Work Station
- 1 Directional Insert

Note: Extra components (tubes, pipettes, capillary pipettes) have been provided for your convenience.

Store the Test Sticks and reagents tightly capped at 15-30 C.

Do not use the Test Sticks or reagents after their expiration dates.

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection containers.
- A timer or watch.

PRECAUTIONS

- **Caution:** Federal Law restricts sales of this device to or on the order of a licensed practitioner.
- All reagents are for in vitro diagnostic use only.
- Follow your laboratory safety guidelines in the collection, handling, storage and disposal of patient specimens and all items exposed to patient specimens.
- The Diluent and Controls contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azide. Large quantities of water must be used to flush discarded Diluent or Controls down a sink.
- Do not interchange or mix components from different kit lots.

HAZARDS AND PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com for Hazards and Precautionary Statements.


SPECIMEN COLLECTION AND PREPARATION

Serum, Plasma, or Whole Blood Sample

Obtain specimens by acceptable medical technique. Collect whole blood samples using a tube containing EDTA or heparin as an anticoagulant. Other anticoagulants have not been tested. Serum and plasma specimens may be refrigerated (2 - 8 C) and tested within 48 hours; serum and plasma specimens held for longer times should be frozen (below -10 C) and tested within 3 months. Test whole blood specimens within 24 hours. Specimens must be at room temperature (15 - 30 C) when tested.

Fingertip Whole Blood


Hold the capillary pipette horizontally, and touch the tip of the pipette to the drop of blood on the patient's finger until it fills completely to the line. **Note:** Filling is automatic; never squeeze the pipette bulb while collecting sample.

TEST PROCEDURE

- Addition of Specimen:

For serum, plasma, or whole blood samples in tubes:
Use the Transfer Pipette provided and add one drop to the Test Tube.

For fingertip blood:
To dispense all the patient sample, place pipette into the test tube and squeeze the bulb. Remove pipette and discard in the appropriate biohazard container. **NOTE:** If pipette does not fully dispense the patient's sample, recollect the sample using a new capillary pipette.

- Slowly add 1 drop of Diluent to the bottom of the Test Tube.


- Mix.

- Remove the Test Stick(s) from the container. Re-cap the container immediately.

- Place the Absorbent End of the Test Stick into the treated sample. Leave the Test Stick in the Test Tube.

- Read results at 5 minutes. Positive results may be read as soon as the red Control Line appears.


INTERPRETATION OF TEST RESULTS

Notes


A blue or red line which appears uneven in color density is considered a valid result.

Positive


A blue Test Line and a red Control Line is a positive result for the detection of infectious mononucleosis heterophile antibody. Note that the blue line can be any shade of color and can be lighter or darker than the line in the picture.

Negative

A red Control Line but no blue Test Line is a negative result. No infectious mononucleosis heterophile antibody has been detected.

Invalid

If after five minutes, no red Control Line appears or background color makes reading the red Control Line impossible, the result is invalid. If this occurs, repeat the test on a new Test Stick or call Meridian Bioscience Technical Services.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

External Quality Control

For external QC testing, use the controls provided in the kit. Add one free falling drop of control to the Test Tube and then proceed in the same manner as with a patient sample. Quality Control requirements should be established in accordance with local, state and federal regulations or accreditation requirements. Minimally, Meridian Bioscience recommends that positive and negative external controls be run with each new lot and with each new untrained operator. Some commercial controls may contain interfering additives. The use of these controls is not recommended.

Internal Quality Controls

The ImmunoCard STAT! Mono Test provides two levels of internal procedural controls with each test procedure.

- The red Control Line is an internal positive procedural control. The Test Stick must absorb the proper amount of test material and be working properly for the red Control Line to appear.
- A clear background is an internal negative procedural control. If the test has been performed correctly and the Test Stick is working properly, the background will clear to give a discernible result.

If the red Control Line does not appear, the test may be invalid. If the background does not clear and interferes with the test result, the test may be invalid. If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

EXPECTED VALUES

A heterophile antibody response is observed in approximately 80-90% of adults and children with EBV-caused IM. This percentage drops to approximately 50% for children under four years of age.¹ While the incidence of IM reflects wide seasonal, ethnic and geographical variation, a large epidemiological study noted that the highest incidence of symptomatic IM occurred during late adolescence (15-24 years of age).²

LIMITATIONS OF THE PROCEDURE

- As with all diagnostic assays, the results obtained by this test yield data that must be used as an adjunct to other information available to the physician.
- The ImmunoCard STAT! Mono Test is a qualitative test for the detection of IM heterophile antibody.
- A negative result may be obtained from patients at the onset of the disease due to heterophile antibody levels below the sensitivity of this test kit. If symptoms persist or intensify, the test should be repeated.
- Some segments of the population with acute IM are heterophile antibody negative.¹

SPECIFIC PERFORMANCE CHARACTERISTICS

A total of 439 specimens (183 serum, 176 plasma and 80 whole blood) were evaluated by two clinical labs in a clinical study. Test results of the ImmunoCard STAT! Mono Test were compared to results obtained with a commercially available latex particle agglutination test for the qualitative determination of infectious mononucleosis heterophile antibodies. Discrepancies between the results given by the ImmunoCard STAT! Mono Test and the latex particle agglutination test were resolved by Epstein-Barr Virus (EBV) specific serological assays. In these assays, the specific antibodies to the EBV capsid antigen (IgM) and EBV nuclear antigen-1 (IgM and IgG) were determined.

Serum Specimens:

		Comparative Test	
		+	-
ImmunoCard STAT! Mono Test	+	74	8*
	-	0	101

*6 out of 8 tested positive by EBV testing

Plasma Specimens:

		Comparative Test	
		+	-
ImmunoCard STAT! Mono Test	+	67	15*
	-	0	94

*8 out of 15 tested positive by EBV testing

Whole Blood Specimens:

		Comparative Test	
		+	-
ImmunoCard STAT! Mono Test	+	30	3*
	-	0	47

*1 out of 3 tested positive by EBV testing

All Specimens:

		Comparative Test	
		+	-
ImmunoCard STAT! Mono Test	+	171	26*
	-	0	242

*15 out of 26 tested positive by EBV testing

When compared to a commercially available latex particle agglutination test for infectious mononucleosis heterophile antibodies, the ImmunoCard STAT! Mono Test showed a sensitivity of 100% and a specificity of 90.3%. The overall agreement was 94.1%.

Fifteen of the 26 discrepant samples were determined to be recent or acute EBV infections by EBV serological testing, in which case the sample was considered positive. Including the samples confirmed positive by EBV serological testing, the overall clinical study specificity of the ImmunoCard STAT! Mono Test is 95.9% and the overall sensitivity is 100%.

POL Studies

An evaluation of the ImmunoCard STAT! Mono Test was conducted at three physicians' offices or clinical laboratories where testing was performed by personnel with diverse educational backgrounds. Each site tested the randomly coded panel consisting of negative (5), low positive (3) and moderate positive (4) specimens for three days. The results obtained had 99.1% agreement (107/108) with the expected results.

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ImmunoCard STAT! is a registered trademark of Meridian Bioscience, Inc.

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Licensed under U.S. Patent Nos. 5,714,389; 5,989,921; 6,485,982 and related non-U.S. patents and patent applications.

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Rev. 3854-6, 11/18

Manufactured For

Meridian Bioscience, Inc.
Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244 USA
Telephone: 513.271.3700
Orders/Customer Service:
800.543.1980
Technical Support Center:
800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124

Authorized Representative

Meridian Bioscience Europe S. r. L
Via dell'Industria, 7
20020 Villa Cortese, Milano
ITALY
Tel: +39 0331 43 36 36
Fax: +39 0331 43 36 16
Email: info@meridianbioscience.eu
WEB: www.meridianbioscience.com/eu

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bnl@meridianbioscience.eu

Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bnl@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/package of this product.

Key Guide to symbols (Guida ai simboli, Guide des symbols, Guía de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfils the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenzubereitung, in dem sich Probenverdünnungspuffer befindet
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Numero de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Eingrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" test / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de température / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / testgerät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
CONJ	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
Rx Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	DET REAG	Detection Reagent / Reagente Diretto / Réactif de Detection / Reactivo de Detección / Nachweis Reagenz
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / ne pas utiliser si le paquet est endommagé / No use si el paquete está dañado / Nicht verwenden, wenn die Verpackung beschädigt ist	TUBE	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.