# curian®

## CAMPY

### A Rapid Fluorescent Immunoassay for the Detection of a *Campylobacter*-specific antigen in Stool Specimens

REF 760730

IVD

**Rx Only**  
For Professional Use Only

#### INTENDED USE

Curian Campy, for use with the Curian Analyzer, is a rapid, qualitative fluorescent immunoassay for the detection of a *Campylobacter*-specific antigen in human fecal specimens. Curian Campy is intended to detect *C. jejuni*, *C. coli*, *C. upsaliensis*, and *C. lari* in human stool from patients with signs and symptoms of gastroenteritis. The test is intended for use with unpreserved fecal specimens or preserved fecal specimens in transport media. Test results are to be used in conjunction with information available from the patient clinical evaluation and other diagnostic procedures. Curian Campy is intended to aid in the diagnosis of *Campylobacter* infection.

#### SUMMARY AND EXPLANATION OF TEST

*Campylobacter* is a gram-negative, microaerophilic bacterium that is carried in the intestine of many wild and domestic animals, especially avian species including poultry.<sup>1</sup> In fact, poultry is considered a major source for *Campylobacter* infections in humans, with transmission occurring through ingestion of contaminated food and water.<sup>2,3</sup> *Campylobacter* infection is referred to as campylobacteriosis, which presents with symptoms of acute watery or bloody diarrhea, abdominal pain, and fever that usually begin within 24 to 72 hours following ingestion, lasting an average of 6 days.<sup>4</sup> In some instances, infection can progress to life-threatening extragastrointestinal diseases.<sup>5</sup>

*Campylobacter* infection is one of the most widespread infectious diseases of the last century, and evidence suggests that the incidence and prevalence of campylobacteriosis have increased in both developed and developing countries over the last decade, especially in children and the elderly.<sup>1,4</sup> In the United States (from 1998 to 2008), the annual number of campylobacteriosis cases was estimated to be 845,024 cases, resulting in 8,463 hospitalizations and 76 deaths. The U.S. Food-Borne Disease Active Surveillance Network (from 1996 to 2012) reported an annual incidence of 14.3 per 100,000 population.<sup>4</sup>

Over 90% of *Campylobacter* infections are caused by *C. jejuni*, about 5% by *C. coli* and the remainder by other *Campylobacter* species, such as *C. lari* or *C. upsaliensis*.<sup>6</sup> There is increasing acknowledgment of the clinical importance of *Campylobacter* species other than *C. jejuni* and *C. coli*, due to the recognition of these species as emerging human and animal pathogens.<sup>4</sup> The conventional approach to the laboratory diagnosis of *Campylobacter* infection is through culture methods, requiring the use of specialized media and a selective incubation temperature of 42 C for several days in an artificially created microaerophilic environment in order to recover organisms from a stool specimen. Recent studies suggest that the culture methodology has its limitations, and that immunoassays are convenient and dependable alternatives for the laboratory diagnosis of *Campylobacter* infection.<sup>7,8</sup> The Curian Campy assay provides rapid, reliable results for the detection of *Campylobacter* specific antigens using a fluorescent analyzer to remove subjectivity and assist in the diagnosis of *Campylobacter* infection.

#### BIOLOGICAL PRINCIPLES

The Curian Campy assay consists of a test strip enclosed in a plastic frame (test card), positive control reagent, and assay-specific Aioprep™ sample preparation device. Curian Campy is a lateral flow-based immunoassay for the direct detection of a *Campylobacter*-specific antigen in human stool. Curian Campy utilizes monoclonal antibodies, specific for an antigen common to *C. jejuni*, *C. coli*, *C. lari* and *C. upsaliensis*, as the capture and detector antibodies.

The Aioprep is pre-filled with blue tinted Sample Diluent/Negative Control and contains a filter and a sample metering device. A sample of the patient's stool specimen is transferred from the collection container to the Aioprep, using either the included transfer pipette (for preserved or unpreserved liquid/semi-solid specimens), or the assay-specific sample collection brush (for preserved or unpreserved formed/solid specimens) available separately as an accessory. When using the sample collection brush, the brush with sample is inserted through the metering device to remove excess stool and then pushed directly into the Sample Diluent/Negative Control. When using the transfer pipette, the sample is added directly into the Sample Diluent/Negative Control.

The diluted sample is mixed and dispensed drop-wise into the sample port of the Curian Campy test card. If present, the *Campylobacter* antigen binds to the monoclonal detector antibody conjugated to fluorescent particles, forming a complex. As the sample moves through the test strip, the anti-*Campylobacter* capture antibody, bound to the assay membrane at the test position of the strip, binds the complex and yields a test line. When antigen is not present, a complex is not formed, and a test line will not form. As the sample continues to move further up the test strip, the polyclonal capture antibody, bound to the assay membrane at the control position of the strip, binds the conjugated antibody and yields a control line. A line at the control position of the test strip should be present each time a sample or external control is tested. If the control line is not generated, adequate sample flow has not occurred, and the Curian Analyzer will consider the test invalid.

#### REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Curian Campy Test Card:** A test strip enclosed in a plastic frame which is in a foil pouch with a desiccant. Supplied ready to use.
2. **Curian Campy Aioprep Sample Preparation Device/ Negative Control:** The Aioprep device is fitted with a metering insert and dropper tip. Supplied ready to use. A buffered protein solution containing blue dye and 0.094% sodium azide.
3. **Curian Campy Positive Control:** Inactivated *Campylobacter jejuni* antigen in a phosphate buffered solution containing 0.094% sodium azide. Supplied ready to use.
4. **Transfer Pipettes**

#### MATERIALS PROVIDED SEPARATELY

1. Curian Campy Stool Collection Brushes, Meridian Bioscience, Inc. Catalog 11516 (To order brushes, please contact Meridian's Technical Services Department at 1-800-343-3858 or your local distributor).

#### MATERIALS NOT PROVIDED

1. Disposable gloves, powder free.
2. Wooden applicator sticks
3. Transport Media (Cary Blair or C&S)

#### EQUIPMENT NOT PROVIDED

1. Vortex mixer
2. Interval timer (Optional)
3. Curian Analyzer System, Meridian Bioscience, Inc. Catalog 610190

## PRECAUTIONS

1. All reagents are for *in vitro* diagnostic use only.
2. Rx Only
3. Store the kit at the temperature indicated on labeling when not in use.
4. Do not interchange transfer pipettes, sample collection brushes, or Aiopreps between assays. Specimens must be sampled and prepared with the transfer pipettes or the Curian Campy sample collection brushes (provided as a separate accessory item) and used with the Curian Campy Aioprep device. These components are not interchangeable with other assays.
5. Handle and dispose of all human specimens as if they are biologically hazardous.
6. Inspect foil pouch before removing the test card. Do not use test cards that have holes in the foil pouch or where the pouch has not been completely sealed.
7. Do not use test cards where the desiccant indicator has changed from blue to pink.
8. Inspect Curian Campy Aioprep Sample Preparation Device/ Negative Control prior to use. Tap the Aioprep to ensure the liquid is in the main chamber of the Aioprep.
9. Do not interchange positive control, Aiopreps, or test cards between kit lots.
10. Do not mark over or near the barcode on the test card.
11. The Curian Campy test card must be incubated outside of a laminar flow hood; this can be done either inside the analyzer or on the lab bench. Sample preparation in the Aioprep can be performed inside a laminar flow hood.
12. Do not use transfer pipettes with formed and solid stool specimens.

## HAZARD and PRECAUTIONARY STATEMENTS


 <b>Aioprep Sample Preparation Device / Negative Control</b>	<b>Signal word</b> Danger <b>Hazard statements</b> H360 - May damage fertility or the unborn child <b>Precautionary Statements - EU (§28, 1272/2008)</b> P201 – Obtain special instructions before use P202 – Do not handle until all safety precautions have been read and understood P280 - Wear eye protection/ face protection P308 + P313 – IF exposed or concerned: Get medical device advice/attention P405 – Store locked up P501 – Dispose of contents/container to an approved waste disposal plant.
--	--

## SHELF LIFE AND STORAGE

The expiration date is indicated on the kit box label. Store the kit at 2 - 8 C as indicated on the labeling. Return all kit components to the indicated storage temperature after use.

## PROCEDURAL NOTES

The Curian Campy transfer pipette is diagrammed below:

## OPERATION OF THE CURIAN ANALYZER SYSTEM

The Curian Analyzer is an easy to use, menu driven analyzer system. Instructions to complete testing are provided on the analyzer touchscreen and in the Curian Operator's Manual.

1. Power ON the Curian Analyzer by pressing the power button located on the left side of the analyzer. The Curian Analyzer will initialize and perform SELF TEST.
2. Instrument Checks, Quality Control testing, and Patient Specimen testing are completed by navigating from the Home screen to the test menus. User ID and Sample ID entries are required for Quality Control and specimen testing.
3. Complete use and operation of the Curian Analyzer System by following the instructions in the Curian Operator's Manual and on-screen prompts. **Refer to operator's manual for detailed instructions for the Curian Analyzer.**

## SPECIMEN COLLECTION, HANDLING, AND STORAGE

This procedure is designed to be used with unpreserved stool or stool preserved in Cary-Blair or C&S transport media. DO NOT USE stool on swabs or mixed with fixatives. The patient specimen should be received in an airtight transport container and stored at 2 - 8 C until tested. The specimen should be tested as soon as possible, but preserved specimens may be held for up to 96 hours (four days) at room temperature (19 - 27 C), or unpreserved and preserved specimens (in Meridian Para-Pak C&S media) may be held for up to 7 days at 2 - 8 C prior to testing. If testing cannot be performed within this time frame, unpreserved specimens or specimens preserved in Meridian Para-Pak C&S media should be frozen immediately upon receipt and stored frozen for up to 59 days (or eight weeks) at ≤ -20 C until tested. Specimens may be frozen and thawed up to 5 times.

Specimen Condition	Storage Temperature	Maximum Storage Time
Unpreserved specimens	2 - 8 C (refrigerated)	7 days
	≤ -20 C (frozen)	8 weeks
Preserved specimens in Cary-Blair or C&S media	19 - 27 C (room temperature)	96 hours
	2 - 8 C (refrigerated) <sup>1</sup>	7 days
	≤ -20 C (frozen) <sup>2</sup>	8 weeks
Prepared specimens in Aioprep / Sample Diluent	19 - 27 C (room temperature)	8 hours
	2 - 8 C (refrigerated)	24 hours

<sup>1</sup> Refrigerated storage validated with Meridian Para-Pak C&S media only. If other C&S or Cary-Blair transport media are used, specimens should be stored according to the recommendations in the respective transport media package insert.

<sup>2</sup> Frozen storage validated with Meridian Para-Pak C&S media only. If other C&S or Cary-Blair transport media are used, specimens should be stored according to the recommendations in the respective transport media package insert.

## SAMPLE PREPARATION AND TEST PROCEDURE

**Note: Handle all waste and specimens as biohazardous.**

**Note: Before testing, mix stool specimen thoroughly by vortexing or inverting for liquid/semi-solid specimens, or with a wooden applicator stick for formed/solid specimens.**


Bring all test components, reagents and specimens to 19 - 27 C before testing. Use one Curian Campy test card for each sample. When ready to perform testing, inspect the pouch for damage. If undamaged, remove test card from its foil pouch. Discard the pouch and desiccant. **Do not use if desiccant is pink. Do not cover or mark on the barcode.**

1. **Sample Preparation:**

Prior to sampling, **thoroughly mix stool specimen** by vortexing for approximately 5 seconds or inverting (e.g., 3 - 4 times), or with a wooden applicator stick for formed/solid specimens. Remove the yellow top cap from the Aioprep device and set aside.

**Note: Prepared specimens in the Aioprep / Sample Diluent may be held at room temperature (19 - 27 C) for up to 8 hours or may be stored at refrigerated temperature (2 - 8 C) for up to 24 hours, prior to proceeding to step 2 for testing.**

**Note: Ensure that the measurement markings on the transfer pipette are not obscured during sample aspiration.**

**Pipettable, liquid/semi-solid specimens:**

a. **Human stool specimens, preserved in Cary-Blair or C&S media:**

i. Using the provided transfer pipette, transfer 300 µL of preserved stool (third mark from pipette tip) directly into the Sample Diluent by inserting the tip of the pipette through the metering insert while expelling specimen to ensure the stool is transferred. **NOTE: Transferring too little specimen, or failure to mix and thoroughly suspend the specimen in the Aioprep/ Sample Diluent, may result in a false-negative result. Refer to the examples below.** Recap the Aioprep device.


ii. Proceed to step 2.

b. **Human stool specimens, unpreserved:**

i. Using the provided transfer pipette, transfer 100 µL of unpreserved stool (first mark from pipette tip) directly into the Sample Diluent by inserting the tip of the pipette through the metering insert while expelling specimen to ensure the stool is transferred. **NOTE: Transferring too little specimen, or failure to mix and thoroughly suspend the specimen in the Aioprep/ Sample Diluent, may result in a false-negative result. Refer to the examples below.** Recap the Aioprep device.

ii. Proceed to step 2.

User may encounter slight resistance when inserting tip of pipette through metering insert.

**Non-pipettable, formed/solid specimens (either preserved in Cary-Blair or C&S media or unpreserved):**

a. Using the Curian Campy sample collection brush, insert the brush into the stool specimen to coat the bristles of the brush. It is important that only the brush is covered with stool. Do not get stool on the shaft or handle of the brush; wipe any excess stool on the inside of the specimen container. Refer to the examples below.

b. Insert the Curian Campy sample collection brush into the Aioprep device, pushing the brush through the metering insert. Snap off the sample collection brush at the marked 'Break' point and discard the handle. Recap the Aioprep device, leaving the Curian Campy sample collection brush in the device.

c. Proceed to step 2.

2. **Do not invert or shake Aioprep.** Thoroughly mix the sample in the Aioprep sample diluent by vortexing for approximately 5 seconds. **Sample must be thoroughly suspended in the Aioprep sample diluent.** Remove the white tip cap from the bottom of the Aioprep and discard. While holding the Aioprep vertically, squeeze in the center of the barrel to add 3 drops into the 'SAMPLE' port of the test card. Discard the Aioprep device immediately.

3. Start the Test Read by navigating to the 'TEST' Menu, select 'ANALYZE NOW' or 'INCUBATE AND ANALYZE', and then follow the on-screen instructions.
 

**For Analyze Now:**

  - a. Incubate the test card at 19 - 27 C on the benchtop for 20 minutes.
  - b. Enter Sample ID and press Ok.
  - c. Within 2 minutes of the end of incubation, insert the test card into the drawer of the Curian Analyzer and close the drawer. Analysis of the reaction will be initiated by the Curian Analyzer.
  - d. The Curian Analyzer will analyze the test card and automatically report the test result.

**For Incubate and Analyze:**

  - a. Enter Sample ID and press Ok.
  - b. Immediately insert the test card into the drawer of the Curian Analyzer and close the drawer. The analyzer will time the incubation for 20 minutes.
  - c. After the end of the incubation, the Curian Analyzer will analyze the test card and automatically report the test result.

### ASSAY EXTERNAL QUALITY CONTROL (QC) TEST PROCEDURE

Note: Good Laboratory Practice (GLP) guidelines<sup>9</sup> recommend the use of control material. Users should follow the appropriate federal, state, and local guidelines concerning the running of external quality controls.

Bring all test components and reagents to 19 - 27 C before testing. Use one Curian Campy test card for each external QC test (Negative and Positive Controls). When ready to perform testing, remove the test card from its foil pouch. Discard the pouch and desiccant. **Do not use if desiccant is pink. Do not mark over the barcode on the test card.**

Start the Test Read by navigating to the 'TEST' Menu, selecting the 'QC TEST', and then selecting 'ANALYZE NOW' or 'INCUBATE AND ANALYZE'. Follow the on-screen instructions.

#### Negative Control

1. Remove the white tip cap from the bottom of the Aioprep and discard. While holding the Aioprep vertically, squeeze in the center of the barrel to add 3 drops of Negative Control reagent (Sample Diluent) into the 'SAMPLE' port of the test card. Discard the Aioprep device immediately.
2. **For Analyze Now:**
  - a. Incubate the test card at 19 - 27 C on the benchtop for 20 minutes.
  - b. Enter Sample ID and press Ok.
  - c. Within 2 minutes of the end of incubation, insert the test card into the drawer of the Curian Analyzer and close the drawer. Analysis of the reaction will be initiated by the Curian Analyzer.
  - d. Select 'Curian Campy QC Negative'.
  - e. The Curian Analyzer will analyze the test card and automatically report the test result.

**For Incubate and Analyze:**

  - a. Enter Sample ID and press Ok.
  - b. Insert test card immediately into the Curian Analyzer.
  - c. Select 'Curian Campy QC Negative'.
  - d. The Curian Analyzer will time the incubation for 20 minutes, then read the test card and automatically report the test result.

#### Positive Control

1. Invert the positive control bottle to mix. Remove the Positive Control tip cap. Holding vertically, squeeze 3 drops of Positive Control into the 'SAMPLE' port of the test card.
2. **For Analyze Now:**
  - a. Incubate the test card at 19 - 27 C on the benchtop for 20 minutes.
  - b. Enter Sample ID and press Ok.
  - c. Within 2 minutes of the end of incubation, insert the test card into the drawer of the Curian Analyzer and close the drawer. Analysis of the reaction will be initiated by the Curian Analyzer.
  - d. Select 'Curian Campy QC Positive'.
  - e. The Curian Analyzer will analyze the test card and automatically report the test result.

**For Incubate and Analyze:**

  - a. Enter Sample ID and press Ok.
  - b. Insert the test card immediately into the drawer of the Curian Analyzer and close the drawer.
  - c. Select 'Curian Campy QC Positive'.
  - d. The Curian Analyzer will time the incubation for 20 minutes, then read the test card and automatically report the test result.

**If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor. Do NOT perform patient testing if either of the QC test results are incorrect.**

### INSTRUMENT CHECK (IC) TEST PROCEDURE

Note: The analyzer's External Control (IC Test) has a default schedule of 30 days. The IC Test will be required to be performed prior to running a QC test or patient test if this has not been performed within this timeframe. The user will be locked-out of running a QC test or patient test if this test is overdue.

1. Power ON the Curian Analyzer by pressing the power button located on the side of the Curian Analyzer. The analyzer will initialize and perform SELF TEST (automatic) checks.
2. From the HOME screen, select 'TEST', and then select 'INSTRUMENT CHECK'.
3. Insert the Fluorescent IC Card into the drawer. A 'PASS' or 'FAIL' result will be displayed once the analysis is complete.

**If the expected result is not observed, repeat the IC Test as the first step in determining the root cause of the failure. Refer to the Curian Analyzer Operator's Manual for detailed instructions for the analyzer.**

### INTERPRETATION OF RESULTS

Result interpretation is completed automatically by the Curian Analyzer system. The result will be shown on screen. Results can be retrieved from analyzer storage, printed and/or exported.

#### Patient Test Results:

One of the following result interpretations will be generated by the Curian Analyzer for patient specimens.

1. **Positive:** *Campylobacter* antigen is present in sample.
2. **Negative:** *Campylobacter* antigen is not present in sample.
3. **Control Invalid:** Result indicates inadequate flow of sample. Repeat the test from the original stool specimen. If repeat testing yields the same results, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor for additional assistance.

#### External Quality Controls:

One of the following result interpretations will be generated by the Curian Analyzer for External Controls.

1. **Pass:** This indicates the test card and reagents are performing as intended.
2. **Fail:** This result indicates the test card and/or reagents are not performing as intended or a user error has occurred. Test should be repeated to assist in trouble-shooting the error. If repeat testing results in a failed output, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor for additional assistance. **Do NOT perform patient testing if either of the QC test results are incorrect.**
3. **Control Invalid:** This result indicates inadequate sample flow. Repeat the test. If repeat testing yields the same results, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor for additional assistance.

#### Instrument Check:

One of the following result interpretations will be generated by the Curian Analyzer for Instrument Check.

1. **Pass:** This indicates that the essential and critical components of the Curian Analyzer are working correctly.
2. **Fail:** This indicates that there may be an issue with the Curian Analyzer. Refer to the Curian Operator's Manual if an error message has occurred. Repeat testing to assist in trouble-shooting the error. If repeat testing results in a failed output, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor for additional assistance.

## QUALITY CONTROL

There are four controls for the Curian Campy test system. There is an external and internal control for both the analyzer and assay.

1. Analyzer Internal Control: Self Test (automatic)
2. Analyzer External Control: Instrument Check (IC) Test
3. Assay External Controls: Positive and Negative Controls (QC test)
4. Assay Internal Control: Control Line

## EXPECTED VALUES

The prevalence of *Campylobacter* infection observed during the 2020 study period for prospectively collected specimens was 1.2% (18/1474). The prevalence of *Campylobacter* spp. by sex assigned at birth for the tested population is provided below. The age demographic for the tested population ranged from less than 1 year old to 99 years old, with an average age of 48 years old.

Sex Assigned at Birth	Samples Positive by Curian Campy and Reference Method	Prevalence in Tested Population
Males (n=578)	10	1.7%
Females (n=896)	8	0.9%
Overall (n=1474)	18	1.2%

## LIMITATIONS OF THE PROCEDURE

1. The Curian Campy assay must only be used with the Curian Analyzer.
2. Curian Campy is a qualitative, *in vitro* diagnostic test. The Curian Analyzer will only provide qualitative results. This test is not intended to provide quantitative results.
3. Curian Campy detects both viable and nonviable *Campylobacter* bacteria and may yield a positive result in the absence of living organisms.
4. A negative test result does not definitively rule-out the presence of *Campylobacter* species in suspected patients. Levels of organism may be present in feces beneath the limit of detection for the Curian Campy assay, and therefore, if *Campylobacter* is suspected, alternative testing should be conducted.
5. A negative result may occur if the sample was collected, transported, or stored improperly.
6. Failure to follow the test procedure may adversely affect test performance and/or invalidate the test result.
7. Test results must be evaluated in conjunction with other clinical data available to the physician.
8. Negative test results do not rule out other possible infections.
9. Positive test results do not rule out co-infections with other pathogens.
10. Curian Campy is designed to be used with unpreserved stool samples or stools preserved in Cary-Blair or C&S media. DO NOT USE stool on swabs or in other transport media (e.g., formalin, polyvinyl alcohol). No data exists on the effects of colonic washes, barium enemas, laxatives, or bowel preparations on the performance of Curian Campy. All of these procedures can result in extensive dilution or the presence of additives that may affect test performance.
11. Transferring too little specimen, or failure to mix and completely suspend the specimen in the Sample Diluent, may result in a false negative test result.
12. The performance of this test has not been evaluated for use in patients without signs and symptoms of gastroenteritis.
13. *Campylobacter helveticus* (strain ATCC 51209) at levels greater than  $3.75 \times 10^6$  CFU/mL in unpreserved stool or greater than  $7.50 \times 10^6$  CFU/mL in preserved stool may cross react or interfere with the performance of the test.
14. Non-pipettable, formed/solid specimens should only be processed with the Curian Campy Sample Collection Brush (Meridian Bioscience, Inc. Catalog 11516).

## SPECIFIC PERFORMANCE CHARACTERISTICS

### CLINICAL PERFORMANCE

#### Prospective Study

The Curian Campy assay was evaluated from July 2020 to December 2020 at five clinical study sites representing geographically distinct regions throughout the United States. There were 1,474 specimens from patients with signs and symptoms of gastroenteritis for whom a diagnostic *Campylobacter* test had been ordered by a practicing physician, prospectively collected and enrolled into the study. All specimens were tested at the study sites with the Curian Campy assay and either had current standard of care *Campylobacter* culture and speciation performed and results available (reference method) or culture and speciation was performed as part of the study. The vast majority of specimens were pipettable and were processed with the Curian Campy pipette included in the kit. A small subset of specimens were non-pipettable and were processed with the Curian Campy Collection Brush (sold separately by Meridian). Clinical performance (sensitivity and specificity) for prospective specimens against the reference method (culture and speciation) are presented in the following table. There were no observable differences in performance of the Curian Campy assay with respect to study site, storage condition, kit lot, or patient gender or age. Prospective specimens with discordant results between the Curian Campy assay and the reference method were further evaluated using standard of care (SoC) results obtained with an FDA-cleared commercial nucleic acid amplification test (NAAT); results of discordant testing are footnoted below.

#### Curian Campy Overall Performance for Prospective Specimens versus Culture and Speciation

		Reference Method: Culture and Speciation			Parameter	Estimate	95% CI
		Positive	Negative	Total			
Curian Campy Assay	Positive	18	28**	46	Sensitivity	85.7%	[65.4% - 95.0%]
	Negative	3*	1425	1428	Specificity	98.1%	[97.2% - 98.7%]
	Total	21	1453	1474			

\* The Standard of Care (SoC) testing of two of the three false negative specimens by a high sensitivity, FDA-cleared nucleic acid amplification test (NAAT) assay showed that one of the specimens produced a negative *Campylobacter* result, whereas one of the specimens produced a positive *Campylobacter* result; the third specimen was not subjected to NAAT testing as part of the SoC.

\*\* Of the 28 false positive specimens, 10 were subjected to SoC testing by a high sensitivity, FDA-cleared NAAT assay. Two of these 10 specimens produced a positive *Campylobacter* result, whereas eight produced a negative *Campylobacter* result. Eighteen specimens were not subjected to NAAT testing as part of the SoC.

#### Archived Study

To further estimate sensitivity and specificity of the Curian Campy assay, 290 archived specimens with culture and speciation results were retrospectively tested for *Campylobacter* antigen using the Curian Campy assay at all five study sites. The clinical performance (sensitivity and specificity) for archived specimens against the reference method (culture and speciation) are presented in the table below.

#### Curian Campy Overall Performance for Archived Specimens versus Culture and Speciation

		Reference Method: Culture and Speciation			Parameter	Estimate	95% CI
		Positive	Negative	Total			
Curian Campy Assay	Positive	28	5	33	Sensitivity	96.6%	[82.8% - 99.4%]
	Negative	1	256	257	Specificity	98.1%	[95.6% - 99.2%]
	Total	29	261	290			

### Contrived Study

Additional testing at each site was conducted with contrived samples at 2xLoD and 8xLoD for *C. coli*, *C. upsaliensis*, and *C. lari* and 3 reagent lots. All 210 samples tested as expected with n=150 positive *Campylobacter* spp. samples and n=60 negative samples yielding 100% correlation with the anticipated results. Both sample matrices were represented with n=105 unpreserved stool samples and n=105 stool samples preserved in C&S media. The overall performance of the Curian Campy assay compared to the anticipated results is presented in the table below.

### Curian Campy Overall Performance for Contrived Samples versus Contrived Anticipated Results

		Contrived: Anticipated Result			Parameter	Estimate	95% CI
		Positive	Negative	Total			
Curian Campy Assay	Positive	150	0	150	PPA	100.0%	[97.5% - 100.0%]
	Negative	0	60	60	NPA	100.0%	[94.0% - 100.0%]
	Total	150	60	210			

### ANALYTICAL SENSITIVITY

Analytical sensitivity studies were performed to determine the analytical limit of detection (LoD) of quantified *Campylobacter* whole cell stocks (*C. jejuni*, *C. coli*, *C. upsaliensis*, and *C. lari*) in human stool matrix for the Curian Campy assay. The LoD is defined as the lowest concentration of the target analyte that produces positive results ≥ 95% of the time.

The LoD values determined for the Curian Campy assay for each intended species in unpreserved and preserved (C&S, Cary Blair) stool matrix are listed below.

<i>C. jejuni</i>		<i>C. coli</i>		<i>C. upsaliensis</i>		<i>C. lari</i>	
CFU/mL	CFU/test	CFU/mL	CFU/test	CFU/mL	CFU/test	CFU/mL	CFU/test
<b>Unpreserved Stool Matrix</b>							
4.00x10 <sup>5</sup>	1818	3.00x10 <sup>6</sup>	13636	1.62x10 <sup>6</sup>	7386	5.00x10 <sup>6</sup>	22727
<b>Preserved (C&amp;S) Stool Matrix</b>							
7.25x10 <sup>5</sup>	2266	1.57x10 <sup>7</sup>	49063	1.18x10 <sup>6</sup>	3681	1.16x10 <sup>7</sup>	36250
<b>Preserved (Cary Blair) Stool Matrix</b>							
7.25x10 <sup>5</sup>	2266	1.57x10 <sup>7</sup>	49063	2.36x10 <sup>6</sup>	7375	1.16x10 <sup>7</sup>	36250

### ASSAY REACTIVITY/ INCLUSIVITY

Several strains of *C. jejuni*, *C. coli*, *C. upsaliensis*, and *C. lari* were used to evaluate the specificity of the Curian Campy assay using both unpreserved and preserved stool (in C&S).

All strains listed below generated positive results when tested.

<i>C. jejuni</i>	<i>C. coli</i>	<i>C. upsaliensis</i>	<i>C. lari</i>
Zeptomatrix Z086	ATCC 43482	ATCC BAA-1059	ATCC BAA-1060
ATCC 33292	ATCC 43476	ATCC 49815	ATCC 35222*
ATCC 49350	ATCC 43478	2017/0506H	ATCC 35223*
ATCC 43442	ATCC 43485	2016/1950*	2013/0823H*
ATCC 33560	ATCC BAA-1061	ATCC 43954	2015/0814*
		ATCC 43953*	2015/2983
		2016/2697*	2016/0235*
		2018/0319H*	2015/0519
		2016/1931	2015/1582
		ATCC 700558*	2015/2189

\*The following *C. upsaliensis* and *C. lari* strains exhibited elevated LoDs in comparison to the reference strains (ATCC 49816 and ATCC 43675) that were used in the LoD determination:

Species	Strain	Unpreserved Stool	Stool Preserved in C&S
<i>C. upsaliensis</i>	2016/1950	4x	8x
	ATCC 43953	2x	8x
	2016/2697	56x	78x
	2018/0319H	40x	55x
	ATCC 700558	24x	33x
<i>C. lari</i>	2013/0823H	8x	4x
	2015/0814	8x	4x
	2016/0235	8x	10x
	ATCC 35222	8x	10x
	ATCC 35223	24x	10x

### REPRODUCIBILITY

Ten manufactured reference panels per sample type (unpreserved stool and stool preserved in C&S transport media) were supplied to three laboratories for this reproducibility study. The panels each contained sixteen blinded samples, including five high-negative samples (just below C<sub>5</sub>), five low-positive samples (~C<sub>95</sub>), five moderate-positive samples (3x higher than the C<sub>95</sub>), and one true negative sample. Panel members were manufactured by spiking *Campylobacter jejuni* organism into pooled negative stool matrix. Positive and negative external controls were run daily.

Different operators, per site, tested each panel type per day for five days, totaling 320 samples (160 unpreserved stools and 160 stools preserved in C&S transport media). Each operator ran the external controls each day of testing. Three lots of Curian Campy kits and four Curian Analyzers were used in this study. The results of the study are provided in the table below.

Sample Type	Sample Category	Site 1 - Percent Agreement		Site 2 - Percent Agreement		Site 3 - Percent Agreement		Total Percent Agreement	
		Rate	%	Rate	%	Rate	%	Rate	%
Controls	Negative Control	10/10	100.0	11/11	100.0	20/20	100.0	41/41	100.0
	Positive Control	10/10	100.0	11/11	100.0	20/20	100.0	41/41	100.0
Preserved in C&S media	True Negative	10/10	100.0	10/10	100.0	10/10	100.0	30/30	100.0
	High Negative	45/50	90.0	48/49 <sup>1</sup>	98.0	49/50	98.0	142/149	95.3
	Low Positive	50/50	100.0	50/50	100.0	50/50	100.0	150/150	100.0
	Moderate Positive	50/50	100.0	50/50	100.0	50/50	100.0	150/150	100.0
Unpreserved Stool	True Negative	10/10	100.0	10/10	100.0	10/10	100.0	30/30	100.0
	High Negative	48/50	96.0	50/50	100.0	50/50	100.0	148/150	98.7
	Low Positive	50/50	100.0	41/50	82.0	33/50	66.0	124/150	82.7
	Moderate Positive	50/50	100.0	50/50	100.0	50/50	100.0	150/150	100.0

<sup>1</sup> A single replicate was not performed because of instrument error and there was insufficient volume for a repeat test.

<sup>2</sup> For the unpreserved low positive panel member, Percent Agreement (PA) was lower than expected; under-sampling with the Curian Campy transfer pipette was the cause.

## CROSS-REACTIVITY/ MICROBIAL INTERFERENCE

The Curian Campy assay was evaluated for cross-reactivity and microbial interference with the organisms listed below. Unless otherwise indicated, each organism was tested at minimum concentrations of  $1.0 \times 10^7$  CFU/mL for bacteria/fungi or  $1.0 \times 10^5$  TCID<sub>50</sub>/mL for viruses. The assay's reactivity with Norovirus was evaluated using clinical Norovirus-positive stool specimens. None of the organisms showed cross-reactivity or microbial interference in the Curian Campy assay, except for *Campylobacter helveticus*, which was found to be positive at concentrations greater than  $3.75 \times 10^6$  CFU/mL in unpreserved stool and  $7.50 \times 10^6$  CFU/mL in preserved stool.

<i>Acinetobacter baumannii</i>	<i>Klebsiella pneumoniae</i>
<i>Aeromonas hydrophila</i>	<i>Lactobacillus acidophilus</i>
<i>Bacillus cereus</i>	<i>Lactococcus lactis</i>
<i>Bacillus subtilis</i>	<i>Listeria monocytogenes</i>
<i>Bacteroides fragilis</i>	<i>Peptostreptococcus anaerobius</i>
<i>Campylobacter concisus</i>	<i>Plesiomonas shigelloides</i>
<i>Campylobacter fetus</i>	<i>Porphyromonas asaccharolytica</i>
<i>Campylobacter helveticus</i>	<i>Prevotella melaninogenica</i>
<i>Campylobacter hyointestinalis</i>	<i>Proteus vulgaris</i>
<i>Candida albicans</i>	<i>Pseudomonas aeruginosa</i>
<i>Citrobacter freundii</i>	<i>Pseudomonas fluorescens</i>
<i>Clostridium bifermens</i>	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Hilversum
<i>Clostridium difficile</i>	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium
<i>Clostridium perfringens</i>	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Minnesota
<i>Edwardsiella tarda</i>	<i>Serratia marcescens</i>
<i>Enterobacter cloacae</i>	<i>Shigella boydii</i>
<i>Enterococcus faecalis</i>	<i>Shigella dysenteriae</i>
<i>Escherichia coli</i>	<i>Shigella flexneri</i>
<i>Escherichia coli</i> EIEC	<i>Shigella sonnei</i>
<i>Escherichia coli</i> EPEC	<i>Staphylococcus aureus</i>
<i>Escherichia coli</i> ETEC	<i>Staphylococcus aureus</i> (Cowan's)
<i>Escherichia coli</i> O157:H7 (non-toxigenic)	<i>Staphylococcus epidermidis</i>
<i>Escherichia coli</i> O157:H7 (toxigenic)	<i>Streptococcus agalactiae</i>
<i>Escherichia fergusonii</i>	<i>Streptococcus dysgalactiae</i> subsp. <i>equisimilis</i>
<i>Escherichia hermanii</i>	<i>Vibrio parahaemolyticus</i>
<i>Helicobacter pylori</i>	<i>Yersinia enterocolitica</i>
Adenovirus Type 1, 2, 3, 5, 40, 41	Human Coronavirus
Coxsackievirus B2, B3, B4, B5	Human Rotavirus
Echovirus 9, 11, 18	Norovirus
Enterovirus 68, 69, 70, 71	Parechovirus 1 (formerly Echovirus 22)

## TESTS FOR INTERFERING SUBSTANCES

The chemical and biological substances listed below were evaluated at the indicated concentrations for interference in the Curian Campy assay. None of the substances showed interference with the Curian Campy assay performance.

Barium Sulfate (5% w/v)	Mylanta® (4.2 mg/mL)
Benzalkonium chloride (1% w/v)	Naproxen sodium (5% w/v)
Ciprofloxacin (0.25% w/v)	Nonoxynol-9 (1% w/v)
Ethanol (1% w/v)	Nystatin (1% w/v)
Hog gastric mucin (3.5% w/v)	Palmitic Acid/Fecal Fat (40% w/v)
Human blood (40% v/v)	Pepto-Bismol® (5% v/v)
Human hemoglobin (10.0% w/v)	Phenylephrine (1% w/v)
Human urine (5% v/v)	Polyethylene glycol 3350 (10% w/v)
Hydrocortisone (1% w/v)	Prilosec OTC® (5 µg/mL)
Imodium® A-D (5% v/v)	Sennosides (1% w/v)
Kaopectate® (5% v/v)	Simethicone (10% w/v)
Leukocytes (0.05% v/v)	Stearic Acid/Fecal Fat (40% w/v)
Mesalazine (10% w/v)	Tagamet® (5 µg/mL)
Metronidazole (0.25% w/v)	TUMS® (50 µg/mL)
Mineral Oil (10% w/v)	Vancomycin (0.25% w/v)

## BRUSH BRIDGING STUDY

The Curian Campy assay is for use with unpreserved human stool specimens and human stool specimens preserved in Cary-Blair and C&S transport media. Most specimens are easily sampled using the transfer pipette provided in the kit; however, some unpreserved stool specimens are non-pipettable and require use of a specific brush (i.e., Curian Campy Stool Collection Brush, sold separately by Meridian) to adequately collect the sample for analysis in the Curian Campy assay. Fifty-three non-pipettable stool specimens (5 *Campylobacter* positives and 48 negatives) were processed with the brush during the prospective and archived clinical studies, combined.

To further support use of the brush with the Curian Campy assay, an analytical bridging study was performed that evaluated a panel consisting of contrived positive and negative samples collected/processed with the brush. Contrived positive samples were generated by spiking *C. jejuni* (strain ATCC #BAA-1234) into negative, non-pipettable stool samples. Ten samples were prepared at 3x LoD and 25 samples were prepared at 5x LoD. Twenty-five negative samples were also tested. Three non-expert operators tested each sample at one internal site. All positive samples gave the expected positive results, and all negative samples were negative. These results indicate that non-pipettable stool specimens can be collected with the brush prior to testing with the Curian Campy assay.

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**REFERENCES**

1. Silva Joana, Daniela Leite, Mariana Fernandes, Cristina Mena, Paul Anthony Gibbs, and Paula Teixeira. Campylobacter spp. as a foodborne pathogen: a review. *Frontiers in microbiology* 2; 2011:200.
2. Natsos G, NK Mouttotou, S Ahmad, Z Kamran, A Ioannidis, and KC Koutoulis. The genus campylobacter: detection and isolation methods, species identification & typing techniques. *J Hell Vet Med Soc* 2019;70(1):1327-1338.
3. Facciola A, Riso R, Avventuroso E, Visalli G, Delia SA, and Laganà P. Campylobacter: from microbiology to prevention. *J Prev Med Hyg* 2017;58(2):E79.
4. Kaakoush Nadeem O, Natalia Castaño-Rodríguez, Hazel M Mitchell, and Si Ming Man. Global epidemiology of campylobacter infection. *Clinical microbiology reviews* 2015;28(3):687-720.
5. Man Si Ming. The clinical importance of emerging campylobacter species. *Nature reviews gastroenterology & hepatology* 2011;8(12):669.
6. Buettner S, Wieland B, Staerk KDC, and Regula G. Risk attribution of campylobacter infection by age group using exposure modelling. *Epidemiology & infection* 2010;138(12):1748-1761.
7. Bessede, Emilie, Delcamp Adline, Sifré Elodie, Buissonnière Alice, and Mégraud Francis. New methods for detection of campylobacters in stool samples in comparison to culture. *J Clin Microbiol* 2011;49(3):941-944.
8. Granato, Paul A, Chen Li, Holiday Iris, Rawling Russell A, Novak-Weekley Susan M, Quinlan Tammy, and Kimberlee A Musser. Comparison of Premier CAMPY enzyme immunoassay (EIA), ProSpecT Campylobacter EIA, and ImmunoCard STAT! CAMPY tests with culture for laboratory diagnosis of Campylobacter enteric infections. *J Clin Microbiol* 2010;48(11):4022-4027.
9. US Department of Health and Human Services PHS/CDC/NIH Biosafety in Microbiology and Biomedical laboratories, Washington, DC. US Government Printing Office;2007.
10. Epps, Sharon VR, et al. "Foodborne Campylobacter: infections, metabolism, pathogenesis and reservoirs." *International journal of environmental research and public health* 10.12 (2013): 6292-6304.

SN11287

REV. 04/24

 <p>Manufactured By</p>	<p><b>Meridian Bioscience, Inc.</b>            3471 River Hills Drive            Cincinnati, OHIO - 45244 USA  <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><u>Contacts:</u>            Main Telephone (+1) 513.271.3700            Customer Service/Orders 800.543.1980            Technical Support Center 800.343.3858            Information Fax: 513.272.5432            Ordering Fax: 513.271.0124            E-mail: <a href="mailto:info@meridianbioscience.com">info@meridianbioscience.com</a></p>
	<p><b>Meridian Bioscience Europe, SRL</b>            Via Dell'Industria 7, 20035 Villa Cortese            (Milano) ITALY  <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><u>Contacts:</u>            Main Telephone (+39) 0331.433636            E-mail: <a href="mailto:info@meridianbioscience.eu">info@meridianbioscience.eu</a>            Technical Support: <a href="mailto:MBE-TechService@meridianbioscience.eu">MBE-TechService@meridianbioscience.eu</a>            Customer Service/Orders:</p> <ul style="list-style-type: none"> <li>• For Italian Customers:  <a href="mailto:ordini@meridianbioscience.com">ordini@meridianbioscience.com</a></li> <li>• For Distributors / International Customers:  <a href="mailto:Export.CustomerService@meridianbioscience.eu">Export.CustomerService@meridianbioscience.eu</a></li> </ul>
	<p><b>MedEnvoy Switzerland</b>            Gotthardstrasse 28            6302 Zug            Switzerland</p>

**INTERNATIONAL SYMBOL USAGE**

You may see one or more of these symbols on the labeling/packaging of this product:

**Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Zeichenerklärung, Legenda, Guia-chave para símbolos)**

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis / Uiterste gebruiksdatum / Prazo de validade	<b>CONTROL +</b>	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle / Positieve controle / Controllo positivo
<b>LOT</b>	Batch Code / Codice del lotto / Code du lot / Código de lote / chargenbezeichnung / Batchcode / Código do lote	<b>CONTROL -</b>	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle / Negatieve controle / Controllo negativo
<b>IVD</b>	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum / Medisch hulpmiddel voor in-vitrodiagnostiek / Dispositivo médico de diagnóstico in vitro	<b>EC REP</b>	Authorized representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft / Gevolmachtigd vertegenwoordiger in de Europese Gemeenschap / Representante autorizado na Comunidade Europeia
<b>CE</b>	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sul dispositivo medico-diagnostico in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG / CE-markering / Dit product voldoet aan de eisen van Richtlijn 98/79/EG met betrekking tot medische hulpmiddelen voor in-vitrodiagnostiek / Este producto cumple os requisitos da Diretiva 98/79/CE relativa a dispositivos médicos de diagnóstico in vitro	<b>SMP PREP DIL SPE</b>	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, In dem sich Probenverdünnungsbuffer befindet / Monsterpreparatieapparaat met monsteroplossingmiddel / Aparelho de preparação de amostras com diluente de amostras
			CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Risiko/Gefahr / LET OP: Risico van gevaar / ATENÇÃO: potencial perigo
<b>REF</b>	Catalogue number / Numero di catalogo / Référence du catalogue / Numero de catálogo / Bestellnummer / catalogusnummer / Referência		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Eingrieren / Niet invriezen / Não congelar
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisungen beachten / Raadpleeg de gebruiksaanwijzing / Consultar as instruções de Utilização	<b>BUF RXN</b>	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer / Reactiebuffer / Tampão de reação
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller / Fabrikant / Fabricante		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung / Uitsluitend voor prestatie-evaluatie IVD / Apenas para avaliação do desempenho IVD
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" test / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen / Bevat voldoende voor <n> testen / Conteúdo suficiente para <n> testes	<b>SOLN STOP</b>	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung / Stoploestof / Solução de paragem
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung / Temperatuurgrens / Limites de temperatura	<b>CONJ ENZ</b>	Enzyme Conjugate / Conjugato enzimatico / Conjugué enzymatique / Conjugado enzimático / enzymkonjugat / Enzymconjugaat / Conjugado Enzimático
<b>SN</b>	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer / Seriennummer / Número de série	<b>CONTROL</b>	Assay Control / Control del test / Test de contrôle / Control de Ensayo / Kontrolltest / Assaycontrol / Controlo de Ensaio
<b>TEST</b>	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / testgerät / Testinstrument / Dispositivo de teste	<b>REAG</b>	Reagent / Reagente / Réactifs / Reactivos / Reagenzien / Reagens / Reagente
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum / Datum van vervaardiging / Data de fabrico	<b>BUF WASH</b>	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer / Concentrate washbuffer: 20X / Tampão de lavagem
<b>BUF</b>	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer / Buffer / Tampão		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise / Waarschuwing / Aviso
<b>CONJ</b>	Conjugate / Conjugato / Conjugué / Conjugado / Konjugat / Conjugaat / Conjugado	<b>DIL SPE</b>	Specimen Diluent (or Sample Diluent) / Diluente del Campione / diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer / Specimenverdünnungsmittel (of monsterverdünningsmiddel) / Diluente de amostra
<b>SUBS</b>	Substrate / Substrato / Substrat / Substrato / Substrat / Substraat / Substrato	<b>BUF WASH 20X</b>	Wash Buffer Concentration 20X / Soluzione dil lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat / Bufferconcentratie 20x wassen / Concentração do tampão de lavagem: 20X
<b>Rx Only</b>	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso por Receta / verschreibungspflichtig / Gebruik uitsluitend op voorschrift / Apenas para utilização mediante prescrição	<b>DET REAG</b>	Detection Reagent / Reagente Diretto / Réactif de Detection / Reactivo de Detección / Nachweis Reagenz / Detectiereagens / Reagente de detecção
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / ne pas utiliser si le paquet est endommagé / No use if el paquete esta dañado / Nicht verwenden, wenn die Verpackung beschädigt ist / Niet gebruiken als de verpakking beschadigd is / Não utilizar se a embalagem estiver danificada	<b>TUBE</b>	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß / Leeg buisje / Tubo vazio
	Single Use Only / Prodotto Monouso / A usage unique / Para Un Solo Uso / nur für die einmalige Anwendung / Uitsluitend voor eenmalig gebruik / Apenas para uma única utilização	<b>CH REP</b>	Swiss Authorized Representative / Mandatario svizzero / Mandataire Suisse / Representante Autorizado Suizo / Schweizer Bevollmächtigter / Gevolmachtigd vertegenwoordiger in Zwitserland / Representante autorizado na Suíça

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.