**Para-Pak® Systems
for use with 10% Formalin**

(Patent #5624554)

REF 900412

IVD

Rx Only

INTENDED USE

Para-Pak based systems provide standardized procedures for the routine collection, transportation, preservation, and examination of stool specimens for intestinal parasites. Kit systems are designed for easy use by individuals not trained in microbiological procedures and afford an excellent means of minimizing the adverse effects of delay in specimen transportation. Bulk products, for preservation of fresh specimens, are also available.

SUMMARY AND EXPLANATION OF THE TEST

Diagnosis of intestinal parasitic disease is confirmed by recovery and identification of helminth eggs and larvae, or protozoan trophozoites and cysts in the clinical parasitology laboratory. Timely collection and transportation of "fresh" stool specimens to the laboratory cannot always be insured. Workload conditions and priorities in clinical laboratories frequently do not permit immediate examination of "fresh" specimens. Procedures such as incubation, refrigeration or freezing of stool specimens will not guarantee the recovery of all diagnostic stages of all parasites.^{2, 3, 5, 6, 8, 9}

Proper use of the Para-Pak systems thus assures the parasitologist that diagnostic stages of intestinal parasites will be preserved if present in the fecal material.

BIOLOGICAL PRINCIPLES

Para-Pak 10% Formalin preserved specimens may be examined directly and concentrated for recovery of eggs, larvae, and protozoan cysts.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

Each kit consists of one vial containing 15 mL 10% formalin (item #9004) preservative; each with built-in collection spoons. Simple directions for patients and nursing personnel are also provided.

MATERIALS NOT PROVIDED

1. Ethyl acetate (suggested) or diethyl ether (optional)
2. Zinc sulfate solution (specific gravity = 1.18)
3. Physiological saline
4. Cotton tipped applicator sticks
5. Microscope slides and coverslips
6. Centrifuge
7. Microscope
8. Transfer pipettes

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Avoid contact of fixative solution with the skin and eyes. Should contact occur, flush with running water. If irritation should develop, see a physician.
3. Fixative solutions are poisonous. If ingested, dilute by drinking milk or water. Then call local poison center or physician immediately.
4. Due to the infectious nature of unpreserved stools, care and handwashing should be employed when the specimen is collected and manipulated.
5. Formalin is a potential cancer hazard.
6. Any serious incident that has occurred in relation to the device should be reported to Meridian Bioscience, Inc., 3471 River Hills Drive, Cincinnati, Ohio 45244 USA or Technical Support Center 800.343-3858 and competent authority of the EU Member State in which the clinician and/or patient is established.
7. **IMPORTANT:** See SDS for additional safety and hazard information.

HAZARD and PRECUATIONARY STATEMENTS

Para-Pak 10% Formalin

Signal Word

Danger

Hazard Statements

H301 - Toxic if swallowed
H311 - Toxic in contact with skin
H317 - May cause an allergic skin reaction
H330 - Fatal if inhaled
H331 - Toxic if inhaled
H341 - Suspected of causing genetic defects
H350 - May cause cancer
H315 - Causes skin irritation
H318 - Causes serious eye damage
H370 - Causes damage to organs
H402 - Harmful to aquatic life
Contains Formaldehyde, Methyl alcohol

Precautionary Statements - EU (§28, 1272/2008)

P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor.
P321 - Specific treatment (see Section 4 and Section 11 on this label).
P280 - Wear protective gloves/protective clothing/eye protection/face protection.
P321 - Specific treatment (see supplemental instructions on the administration of antidotes on this label).
P361 + P364 - Take off immediately all contaminated clothing and wash it before reuse.
P260 - Do not breathe dust/fume/gas/mist/vapors/spray.
P304 + P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P310 - Immediately call a POISON CENTER or doctor.
P320 - Specific treatment is urgent (see Section 4 and Section 11 on this label).
P403 + P233 - Store in a well-ventilated place. Keep container tightly closed.
P280 - Wear eye protection/ face protection.
P201 - Obtain special instructions before use.
P308 + P313 - IF exposed or concerned: Get medical advice/attention.
P281 - Use personal protective equipment as required.
P202 - Do not handle until all safety precautions have been read and understood.
P270 - Do not eat, drink or smoke when using this product.
P271 - Use only outdoors or in a well-ventilated area.
P284 - Wear respiratory protection.
P272 - Contaminated work clothing should not be allowed out of the workplace.
P307 + P311 - IF exposed: Call a POISON CENTER or doctor/physician.
P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P302 + P352 - IF ON SKIN: Wash with plenty of soap and water.
P312 - Call a POISON CENTER or doctor/physician if you feel unwell.
P361 - Remove/Take off immediately all contaminated clothing.
P332 + P313 - If skin irritation occurs: Get medical advice/attention.
P363 - Wash contaminated clothing before reuse.
P304 + P340 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P405 - Store locked up.
P501 - Dispose of contents/ container to an approved waste disposal plant.

SHELF LIFE AND STORAGE

Expiration dating for Para-Pak System is indicated on the outer packaging. Store at room temperature (15-30 C). Excessive heat and cold should be avoided.

SPECIMEN COLLECTION AND PREPARATION

1. The patient should be cautioned against the use of antacids, barium, bismuth, antidiarrheal medication, or oily laxatives prior to collection of the specimen.
2. To assure recover of parasitic elements which are passed intermittently and in fluctuating numbers, three specimens spaced a few days apart must be examined. In the case of hospitalized patients, it is suggested that all fecal passages be collected for a designated length of time to avoid prolonging the hospital stay.^{2, 6}
3. The specimen is ideally passed into a bedpan but must not be contaminated with urine. Alternatively, a large plastic bag or "saran wrap" may be placed over the toilet seat opening and the specimen passed into the bag. A thoroughly cleaned and dried milk carton, cut so as to remove the upper two thirds of the carton may also be used. It will be easier to collect the specimen if the water supply to the toilet is shut off and water drained from the bowl by "flushing twice".
4. An appropriate (i.e., bloody, slimy, watery) area of stool should be selected and sampled with the collection spoons provided in the caps of the containers. Sufficient stool is added to each container to bring the liquid level up to the "Fill to Here" line. This will result in approximately 5 mL of sample. To insure adequate sampling of the formed stool, material should be removed from the sides, ends, and middle of the bolus.
5. Agitate each specimen with the spoon along the sides of the container, tighten the cap and shake firmly to insure that the specimen is adequately mixed. When mixing is completed, the specimen should appear homogenous.
6. Return the vials to their container, seal the container, and label appropriately.

TEST PROCEDURE

The Para-Pak Systems lend themselves to a wide variety of procedures in common use. The following discussion is not exhaustive and alternatives may be found in the literature cited. While variations may exist from lab to lab, a thorough examination should include at least three steps:

1. Gross examination:
Record the presence of blood, worms, mucus, or proglottids.
2. Direct microscopic examination from the 10% formalin preserved specimen:
 - A. Place a clean glass slide on a sheet of newsprint.
 - B. Add a drop of saline (iodine may be substituted if desired) to the slide.
 - C. Add a representative sample of formalin preserved specimen to the drop of saline and mix thoroughly with the collecting spoon. The newsprint must be just legible through the slide.
 - D. Place a double width coverslip on the suspension and examine immediately.
3. Concentration procedures:
One or more concentration procedures should be employed. No one concentration procedure works equally well for all parasites,^{1, 2, 8} however, two in common usage that lend themselves well to use with the Para-Pak Systems are:
 - A. **Formalin-ether (ethyl acetate)^{4, 10} sedimentation:**
 1. Mix the 10% Formalin specimen thoroughly. The specimen is now ready for processing with the Para-Pak Macro-Con or Para-Pak CON-Trate Stool Concentration Systems. See the appropriate package insert for further directions.
If Para-Pak Macro-Con or Para-Pak CON-Trate[®] are not available, a sufficient quantity of specimen must be strained into a 15 mL conical centrifuge tube through one layer of narrow mesh or two layers of wide mesh gauze to provide the 1.0 mL of sediment. This amount will vary with the size and density of the specimen.
 2. Add approximately 6-8 mL of 10% Formalin (or saline), mix thoroughly and allow to stand five minutes.
 3. Add 3 mL of ethyl acetate or ether then stopper and shake the tube vigorously for at least 30 seconds. Carefully remove the stopper.
 4. Centrifuge at 500 xg for 10 minutes (1800-2200 rpm for most tabletop centrifuges).
 5. Four layers will be apparent:
 - a. Top layer: ethyl acetate or ether
 - b. Second layer: plug of debris
 - c. Third layer: formalin
 - d. Bottom layer: sediment
 6. After ringing the plug of debris from the sides of the tube with an applicator stick, carefully decant the top 3 layers. While keeping the tube inverted, a cotton swab may be used to remove debris sticking to the sides of the tube. This is particularly important for obtaining suitable results with ethyl acetate and avoids solvent bubbles in the wet mount.
 7. Add a few drops of physiological saline or 10% Formalin to resuspend the remaining sediment. If the resulting slides are too dense (newsprint should be legible through them) more saline or formalin may be added.
 8. Iodine and saline mounts are suggested for microscopic examination.
NOTE: Rarely, a specimen may be extremely thick or mucoid, necessitating a preliminary wash. The wash, in physiological saline or tap water, may be implemented between Steps 1 and 2. The washed material should be centrifuged as in Step 4.
 - B. **Zinc Sulfate flotation:**
 1. Thoroughly mix a representative portion of the 10% formalin stool suspension or **fresh** unpreserved specimen in a 15 mL centrifuge tube and q.s. with tap water to approximately 10-12 mL. The amount of specimen to use will vary with its size and density.
 2. Centrifuge one minute at 1000-1200 xg.
 3. If sediment is about one milliliter in volume, decant the supernatant fluid. Otherwise adjust the density of the suspension by adding material from the 10% Formalin suspension or diluting with more water. If adjustment of the sediment was necessary or if the stool is very oily, repeat the wash procedure.
 4. When using formalinized specimens, the specific gravity of the zinc sulfate solution must be adjusted to 1.2.^{1, 2} Fill the tube about half full with zinc sulfate solution and resuspend the sediment by mixing thoroughly with applicator sticks.
 5. Add additional zinc sulfate solution to within a thumb-width of the top.
 6. Centrifuge one minute at 1000-1200 xg.
 7. Carefully remove the tube from the centrifuge and, avoiding agitation, place it upright in a test tube rack or other suitable holder.
 8. Carefully fill the tube to the brim with zinc sulfate solution. Do not allow any overflow.
 9. A clean coverslip may now be placed on top of the tube. If the coverslip does not contact the meniscus of the liquid, carefully add more zinc sulfate solution until it does.
 10. Do not disturb the tube or coverslip for ten minutes.
 11. With a quick, deft motion remove the coverslip straight upward so that a drop of liquid containing eggs and cysts adheres to the center of the coverslip.
 12. The coverslip may now be placed on a clean glass slide. If an iodine mount is desired, place a small drop of iodine on the slide prior to adding the coverslip. Sealing the edge of the slide with Vaspar (Vaseline/paraffin (1:1) mixture) will prevent drying and distortion of the larger eggs.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. Visual inspection: Vials should contain approximately 15 mL of clear fluid and be free from gelling.
2. If gelled, the fixative may be cleared by placing in a 50 C water bath until clear and fluid. Fixative vials containing patient specimens may be cleared once or twice in this manner also, provided the product is used cooled.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

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REFERENCES

1. Bartlett, Marilyn, et al. Comparative evaluation of a modified zinc sulfate floatation technique. J Clin Microbiol. 1978;7:524-528.
2. Brooke MM. Intestinal and urogenital protozoa. Manual of Clinical Microbiology, ASM, Washington, D.C. 2nd ed. 1974;582-601.
3. Burrows RB. Microscopic diagnosis of the parasites of man. Yale University Press, New Haven. 1965.
4. Erdman Dean. Clinical comparison of ethyl acetate and diethyl ether in the formalin-ether sedimentation technique. J Clin Microbiol. 1981;14:483-385.
5. Garcia LS and Ash Lawrence. Diagnostic Parasitology, C. V. Mosby Co., St. Louis, 1979;pp 9-24.
6. Harper K. et al. Advantages of the PVA fixative two-bottle stool collection technique in the detection and identification of intestinal parasites. Pub Health Lab 1957;15:96-108.
7. Melvin DM and Brooke MM. Laboratory procedures for the diagnosis of intestinal parasites, U. S. Department of Health, Education and Welfare Publication (CDC) – 1975;75-82821, Atlanta.
8. Scholten TH. and Yang J. Evaluation of unpreserved and preserved stools for the detection and identification of intestinal parasites. Am J Clin Path. 1974;62:563-567.
9. Swartzwalder J, Clyde Hunter GW, Frye WW. Manual of Tropical Medicine, W. B. Saunders Co., Philadelphia. 1966.
10. Young Kirk H. et al. Ethyl acetate as a substitute for diethyl ether in the formalin-ether sedimentation technique. J Clin Microbiol. 1979;10:852-853.

SN11000

REV. 11/25

 Manufactured By	Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com
EC REP	Meridian Bioscience Europe, SRL Via Dell'Industria 7, 20035 Villa Cortese (Milano) ITALY www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+39) 0331.433636 E-mail: info@meridianbioscience.eu Technical Support: MBE-TechService@meridianbioscience.eu Customer Service/Orders: • For Italian Customers: ordini@meridianbioscience.com • For Distributors / International Customers: Export.CustomerService@meridianbioscience.eu
UK Authorised Representative	Launch Diagnostics Ash House Ash Road Longfield DA3 8JD UK
CH REP	MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland

INTERNATIONAL SYMBOL CHART

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.