For use with the Revogene®

REF 410200

IVD For *in vitro* Diagnostic Use

Rx Only

INTENDED USE

The Revogene® GBS LB assay performed on the Revogene instrument is a qualitative *in vitro* diagnostic test designed to detect Group B *Streptococcus* (GBS) DNA from 18-24 hour LIM broth enrichments of vaginal/rectal specimen swabs obtained from pregnant women. The Revogene GBS LB assay utilizes automated sample processing and real-time polymerase chain reaction (PCR) to detect a *cbf* gene sequence specific to the *Streptococcus agalactiae* genome.

The Revogene GBS LB assay is indicated for the identification of antepartum GBS colonization and does not provide susceptibility results. It is not intended to diagnose or monitor treatment of GBS infection. Culture isolates are needed for performing susceptibility testing as recommended for penicillin-allergic women.

SUMMARY AND EXPLANATION

Group B *Streptococcus* is a Gram-positive bacterium commonly found in the gastrointestinal, genital, and urinary tract of healthy adults. Approximately 10-30% of all pregnant women are colonized with GBS in the vagina or rectum.¹

GBS colonized mothers typically show no symptoms or health effects. However, this type of bacterial infection is associated with serious illnesses in newborns from women who are colonized with the microorganism. This type of disease in newborns frequently causes sepsis (blood infection), pneumonia (lungs infection), and sometimes meningitis (fluid and lining around the brain infection). The fatality rate for infants with early on-set Group B Streptococcal disease (during the first week of life) is currently estimated at 4% to 6%.¹ Surviving infants may experience long-term disabilities including hearing loss, vision loss or mental retardation.²

Transmission of GBS mainly occurs from the mother to her newborn during labor and birth (*intrapartum*).² Screening for GBS colonization in antepartum women between 35 and 37 weeks' gestation, followed by *intrapartum* prophylactic antibiotic treatment for women with positive colonization status, has proven to be an effective mechanism for prevention of perinatal Group B Streptococcal disease. As colonization may be transient, intermittent or persistent throughout pregnancy, screening is most effective when specimens are collected no more than five weeks (35-37 weeks' gestation) prior to delivery and performed after enrichment with selective broth medium.¹

Most GBS testing is performed by culture and can take up to 48 hours for definitive identification of GBS following the initial 18 to 24 hours of incubation of vaginal/rectal swab in a selective broth medium.

The GBS LB assay can provide results from up to eight specimens in approximately one hour once the LIM broth culture sample is available. The assay minimizes operator intervention from the time the single-use disposable microfluidic cartridge (PIE) containing the sample is placed into the Revogene carousel until results are available. The assay minimizes operator intervention from the time the single-use disposable microfluidic cartridge (PIE) containing the sample is placed into the Revogene carousel until results are available.

PRINCIPLES OF THE PROCEDURE

The Revogene automates sample homogenization, sample dilution, cell lysis, DNA amplification, and detection of the amplified PCR products. User intervention is only required for discharging the patient specimen into the Sample Buffer Tube (SBT), transferring the sample into the PIE, and loading/unloading the PIEs into the Revogene carousel.

Each PIE is a completely integrated closed device in which a sample is dispensed and processed through different microfluidic chambers and channels which allow for the sample processing (i.e. homogenization, sample dilution and cell lysis). The liquid from a single sample is transferred by centrifugation from one chamber to the next in sequence and all reagents specific for the PCR reaction are incorporated and dried within the PCR well(s). A Process Control (PrC) is also incorporated into each PIE to control sample processing and amplification steps (Figure 1). The PrC allows for the monitoring of potentially inhibitory substances as well as microfluidic, instrument or reagent failure. No operator intervention is necessary once a PIE is loaded into the Revogene. The amplified products are detected in real time using target-specific TaqMan® chemistry- based probes.


Figure 1 Top View of an example of PIE
1: Sample Loading chamber, 2: Overflow chamber,
3: Homogenization chamber, 4: Dilution/Lysis chamber,
5: Three PCR wells (A to C from left to right) and one Waste Chamber (at the right end).

The Revogene can process from one up to a maximum of eight samples simultaneously in the same run. The carousel must contain eight PIEs to maintain thermodynamic balance within the run. On completion of a run, the results are interpolated by the system from measured fluorescent signals and embedded calculation algorithms. Results that are displayed on the touchscreen may be printed, transferred and/or stored by the user using the USB port or the connectivity option.

REAGENTS AND MATERIALS

The GBS LB kit contains sufficient reagents and material to process 24 specimens. The kit contains the following materials:

- 24 **Sample Buffer Tubes (SBT)**: Barcode-labeled tube containing TE 1X buffered solution (Tris-HCl pH 8.0/EDTA.Na₂) as a dilution and preservation buffer for sample.
- 24 **Disposable Transfer Tools (DTT)**: DTT consists of a single-use transfer pipette for transferring the sample from the SBT to the PIE.
- 24 individual pouches containing one (1) **GBS LB PIE**: Integrated device which comprises dried reagents allowing sample process and real-time PCR steps for GBS amplification and detection. Each PIE contains PrC, PrC-specific primers and probe, GBS *cbf* gene-specific primers and probe, dNTPs, buffer and DNA polymerase.

MATERIALS/EQUIPMENTS REQUIRED BUT NOT PROVIDED

1. Revogene® (cat# 610210)
2. Disposable gloves; powderless
3. Vortex Mixer (cat# 131132 or an equivalent)
4. Sample rack (cat# 132539; optional)
5. Swab compatible with vaginal/rectal specimen collection and recommended transport media (e.g. Liquid Stuart or Liquid Amies)
6. Selective enrichment LIM broth (Todd Hewitt broth supplemented with colistin (10 µg/mL) and nalidixic acid (15 µg/mL) (Alere cat# L57 or equivalent)
7. Calibrated micropipette (P20 recommended, VWR cat# 89079-964 or equivalent)
8. DNase/RNase-free, aerosol resistant extended length micropipette tips (Sarstedt cat# 70.1189.215 or equivalent)
9. MOCK PIE (cat# 610208; optional)

WARNING AND PRECAUTIONS

1. This product can only be used with the Revogene.
2. Do not use the kit if the label that seals the outer box is broken upon arrival.
3. Do not use GBS LB reagents if the protective pouches are open or broken upon arrival.
4. Do not interchange DTT, SBT and PIE between kit lots.
5. Each single-use GBS LB PIE and DTT are used to process one sample at a time. Do not reuse PIE and DTT.
6. Always handle specimens as if they are infectious and in accordance with safe laboratory procedures such as those described in Biosafety in Microbiological and Biomedical Laboratories⁴ and in CLSI Document M29-A4.⁵
7. Wear disposable gloves while handling patient specimens and thoroughly wash hands afterwards.
8. The GBS LB PIE contains dried reagents. The protective pouch should not be opened until ready to perform the test.
9. Dispose of unused reagents and waste in accordance with country, federal, provincial, state and local regulations.
10. Do not open or break apart the PIE after use. The cap and the seals in the PIE prevent contamination with amplification products and/or infectious particles.
11. Do not use a PIE that has been dropped, shaken or inverted after the sample has been loaded as this may cause invalid results.
12. The GBS LB assay does not provide susceptibility results. Additional time is required to culture isolates and perform susceptibility testing as recommended for penicillin-allergic women.
13. Do not use kits or kit components that have passed their stated expiration date(s).
14. Do not refrigerate the loaded PIE.

HAZARD AND PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

STORAGE AND STABILITY

1. Vaginal/rectal swab specimens should be stored according to established guidelines.¹
2. Enriched LIM broth can be kept at 25 C for up to 2 days, or at 2-8 C for up to 3 days.
3. Store the GBS LB kit at 2-25 C. The expiration date is indicated on the kit box label. Do not open a pouch until ready to perform testing. Use the PIE within 1 hour after opening the pouch.
4. Inoculated SBT can be kept at 25 C for up to 3 days, or at 2-8 C for up to 5 days.

INSTRUCTION FOR USE

SPECIMEN COLLECTION AND TRANSPORT

Specimen type: Vaginal/rectal swab taken from antepartum woman at 35-37 weeks of gestation.

Vaginal/rectal sample collection from the same patient should be performed in accordance with published guidelines for collection of clinical specimens for culture of GBS.¹

1. Collect the vaginal-rectal swab according to CDC recommendations.¹ Briefly:
 - a. Swab the lower vagina (vaginal introitus), followed by the rectum (i.e. insert swab through the anal sphincter) using the same swab or two different swabs. Take care not to mix stool, water, urine or soap in the sample.
 - b. Place the swab(s) into a non-nutritive transport medium (e.g. Liquid Stuart or Liquid Amies) right after sample collection. If vaginal-rectal swabs are collected separately from the same patient, both swabs can be placed in the same transport container.
2. Label the transport container with the sample or patient ID and transport to the laboratory for testing (refer to the "Storage and Stability" section).

SPECIMEN ENRICHMENT

3. Remove vaginal/rectal swab(s) from the transport medium and inoculate into selective enrichment LIM broth.
4. Incubate inoculated LIM broth at 35-37 C for 18-24 hours at ambient air or 5% CO₂.
5. Enriched LIM broth can be kept at 25 C for up to 2 days or at 2-8 C for up to 3 days.
6. Proceed to Specimen Preparation.

SPECIMEN PREPARATION

NOTE 1: Start the test within 1 hour after opening the pouch containing the PIE.

1. For each specimen to be tested, obtain one SBT from the kit box.
2. Label the SBT with the appropriate specimen identification without obscuring or writing over the label barcodes.
3. Vortex the enriched LIM broth specimen at high speed for 15 seconds.
4. Using a calibrated micropipette and an aerosol resistant extended length tip, aspirate 15 µL of the enriched LIM broth specimen.
5. Remove the cap from the SBT and dispense the enriched specimen aliquot into the SBT, taking care not to aerosolize the specimen. Pipet liquid up and down to ensure complete transfer of the specimen aliquot.
6. Close the SBT cap tightly, and place it on the Sample rack, if used. Make sure that only one SBT is open at once.
7. Prepare any additional specimens for testing by repeating steps 1 to 6 then proceed to step 8.
8. Mix the SBT for 15 seconds at maximum speed using a vortex mixer.
9. Unseal the PIE pouch and remove the PIE.
10. Obtain one DTT from the kit box and aspirate the SB by squeezing the entire bulb. The liquid level into the DTT must be anywhere between the two marks (Figure 2). If the liquid level is not between the two marks, discharge completely the SB in the SBT and repeat step 10.
11. Discharge completely the SB into the sample injection chamber of the PIE (Figure 1).
12. Close the cap of the PIE tightly, making sure the cap seals the chamber completely. Place the PIE on the Sample rack, if used. Do not refrigerate the loaded PIE.
13. Prepare all additional samples for testing by repeating steps 8 to 11, then proceed to step 13, making sure only one PIE is open at once.
14. Store the enriched LIM broth specimens and the inoculated SBT according to the "Storage and Stability" section.


Figure 2.
Representation of an appropriate Sample Buffer (SB) level using the Disposable Transfer Tool (DTT)

REVOGENE SYSTEM OPERATION

NOTE 1: A maximum of eight samples can be processed simultaneously in a single run (including external controls).

NOTE 2: Each run must be performed with eight PIEs in the Revogene. When less than eight samples are processed, fill the empty places with MOCK PIEs*.

1. Ensure that the Revogene is powered on. The software will launch automatically. Refer to the Revogene Operator's Manual³ for further information regarding instrument set-up and operation.
2. Enter Username and Password and press **<Login>**.
3. Press the **<Setup Run>** menu.
4. Enter the sample identification number using either the barcode scanner or manual entry.
5. Enter the SBT and PIE barcodes manually or using the barcode scanner by gently positioning the PIE almost vertically in front of the scanner. Handle the PIE carefully without shaking or dropping it.
6. (Optional) Enter the patient identification number using either the barcode scanner or manual entry.
7. (Optional) Press **<Add Comments>** and type a comment.
8. Insert the PE into the instrument, at any position of the carousel, and press **<OK>**. The software will automatically associated sample and SBT to the correct PIE.
9. To add another sample, press **<+Add>** and repeat steps 4 to 8.
10. When all samples are placed on the carousel, press **<Next>**.
11. Scan the retention ring and proceed to its installation. Close the instrument lid.
12. Initiate the test run by pressing **<Start>**.

*If MOCK PIEs are not available, use assay PIEs filled with SB or with External Controls.

VIEWING AND EXPORTING RESULTS

1. Once the run is completed, the lid unlocks automatically.
2. If the Revogene has logged out, re-enter Username and Password and press **<Login>**.
3. Press the **<Results>** menu to access test results.
4. Press **<Last Run>** to see the latest test results.
5. From the **<Last Run>** menu, select samples for which results report(s) has (have) to be exported. All samples can be selected at one time by clicking the first box above.
6. Press **<Export>** and save where appropriate (e.g. USB key).
7. Remove the retention ring and PIEs from the Revogene.

NOTE: When Indeterminate (IND) or Unresolved (UNR) results are obtained, or when an external control failure occurs, a repeat test from the prepared SBT must be performed within the specified timeframe (refer to the Repeat Testing Procedure section).

REPEAT TESTING PROCEDURE

NOTE: Only one repeat testing from the original corresponding enriched LIM broth specimen or the corresponding inoculated SBT is allowed.

1. If retesting is performed directly from the remaining enriched LIM broth specimen stored at 25 C for up to 2 days, or at 2-8 C for up to 3 days, refer to steps 1 to 11 of the Specimen Preparation section. Then follow the Revogene System Operation section.
2. If retesting is performed from the inoculated SBT stored at 25 C for up to 3 days, or at 2-8 C for up to 5 days following sample preparation, vortex the sample for 15 seconds at maximum speed using a vortex mixer, and follow steps 10 and 11 of the Specimen Preparation section for each sample. Then follow the Revogene System Operation section.



QUALITY CONTROL

Quality control procedures monitor the accuracy and precision of the analytical process. Each laboratory must establish the number, type and frequency of testing control materials per applicable regulations or accrediting agencies. The procedure described below may be employed, if appropriate, based on local policies and procedures.

1. Each PIE contains a Process Control (PrC) that monitors for sample homogenization, sample dilution, cell lysis, DNA amplification inhibition and assay reagents failure.
2. Good laboratory practice recommends the use of control materials. Users should follow the appropriate guidelines concerning running of external quality controls. It is recommended that one Positive External Control and one Negative External Control should be run at least on a daily basis until adequate process validation is achieved on the Revogene in each laboratory setting. External control materials are not provided by Meridian Bioscience, Inc.
3. Commercial control material (e.g. *Streptococcus agalactiae* ATCC® 13813) can be used as a Positive External Control. It is recommended that bacterial strains be freshly prepared in LIM broth. Transfer a 15 µL aliquot of a 18-24 hours culture in LIM broth into a SBT and proceed from step 6 of the Specimen Preparation section.
4. Pure LIM broth is recommended for use as a Negative External Control. Transfer a 15 µL aliquot of LIM broth into an SBT and proceed from step 6 of the Specimen Preparation section.
5. Separate DTT, SBT and PIE must be used for each external control reagent.

RESULTS INTERPRETATION

The results are computed by the Revogene from measured fluorescent signals and embedded calculation algorithms and are available on the Results window of the Revogene. Possible results are:

Sample	Symbol displayed on user screen	Overall reported result	Interpretation
Patient Specimen		Positive	Sample contains GBS target DNA.
		Negative	No GBS target DNA detected or number of organisms may be below the limit of detection of the assay.
		Indeterminate	No reportable result due to possible Revogene detection error during the assay treatment, or the data analysis, or if the run is interrupted by the user. Repeat testing must be performed using the original corresponding enriched LIM broth specimen or corresponding inoculated SBT within the timeframe defined above.
		Unresolved	Amplification/detection failure for the Process Control as well as for the GBS target. Could be caused by inhibitory specimens, microfluidic or reagent failure. Repeat testing must be performed using the original corresponding enriched LIM broth specimen or corresponding inoculated SBT within the timeframe defined above.
Positive External Control		Positive	Valid Positive External Control result.
		Negative	Incorrect Positive External Control result. The run is invalid. Repeat testing must be performed for all samples of the run, using the original corresponding enriched LIM broth specimen or corresponding inoculated SBT with a new set of external controls within the timeframe defined above.
		Indeterminate	Incorrect Positive External Control result. No reportable result due to possible Revogene detection error during the assay treatment, or the data analysis, or if the run is interrupted by the user. The run is invalid. Repeat testing must be performed for all samples of the run, using the original corresponding enriched LIM broth specimens or corresponding inoculated SBTs with a new set of external controls within the timeframe defined above.
		Unresolved	Incorrect Positive External Control result. Amplification/detection failure for the Process Control as well as for the GBS target. Could be caused by inhibitory specimens, microfluidic or reagent failure. The run is invalid. Repeat testing must be performed for all samples of the run, using the original corresponding enriched LIM broth specimens or corresponding inoculated SBTs with a new set of external controls within the timeframe defined above.
Negative External Control		Positive	Incorrect Negative External Control result due to specimen handling and/or a contamination problem. The run is invalid. Repeat testing must be performed for all samples of the run, using the original corresponding enriched LIM broth specimens or corresponding inoculated SBTs with a new set of external controls within the timeframe defined above.
		Negative	Valid Negative External Control result.
		Indeterminate	Incorrect Negative External Control result. No reportable result due to possible Revogene detection error during the assay treatment, or the data analysis, or if the run is interrupted by the user. The run is invalid. Repeat testing must be performed for all samples of the run, using the original corresponding enriched LIM broth specimens or corresponding inoculated SBTs with a new set of external controls within the timeframe defined above.
		Unresolved	Incorrect Negative External Control result. Amplification/detection failure for the Process Control. Could be caused by inhibitory specimens, microfluidic or reagent failure. The run is invalid. Repeat testing must be performed for all samples of the run, using the original corresponding enriched LIM broth specimens or corresponding inoculated SBTs with a new set of external controls within the timeframe defined above.

LIMITATIONS OF THE PROCEDURE

- The GBS LB assay can only be used with the Revogene by trained personnel.
- A positive test result does not necessarily indicate the presence of viable organisms. It is, however, indicative for the presence of GBS DNA.
- Performance of the GBS LB assay was established with vaginal/rectal swab specimens collected from antepartum women, transported in non-nutritive medium (e.g., Liquid Stuart) and enriched in selective LIM broth. Use of the GBS LB assay for clinical specimen types other than those specified has not been evaluated and performance characteristics are not established.
- Cervical, perianal, perirectal or perineal specimens are not acceptable, and a speculum should not be used for culture collection.
- The GBS LB assay does not provide susceptibility results. Antimicrobial susceptibility testing should be performed on antenatal GBS isolates from penicillin-allergic women at high risk for anaphylaxis.
- The use of enrichment broths different than LIM broth was not evaluated.
- Erroneous test results may occur from improper specimen collection, handling or storage, technical error, or sample mix-up. Careful compliance to the instructions in this insert and to established guidelines is necessary to avoid erroneous results.
- Erroneous assay results may occur if a PIE cap is incorrectly closed.
- While there are no known strains/isolates of GBS lacking the *cfb* gene, the occurrence of such a strain would lead to an erroneous result (false negative) using the GBS LB assay.
- Enterococcus faecalis*, *Enterococcus faecium* and *Lactobacillus acidophilus* can potentially inhibit the detection of GBS if individually present at > 10⁴ CFU/mL of SB.
- Potentially inhibitory effect of fecal fat has not been evaluated.
- A combination of loperamide hydrochloride (e.g. Imodium®), bismuth subsalicylate (e.g. Pepto-Bismol™) and sennosides (e.g. Senokot®), while highly unlikely, may potentially inhibit the detection of GBS when these substances are present in SBT at concentrations of 0.023 µg/mL, 0.675 µg/mL and 0.011 µg/mL respectively.
- Mutations or polymorphisms in primer- or probe-binding regions may affect detection of GBS *cfb* gene variants, resulting in a false negative result with the GBS LB assay.
- Results from the GBS LB assay should be used as an adjunct to clinical observations and other information available to the physician.
- A negative result does not rule out the possibility of GBS colonization. False negative results may occur when the GBS concentration is below limit of detection of the assay.
- The test is not intended to differentiate carriers of GBS from those with streptococcal disease.
- Test results may be affected by concurrent antimicrobial therapy as GBS DNA may continue to be detected.

EXPECTED VALUES

In the investigational study for the Revogene GBS LB assay, a total of 771 reportable results from specimens compliant at the culture and PCR levels were obtained from 4 geographically diverse regions. Overall, GBS prevalence rate was 22.0 % (170/771) with a 95 % CI of 19.3 – 25.1 %.

PERFORMANCE CHARACTERISTICS

CLINICAL PERFORMANCE

Performance of the GBS LB assay was established during a prospective investigational study involving two Canadian and two American sites using the Revogene.

Prospectively collected and tested specimens consisted of excess leftover deidentified enriched LIM broth specimens collected for standard of care antepartum GBS testing from women at 35-37 weeks of gestation. The leftover LIM broth was used for the GBS LB assay and for the reference culture testing performed on each site.

The GBS LB assay performance was compared to a composite reference method consisting of LIM broth enrichment of the vaginal/rectal swab followed by a subculture on sheep blood agar plates and biochemical identification of GBS. All specimens with a negative bacterial culture result were re-tested at a single independent core laboratory using the same standard operating procedures required for bacterial culture.

A total of 839 specimens were enrolled in the study. Of those, 63 specimens were regarded as noncompliant per reference method and/or GBS LB assay protocol criteria, 1 specimen was regarded as noncompliant in terms of transport and storage conditions and 4 specimens were regarded as noncompliant due to missing test results. Two (2) fully compliant specimens gave final non-reportable PCR results. A total of 771 specimen results were used to determine the clinical performance of the GBS LB assay in comparison to the composite reference method (Tables 1 and 2).

Overall performance characteristics of the GBS LB assay demonstrated 95.9% sensitivity (163/170, 95%CI 91.7%- 98.0%) and 95.5% specificity (574/601), 95%CI 93.5% - 96.9%) in comparison to the composite reference method. Of the 812 specimens tested with the GBS LB assay that were compliant at the specimen and PCR level, 6 (0.74%) were reported as Unresolved after initial testing. The unresolved rate after repeat testing was 0.25% (2/812).

Out of the same 812 specimens tested with the GBS LB assay, 12 (1.48%) were initially reported as Indeterminate.; 6 from a run with an external control failure and 6 from an incomplete run. No result remained Indeterminate upon repeat testing.

Table 1. Overall performance characteristics (all sites combined) of the GBS LB assay in comparison to the Composite Reference Method.

Overall performance		Composite Reference Method		Total
GBS LB assay	Positive	Positive	Negative	
	Negative	163	27 ^B	190
	Total	7A	574	581
		170	601	771

Sensitivity: 95.9% (95%CI: 91.7 - 98.0 %)
Specificity: 95.5% (95%CI: 93.5 - 96.9 %)

^A 5 out of 7 false negative Revogene GBS LB results were tested on a FDA-cleared molecular device and yielded negative results

^B 10 out 27 false positive Revogene GBS LB results were tested on a FDA-cleared molecular device and yielded positive results.

Table 2. Summary of performance characteristics of the GBS LB assay per site

Site	Sensitivity	Specificity	Prevalence ¹
1	87.8% (36/41)	92.8% (142/153)	21.1% (41/194)
2	98.0% (48/49)	97.5% (194/199)	19.8% (49/248)
3	100% (43/43)	97.0% (164/169)	20.3% (43/212)
4	97.3% (36/37)	92.5% (74/80)	31.6% (37/117)
Total	95.9% (163/170)	95.5% (574/601)	22.0% (170/771)

¹Prevalence was calculated with samples compliant at culture reference method and Revogene GBS LB assay levels.

ANALYTICAL SENSITIVITY

The analytical sensitivity (Limit of Detection or LoD) of the GBS LB assay was determined using clinical LIM broth matrix previously tested negative for GBS and spiked with different concentrations of GBS bacterial suspensions. Two strains of GBS (ATCC 12403 and ATCC 13813) were tested in 24 replicates per concentration by 2 operators using 3 different lots of GBS LB kits. The LoD is defined as the lowest concentration at which 95% of all replicates tested positive. The LoD of the GBS LB assay with the two strains tested ranged from 200 to 375 CFU/mL of SB (Table 3).

Table 3. LoD of the GBS LB assay

GBS strain (ATCC number)	LoD (CFU/mL of SB)
Serotype III (ATCC 12403)	375
Non-hemolytic (ATCC 13813)	200

INCLUSIVITY

Inclusivity of the GBS LB assay was determined for 12 strains of GBS representing 11 known serotypes and 1 non-hemolytic strain. Each GBS strain was tested using a quantitated culture spiked in a GBS-negative clinical LIM broth matrix. Twenty-four replicates per concentration per strain were tested using 3 different lots of GBS LB kits. The lowest concentrations at which strains reached a positivity rate of 100% are summarized in Table 4.

Table 4. Organisms tested for inclusivity with the Revogene GBS LB assay.

GBS strain	Test concentration	Test concentration (CFU/mL of SB)	Positivity rate
Serotype Ia (ATCC 12400)	< 2	735	100%
Serotype Ib (ATCC 51487)	15	5,625	100%
Serotype Ic (ATCC 27591)	5	1,875	100%
Serotype II (ATCC 12973)	< 2	735	100%
Serotype III (ATCC 12403)	1	375	100%
Serotype IV (ATCC 49446)	7	2,625	100%
Serotype V (ATCC BAA-611)	5	1,875	100%
Serotype VI (ATCC BAA-2671)	3	1,125	100%
Serotype VII (ATCC BAA-2670)	7	2,625	100%
Serotype VIII (ATCC BAA-2669)	10	3,750	100%
Serotype IX (ATCC BAA-2668)	3	1,125	100%
Non-hemolytic (ATCC 13813)	1.33	500	100%

*LoD of strain ATCC 12403 determined to be 375 CFU/mL of SB.

In addition, an *in silico* analysis was performed to assess inclusivity of the GBS LB assay regarding 36 additional GBS strains listed in the National Center for Biotechnology Information (NCBI) database of October 20th, 2016. The alignment results showed no mismatch with the selected sequences. The analysis predicted the detection of all these GBS strains.

CROSS-REACTIVITY

A cross-reactivity study was performed to assess unintended reactivity of the GBS LB assay with high loads of organisms phylogenetically related to GBS or present in urogenital flora and intestinal tract and not targeted by the test. The study included 64 bacteria, 4 yeasts, 4 viruses, 2 parasites and human DNA (Table 5). Bacteria and yeasts were tested at a load of $\geq 10^6$ CFU/mL of SB. Nucleic acids from viruses, parasites and human DNA were tested at a load of $\geq 10^5$ DNA or RNA cp/mL of SB. Each organism was tested using a quantitated culture or nucleic acid solution spiked in a GBS-negative clinical LIM broth matrix. Three replicates per organism were tested using 3 different lots of GBS LB kits.

Under the conditions of the study, none of the organisms or nucleic acids tested were detected by the GBS LB assay.

Table 5. List of organisms tested *in vitro* for cross-reactivity with the GBS LB assay

Bacteria	
<i>Acinetobacter baumannii</i>	<i>Mycoplasma genitalium</i> (gDNA)
<i>Aerococcus viridans</i>	<i>Mycoplasma hominis</i> (gDNA)
<i>Aeromonas hydrophila</i>	<i>Neisseria gonorrhoeae</i>
<i>Bacillus cereus</i>	<i>Peptostreptococcus anaerobius</i>
<i>Bacillus subtilis</i>	<i>Porphyromonas asaccharolytica</i>
<i>Bacteroides fragilis</i>	<i>Prevotella melaninogenica</i>
<i>Bifidobacterium adolescentis</i>	<i>Propionibacterium acnes</i>
<i>Bifidobacterium breve</i>	<i>Proteus mirabilis</i>
<i>Brevibacterium linens</i>	<i>Pseudomonas aeruginosa</i>
<i>Campylobacter jejuni</i>	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Dublin</i>
<i>Chlamydia trachomatis</i> (gDNA)	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Minneapolis</i>
<i>Citrobacter freundii</i>	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>typhimurium</i>
<i>Clostridium difficile</i>	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Newport</i>
<i>Clostridium perfringens</i>	<i>Serratia marcescens</i>
<i>Corynebacterium genitalium</i>	<i>Shigella sonnei</i>
<i>Enterobacter aerogenes</i>	<i>Staphylococcus aureus</i>
<i>Enterobacter cloacae</i>	<i>Staphylococcus aureus</i> (Cowan)
<i>Enterococcus avium</i>	<i>Staphylococcus epidermidis</i>
<i>Enterococcus dispar</i>	<i>Staphylococcus saprophyticus</i>
<i>Enterococcus durans</i>	<i>Streptococcus anginosus</i>
<i>Enterococcus faecalis</i>	<i>Streptococcus bovis</i>
<i>Enterococcus faecium</i>	<i>Streptococcus dysgalactiae</i> subsp. <i>disgalactiae</i>
<i>Escherichia coli</i>	<i>Streptococcus dysgalactiae</i> <i>equisimilis</i>
<i>Escherichia fergusonii</i>	<i>Streptococcus intermedius</i>
<i>Gardnerella vaginalis</i>	<i>Streptococcus oralis</i>
<i>Klebsiella oxytoca</i>	<i>Streptococcus pneumoniae</i>
<i>Klebsiella pneumoniae</i>	<i>Streptococcus pyogenes</i>
<i>Lactobacillus acidophilus</i>	<i>Streptococcus salivarius</i>
<i>Lactobacillus brevis</i>	<i>Streptococcus sanguinis</i>
<i>Lactobacillus casei</i>	<i>Streptococcus uberis</i>
<i>Lactobacillus delbreuckii</i>	<i>Ureaplasma urealyticum</i> (gDNA)
<i>Lactobacillus jensenii</i>	<i>Yersinia enterocolitica</i>
Yeasts	
<i>Candida albicans</i>	<i>Candida parapsilosis</i>
<i>Candida glabrata</i>	<i>Candida tropicalis</i>
Viruses	
<i>HerpesSimplexVirus-1</i> (gDNA)	<i>Norovirus</i> GII (RNA)
<i>HerpesSimplexVirus-2</i> (gDNA)	Human Papillomavirus (HPV) (gDNA)
Parasites	
<i>Blastocystis hominis</i> (gDNA)	<i>Trichomonas vaginalis</i> (gDNA)
Human DNA	
Human DNA (gDNA)	

In addition, an *in silico* analysis was performed to assess the risk of cross-reactivity of the GBS LB assay towards 62 additional organisms potentially found in urogenital flora and intestinal tract (Table 6). None of the organisms are presumed to cross-react with the GBS LB assay based on the high degree of mismatches between the selected NCBI gene sequences and the test primers and probes.

Table 6. List of organisms tested *in silico* for cross-reactivity with the GBS LB assay

Bacteria	
<i>Abiotrophia defectiva</i>	<i>Listeria monocytogenes</i>
<i>Achromobacter xylosoxidans</i>	<i>Micrococcus luteus</i>
<i>Alcaligenes faecalis</i>	<i>Mobiluncus mulieris</i>
<i>Campylobacter fetus</i>	<i>Moraxella catarrhalis</i>
<i>Chromobacterium violaceum</i>	<i>Moraxella lacunata</i>
<i>Clostridium bifementans</i>	<i>Moraxella osloensis</i>
<i>Clostridium butyricum</i>	<i>Morganella morganii</i>
<i>Clostridium novyi</i>	<i>Mycobacterium smegmatis</i>
<i>Clostridium septicum</i>	<i>Neisseria flava</i>
<i>Clostridium sordelli</i>	<i>Neisseria flavescens</i>
<i>Clostridium sporogenes</i>	<i>Neisseria lactamica</i>
<i>Corynebacterium urealyticum</i>	<i>Neisseria meningitidis</i>
<i>Corynebacterium xerosis</i>	<i>Neisseria perflava</i>
<i>Deinococcus radiodurans</i>	<i>Plesiomonas shigelloides</i>
<i>Derrxia gummosa</i>	<i>Proteus vulgaris</i>
<i>Eikenella corrodens</i>	<i>Providencia stuartii</i>
<i>Enterococcus gallinarum</i>	<i>Pseudomonas fluorescens</i>
<i>Gemella haemolysans</i>	<i>Pseudomonas putida</i>
<i>Haemophilus ducreyi</i>	<i>Rhodospirillum rubrum</i>
<i>Haemophilus influenza</i>	<i>Serratia liquefaciens</i>
<i>Haemophilus influenza type B</i>	<i>Shigella boydii</i>
<i>Helicobacter pylori</i>	<i>Shigella flexneri</i>
<i>Kingella denitrificans</i>	<i>Staphylococcus simulans</i>
<i>Kingella kingae</i>	<i>Stenotrophomonas maltophilia</i>
<i>Legionella pneumophila</i>	<i>Streptococcus gordonii</i>
<i>Streptococcus mitis</i>	<i>Vibrio parahaemolyticus</i>
<i>Streptococcus mutans</i>	
Mold/Yeast	
<i>Cryptococcus neoformans</i>	<i>Saccharomyces cerevisiae</i>
Viruses	
<i>Adenovirus 40</i>	<i>Echovirus</i>
<i>Adenovirus 41</i>	<i>Epstein-Barr Virus</i>
<i>BK virus</i>	<i>Rotavirus</i>
<i>Coxsackievirus</i>	

INTERFERING ORGANISMS

A study was conducted to assess the potentially inhibitory effect of 29 organisms that may be present in urogenital flora and intestinal tract and which are not targeted by the test. Selection of organisms was based on their prevalence in the vaginal/rectal flora. The *Enterococcus* strains represent the highest interference risk due to their prevalence in the rectal flora as well as their capacity to grow in LIM broth. Other organisms present a low risk of being found at high loads in the specimen since selective LIM broth enrichment inhibits their growth and the loads of those organisms cannot increase much following sampling. Each organism or nucleic acid solution was diluted to reach a final concentration of 10⁶ CFU/mL or cp/mL of SB, except *Escherichia coli* ATCC 11775, *Lactobacillus acidophilus* ATCC 4356, *Trichomonas vaginalis* ATCC 30001 and *Mycoplasma genitalium* ATCC BAA-2641S, which were tested at 10⁵ CFU/mL or cp/mL of SB. Fifteen (15) µL of a GBS dilution prepared in a negative clinical LIM broth matrix to reach a final concentration of 735 CFU/mL of SB was added to the SBT. The 29 organisms included in the study are presented in Table 7.

Table 7. List of organisms tested for interference with the GBS LB assay

Organisms	
<i>Enterococcus faecalis</i>	<i>Neisseria gonorrhoeae</i>
<i>Enterococcus faecium</i>	<i>Peptostreptococcus anaerobius</i>
<i>Acinetobacter baumannii</i>	<i>Propionibacterium acnes</i>
<i>Bacillus subtilis</i>	<i>Pseudomonas aeruginosa</i>
<i>Bacteroides fragilis</i>	<i>Salmonella enterica subsp. enterica serovar Dublin</i>
<i>Campylobacter jejuni</i>	<i>Serratia marcescens</i>
<i>Clostridium difficile</i>	<i>Shigella sonnei</i>
<i>Corynebacterium genitalium</i>	<i>Staphylococcus aureus</i>
<i>Escherichia coli</i>	<i>Staphylococcus epidermidis</i>
<i>Gardnerella vaginalis</i>	<i>Streptococcus pyogenes</i>
<i>Klebsiella oxytoca</i>	<i>Yersinia enterocolitica</i>
<i>Lactobacillus acidophilus</i>	<i>Candida albicans</i>
<i>Chlamydia trachomatis</i> (gDNA)	<i>Trichomonas vaginalis</i> (gDNA)
<i>Herpes Simplex Virus-1</i> (gDNA)	<i>Mycoplasma genitalium</i> (gDNA)
Human Papillomavirus (HPV) (gDNA)	

None of the 29 organisms present at ≥ 10⁵ CFU/mL or cp/mL of SB in the specimen interfered with detection of PrC. Among the 29 organisms tested, only 3 showed a potentially inhibitory effect on detection of GBS (Table 8). Inhibition from *E. faecalis*, *E. faecium* and *L. acidophilus* was observed at concentrations of > 10⁵, > 10⁴ and > 10⁴ CFU/mL of SB respectively.

Table 8. Organisms interfering with the GBS LB assay.

Organism (ATCC)	Concentration tested (CFU/mL of SB)	GBS positivity rate (positive / total)	
		Non-hemolytic	Serotype III
<i>Enterococcus faecalis</i> (ATCC 19433)	106	66.7% (2/3)	33.3% (1/3)
	105	100% (3/3)	100% (3/3)
	106	100% (3/3)	66.7% (2/3)
<i>Enterococcus faecium</i> (ATCC 19434)	105	N/D*	33.3% (1/3)
	104	N/D	100% (3/3)
	105	100% (3/3)	66.7% (2/3)
<i>Lactobacillus acidophilus</i> (ATCC 4356)	104	N/D	100% (3/3)

*N/D: Not determined because there was no interference at the highest concentration.

INTERFERING SUBSTANCES

A study was conducted to assess the potentially inhibitory effect of 22 exogenous and 9 endogenous substances that may be present in urogenital flora and intestinal tract. GBS negative and GBS positive samples at 735 CFU/mL of SB were tested in presence of a negative clinical LIM broth matrix with the highest physiologically relevant concentration of substance. The 31 substances tested in 11 groups are presented in Table 9.

Results demonstrated no reportable interference on PrC. Only the combination of loperamide hydrochloride (e.g. Imodium®), bismuth subsalicylate (e.g. Pepto-Bismol™) and sennosides (e.g. Senokot®) at concentrations of 0.023 µg/mL, 0.675 µg/mL and 0.011 µg/mL SB respectively showed a potentially inhibitory effect on detection of GBS. When tested individually, these substances showed no reportable interference with the GBS LB assay.

Table 9. List of exogenous and endogenous substances tested with the GBS LB assay.

Exogenous substances			
Group	Substance (commercial name)	Concentration in SBT ¹	Results ²
A	Fungicide (Micatin [®])	0.023% W/V	NI
	Hemorrhoid cooling gel (Preparation H [®])	0.023% W/V	
	Lubricating gel (K-Y [®] Personal Lubricant)	0.023% W/V	
	Body powder (Vagisil [®] deodorant powder)	0.023% W/V	
	Moisturizing lotion (Aveeno [®] moisturizing lotion)	0.023% W/V	
	Body oil (Neutrogena [®] Body oil)	0.023% V/V	
B	Deodorant spray (Summer's Eve [®] Spray)	0.023% V/V	NI
	Enemas (Life BRAND [™] Heavy Mineral Oil USP)	0.023% V/V	
	Antimicrobials (Canesten [®])	0.023% W/V	
	Enemas (Pentasa)	0.495 µg/mL	
C	Radiology oral compounds (Barium sulfate)	0.113 µg/mL	NI
	Gastritis medications (Nexium)	0.011 µg/mL	
	Antimicrobials (Flagyl)	0.016 µg/mL	
D	Non-steroidal anti-inflammatory medications (Aleve [®])	0.071 µg/mL	NI
	Antimicrobials (DIFLUCAN [®] One)	0.020 µg/mL	
	Gastritis medications (Tums [®])	1.200 µg/mL	
E	Anti-diarrheal medication (Imodium [®])	0.023 µg/mL	I
	Anti-diarrheal medication (Pepto-Bismol [®])	0.675 µg/mL	
	Laxatives (Senokot [®])	0.011 µg/mL	
F	Spermicide (Trojan [®] with spermicidal lubricant)	0.023% V/V	NI
	Moist towelettes (Equate [™] flushable moist wipes)	0.023% V/V	
	Moist towelettes (Wet Ones [®])	0.023% V/V	
Endogenous substances			
Group	Substance	Concentration in SBT ¹	Results ²
G	Whole blood	0.023% V/V	NI
	Leukocytes	15,000 cells/mL	
H	Amniotic fluid	0.023% V/V	NI
	Mucous	0.023% V/V	
I	Seminal fluid	0.023% V/V	NI
	Urine	0.023% V/V	
J	Feces	0.023% V/V	NI
	Meconium	0.023% V/V	
K	Human DNA	4.65 ng/mL	NI

¹W/V: Weight/Volume; V/V: Volume/Volume

²I: interference with the GBS

LB assay;

NI: no interference with the GBS LB assay

PRECISION/REPRODUCIBILITY

A reproducibility study was conducted to evaluate the within- site (precision), between-site and between-day variance with multiple operators. The panel comprised three GBS strains (ATCC 13813 (non-hemolytic), ATCC 12400 (serotype Ia) and ATCC 12403 (serotype III)) tested at two concentrations and two true negative samples. Each sample was tested by two operators three times per day over five different days on three sites. A total of 90 replicates for each positive sample and of 180 replicates for negative samples were tested (8 panel members x 2 operators x 3 times per day x 5 days x 3 sites). The study was performed using three lots of GBS LB assay.

The panel consisted of three contrived samples categories as follow:

1. Low positive (LP): $1 < 2 \times \text{LoD}$
2. Moderate positive (MP): $> 2 - 3 \times \text{LoD}$
3. True negative (TN): samples without target

Precision study results within sites are shown in Table 10.

Table 10. Summary of the precision study within sites

Site 1				
Panel ID	Results/Total	% Agreement	Ct Mean	Ct %CV
LP	87/90	96.7%	37.3	5.4
MP	87/90	96.7%	36.6	5.0
TN	60/60	100%	32.1	3.7
Site 2				
Panel ID	Results/Total	% Agreement	Ct Mean	Ct %CV
LP	90/90	100%	36.3	4.2
MP	89/90	98.9%	35.7	3.8
TN	60/60	100%	32.7	3.2
Site 3				
Panel ID	Results/Total	% Agreement	Ct Mean	Ct %CV
LP	90/90	100%	35.5	4.1
MP	90/90	100%	34.9	4.3
TN	60/60	100%	30.8	2.4

Reproducibility study results between sites are shown in Table 11.

Table 11. Summary of the reproducibility study results between sites

All sites combined				
Panel ID	Results/Total	% Agreement	Ct Mean	Ct %CV
LP	267/270	98.9%	36.4	2.3
MP	266/270	98.5%	35.7	2.3
TN	180/180	100%	31.9	3.1

CARRY-OVER AND CROSS-CONTAMINATION

A study was conducted to investigate within run and across run carry-over and cross-contamination. Positive samples contained 15 µL of a GBS dilution prepared in a clinical negative LIM broth matrix to reach a final concentration of $> 10^6$ CFU/mL of SB. Negative samples contained 15 µL of clinical negative LIM broth matrix per SBT. For the within run study, 4 replicates of high positive samples and 4 replicates of negative samples were tested in each run by alternating positive and negative replicates. For the across run study, a run of 8 replicates of high positive samples followed by a run of 8 replicates of negative samples were repeated 5 times. Each operator performed 10 runs of 8 samples, for a total of 20 runs.

Overall, all high positive samples were accurately detected as positive while all of the true negative samples were accurately detected as negative. Absence of carry-over and cross-contamination was demonstrated when using the GBS LB assay.

E-LABELING

Documentation related to this product can be accessed online at www.meridianbioscience.com/pi. Additionally, paper copies are available upon request by contacting your local distributor or via the phone number listed on the kit box.

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REFERENCES

1. MMWR Morbidity and Mortality Weekly Report, Recommendations and Reports, Prevention of Perinatal Group B Streptococcal Disease; Revised Guidelines from CDC, 2010 Nov 19; 59(RR-10):1-36
2. A population-based comparison of strategies to prevent early-onset Group B Streptococcal disease in neonates. N Engl J Med 2002; 347:233-239
3. SN134822 Revogene® Operator's Manual
4. Biosafety in microbiological and biomedical laboratories, 5th ed. U.S Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institute of Health, HHS Publication No. (CDC) 21-1112, Revised 2009 Dec
5. M29-A4 Protection of laboratory workers from occupationally acquired infections; Approved Guideline 4th ed. Clinical and Laboratory Standards Institute (CLSI) 2014 May.

SN134763

REV. 04/26

 <p>Manufactured By</p>	<p>Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com</p> <p><u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com</p>
--	--

INTERNATIONAL SYMBOLS USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Zeichenerklärung)

	Use-by date / Data di scadenza / Date de péremption / Fecha de caducidad / Verwendbar bis
	Batch code / Codice di lotto / Code de lot / Código de lote / Chargennummer
	In vitro diagnostic medical device / Dispositivo medico diagnostico <i>in vitro</i> / Instrument de test diagnostique <i>in vitro</i> / Producto sanitario para diagnóstico <i>In vitro</i> / <i>In-vitro</i> -Diagnostikum
	CE Mark / Marcatura CE / Symbole CE / Marcado CE / CE-Kennzeichen
	Catalog number / Numero di catalogo / Référence catalogue / Referencia / Bestellnummer
	Consult instructions for use / Consultare le istruzioni per l'uso / Consulter le mode d'emploi / Consultar las instrucciones de uso / Gebrauchsanweisung beachten
	Manufacturer / Produttore / Fabricant / Fabricante / Hersteller
	Contains sufficient for <n> tests / Contiene una quantità sufficiente per <n> test / Contient le matériel suffisant pour <n> tests / Contiene la cantidad suficiente para <n> ensayos / Inhalt ausreichend für <n> Tests
	Temperature limit / Limite di temperatura / Limite de température / Limite de temperatura / Temperaturgrenze
	Authorized representative in the European Community / Rappresentante autorizzato nella Comunità Europea / Représentant agréé dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter EU-Repräsentant
	Do not reuse / Non riutilizzare / Ne pas réutiliser / No reutilizar / Nicht wiederverwenden
	Keep dry / Conservare all'asciutto / Conserver au sec / Mantener seco / Vor Feuchtigkeit schützen
	Contains # pouches: 1 Disposable Transfer Tool (DTT), 1 Sample Buffer Tube (SBT), 1 PIE / Contiene # buste: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT) / Contient # sachets: Une cartouche PIE, Un tube échantillon, Un outil de transfert jetable (OTJ) / Incluye # bolsitas: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT) / Enthält # Beutel: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT)
	Humidity Limitation / Limitazione dell'umidità / Limite d'humidité / Limitación de humedad / Feuchtigkeitsbegrenzung
	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Ne pas utiliser si l'emballage est endommagé / No utilizar si el envase está dañado / Nicht verwenden, wenn die Verpackung beschädigt ist
	Keep away from sunlight / Tenere lontano dalla luce del sole / Tenir à l'écart de la lumière du soleil / Mantener alejado de la luz solar / Vor Sonnenlicht schützen
	Contains 24 Disposable Transfer Loops (DTL) / Contiene 24 loop di trasferimento monouso (DTL) / Contient 24 boucles de transfert jetables (DTL) / Contiene 24 bucles de transferencia desechables (DTL) / Enthält 24 Einwegübertragungsschleifen (DTL)
	Disposable Transfer Tool (DTT) / Dispositivo Transfer Tool 9DTT / Un outil de transfert jetable (OTJ) / Disposable Transfer Tool (DTT) / Disposable Transfer Tool (DTT)
	Revogene Test Device / Dispositivo Test Revogene / Dispositif de test Revogene / Dispositivo para la Prueba Revogene / Revogene-Analysegefäß
	For Emergency Use Authorization only / Solo per l'autorizzazione per l'uso di emergenza / pour autorisation d'utilisation d'urgence uniquement / para autorización de uso de emergencia solamente / nur für Notfallverwendungsautorisierung
	Disposable Transfer Loop / Loop di trasferimento monouso / Boucles de transfert jetables / Bucles de transferencia desechables / Einwegübertragungsschleifen
	Sample Buffer Tube / Sample Buffer Tube / Un tube échantillon / Sample Buffer Tube / Sample Buffer Tube
	Revogene Test Device / Dispositivo Test Revogene / Dispositif de test Revogene / Dispositivo para la Prueba Revogene / Revogene-Analysegefäß

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order call Customer Service Department at 800-543-1980.

TaqMan is a registered trademark of Roche Molecular Systems, Inc.
ATCC is a trademark of the American Type Culture Collection.
Revogene and associated logos are registered trademarks of Meridian Bioscience, Inc.
© 2020-12 Meridian Bioscience, Inc.

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TaqMan è un marchio registrato di Roche Molecular Systems, Inc. ATCC
Revogene e i relativi loghi sono marchi registrati di Meridian Bioscience, Inc.
TaqMan est une marque commerciale déposée de Roche Molecular Systems, Inc.
Revogene et les logos connexes sont des marques déposées de Meridian Bioscience, Inc.
TaqMan es una marca registrada de Roche Molecular Systems, Inc.
Revogene y los logotipos asociados son marcas registradas de Meridian Bioscience, Inc.