**Para-Pak® SYSTEMS**  
**For use with Zn-PVA and Formalin**

US Patent No. 5624554

REF 901012

IVD

Rx Only

**INTENDED USE**

Para-Pak Zn-PVA based systems provide standardized procedures for the routine collection, transportation, preservation, and examination of stool specimens for intestinal parasites. Kit systems are designed for easy use by individuals not trained in microbiological procedures and afford an excellent means of minimizing the adverse effects of delay in specimen transportation.

**SUMMARY AND EXPLANATION OF THE TEST**

Diagnosis of intestinal parasitic disease is confirmed by recovery and identification of helminth eggs and larvae, or protozoan trophozoites and cysts in the clinical parasitology laboratory. Timely collection and transportation of "fresh" stool specimens to the laboratory cannot always be insured. Workload conditions and priorities in clinical laboratories frequently do not permit immediate examination of "fresh" specimens. Procedures such as incubation, refrigeration, or freezing of stool specimens will not guarantee the recovery of all diagnostic stages of all parasites.<sup>5, 6, 8, 9, 11, 12</sup> In 1949, Brooke and Goldman described the use of PVA-fixative for the preservation of intestinal protozoa. The efficacy of the two vial PVA/10% Formalin method in the transportation and preservation of stool specimens has since been confirmed by numerous authors.<sup>2-6, 8, 9</sup> Proper use of the Para-Pak Zn-PVA/Formalin systems assures the parasitologist that the diagnostic stages of intestinal parasites will be preserved if present in the fecal material.

**BIOLOGICAL PRINCIPLES**

Reagent alcohol (28%) in a stabilizing solution containing zinc sulfate and polyvinyl alcohol provides an acceptable preservative-fixative for protozoan trophozoites and cysts. The resulting permanent slide lends itself to commonly used staining procedures such as Wheatley's trichrome or iron hematoxylin.<sup>2,4, 13</sup> Ten-percent formalin preserved specimens may be examined directly or concentrated for recovery of eggs, larvae, and protozoan cysts.

**REAGENTS/MATERIALS PROVIDED**

***The maximum number of tests obtained from this test kit is listed on the outer box.***

Each kit consists of one vial containing 15 mL of Zn-PVA fixative and one vial containing 15 mL of 10% Formalin preservative. Single vial cases are also available. Simple directions for patients and nursing personnel are also provided in ten languages.

**Catalog# 901012**

Zn-PVA (item# 9012)

**MATERIALS NOT PROVIDED**

1. Ethyl acetate (suggested) or ether (optional)
2. Physiological saline
3. Zinc Sulfate Solution
4. Cotton tipped applicator sticks
5. Microscope slides and coverslips
6. Centrifuge
7. Microscope
8. Transfer pipettes

**PRECAUTIONS**

1. All reagents are for in vitro diagnostic use only.
2. Avoid contact of fixative solutions with the skin and eyes. Should contact occur, flush with running water. If irritation should develop, see a physician.
3. Fixative solutions are poisonous. If ingested, call a physician immediately.
4. Due to the infectious nature of unpreserved stools, care and handwashing should be employed when the specimen is collected and manipulated.
5. Formalin is a potential cancer hazard.
6. Zn-PVA is flammable.
7. Any serious incident that has occurred in relation to the device should be reported to Meridian Bioscience, Inc., 3471 River Hills Drive, Cincinnati, Ohio 45244 USA or Technical Support Center 800-343-3858 and competent authority of the EU Member State in which the clinician and/or patient is established.
8. IMPORTANT: See SDS for additional safety and hazard information.

**HAZARD and PRECAUTIONARY STATEMENTS**

|                                                                                                                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  <p style="text-align: center;">Para-Pak Zn-PVA Fixative</p> | <p><b>Signal Word</b><br/>Danger</p> <p><b>Hazard Statements</b><br/>H302 - Harmful if swallowed<br/>H411 - Toxic to aquatic life with long lasting effects<br/>H314 - Causes severe skin burns and eye damage<br/>H330 - Fatal if inhaled<br/>H350 - May cause cancer<br/>H225 - Highly flammable liquid and vapor<br/>Contains Zinc sulfate, Methyl alcohol, Isopropyl alcohol, Acetic acid</p> <p><b>Precautionary Statements - EU (§28, 1272/2008)</b><br/>P280 - Wear eye protection/ face protection<br/>P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing<br/>P310 - Immediately call a POISON CENTER or doctor/ physician<br/>P210 - Keep away from heat/sparks/open flames/hot surfaces. - No smoking<br/>P370 + P378 - In case of fire: Use dry sand, dry chemical or alcohol-resistant foam for extinction<br/>P304 + P340 - IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing<br/>P403 + P233 - Store in a well-ventilated place. Keep container tightly closed<br/>P260 - Do not breathe dust/fume/gas/mist/vapors/spray<br/>P303 + P361 + P353 - IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/shower<br/>P201 - Obtain special instructions before use<br/>P281 - Use personal protective equipment as required<br/>P308 + P313 - IF exposed or concerned: Get medical advice/ attention</p> |
|-------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## SHELF LIFE AND STORAGE

Shelf life of the Para-Pak Systems is indicated on the outer package label. Store at room temperature (15-30 C). Excessive heat or cold should be avoided.

## SPECIMEN COLLECTION AND PREPARATION

1. The patient should be cautioned against the use of antacids, barium, bismuth, antidiarrheal medication, or oily laxatives prior to collection of the specimen.
2. To assure recovery of parasitic elements which are passed intermittently and in fluctuating numbers, three specimens spaced a few days apart must be examined. In the case of hospitalized patients it is suggested that all fecal passages be collected for a designated length of time to avoid prolonging the hospital stay.<sup>5,9</sup>
3. The specimen is ideally passed into a bedpan but must not be contaminated with urine. Alternatively, plastic wrap may be placed in the toilet seat opening and the specimen passed into the bag. A thoroughly cleaned and dried milk carton, cut so as to remove the upper two thirds of the carton may also be used. It will be easier to collect the specimen if the water supply to the toilet is shut off and water drained from the bowl.
4. An appropriate (i.e. bloody, slimy, watery) area of stool should be selected and sampled with the collection spoons provided in the caps of the vials. Sufficient stool is added to each vial to bring the liquid level up to the "Fill to Here" line. This will result in approximately 5 mL of sample. To insure ideal sampling of a formed stool, material should be removed from the sides, ends, and middle of the bolus.
5. Agitate each specimen with the spoon along the sides of the vial, tighten the cap and shake firmly to insure that the specimen is adequately mixed. When mixing is completed the specimen should appear homogeneous.
6. Label each vial appropriately. Return the vials to the ziplock bag, seal the bag.

## TEST PROCEDURE

The Para-Pak Systems lend themselves to a wide variety of procedures in common use. The following discussion is not exhaustive and alternatives may be found in the literature cited. While variations exist from lab to lab, a thorough examination should include at least four steps:

1. Gross examination: record the presence of blood, worms, mucus, or proglottids.
2. Direct microscopic examination from the 10% formalin preserved specimen:
  - a. Place a clean glass slide on a sheet of newsprint.
  - b. Add a drop of saline (iodine may be substituted if desired) to the slide.
  - c. Add a representative sample of formalin preserved specimen to the drop of saline and mix thoroughly with the collecting spoon. The newsprint must be just legible through the slide.
  - d. Place a double width coverslip on the suspension and examine immediately.
3. Permanent slides for staining with Wheatley's Trichrome Stain (Meridian Cat. No. 400101), iron-hematoxylin, etc.:
  - a. Pour some of the Zn-PVA-fixed material onto a paper towel and allow to stand for approximately three minutes. This will absorb excess PVA and is considered critical to obtain the best possible staining.<sup>8</sup>  
**NOTE:** Do not allow the fixed material to absorb completely into the paper towel.
  - b. Using an applicator stick or brush, spread (avoid smearing) some of the specimen from the paper towel onto one or more clean glass slides. For best adherence, spread the material to the edge of the slide.
  - c. The slides are dried overnight at room temperature or for several hours in a 37 C incubator or slide warmer. Accelerated drying may cause some morphological distortion. The slides must be dried thoroughly to avoid washing the film off during staining.
  - d. Once slides prepared in Zn-PVA have dried completely, staining may begin by placing the slides directly into the trichrome stain.
  - e. Staining Procedure: \*

|                                                             |              |
|-------------------------------------------------------------|--------------|
| 1. Trichrome Stain                                          | 6-8 minutes  |
| 2. 90% Alcohol Acidified                                    | 2 quick dips |
| 3. 95% Alcohol                                              | Dip twice    |
| 4. 95% Alcohol                                              | 5 minutes    |
| 5. 100% Ethanol or Carbol-xylene                            | 10 minutes   |
| 6. Xylene                                                   | 10 minutes   |
| 7. Mount with coverslip using the mounting medium of choice |              |

  
\*Manual of Clinical Microbiology ASM publication. (**Note:** Other Trichrome staining procedures have shown to be compatible with Para-Pak)
4. Concentration procedures for 10% Formalin: One or more concentration procedures should be employed. No one concentration procedure works equally well for all parasites;<sup>1, 5, 11</sup> however, two in common usage that lend themselves well with to use with the Para-Pak/Para-Pak ULTRA systems are:
  - A. Formalin-ether (ethyl acetate)<sup>7, 14</sup> sedimentation:
    1. Mix the 10% formalin specimen thoroughly. The specimen is now ready for processing with the Para-Pak Macro-Con or CON-Trate Stool Concentration Systems. See the appropriate package insert for further directions.  
If Para-Pak Macro-Con or CON-Trate is not available, a sufficient quantity of specimen must be strained into a 15 mL conical centrifuge tube through one layer of narrow mesh or two layers of wide mesh gauze to provide 1.0 mL of sediment. This amount will vary with the size and density of the specimen.
    2. Add approximately 6-8 mL of 10% formalin (or saline), mix thoroughly and allow to stand five minutes.
    3. Add 3 mL of ethyl acetate or ether then stopper and shake vigorously for at least 30 seconds. Carefully remove the stopper.
    4. Centrifuge at 500 xg for 10 minutes (1800-2200 mm for most tabletop centrifuges).
    5. Four layers will be apparent:
      - a. Top layer: ethyl acetate or ether
      - b. Second layer: plug of debris
      - c. Third layer: formalin
      - d. Bottom layer sediment
    6. After ringing the plug of debris from the sides of the tube with an applicator stick, carefully decant the top three layers. While keeping the tube inverted, a cotton swab may be used to remove debris sticking to the sides of the tube. This is particularly important for obtaining suitable results with ethyl acetate and avoids solvent bubbles in the wet mount.
    7. Add a few drops of physiological saline or 10% formalin to resuspend the remaining sediment. If the resulting slides are too dense (newsprint should be legible through them) more saline or formalin may be added.
    8. Iodine and saline mounts are suggested for microscopic examination.  
**Note:** Rarely, a specimen may be extremely thick or mucoid, necessitating a preliminary wash. The wash, in physiological saline or tap water, may be implemented between steps 1e and 2. The washed material should be centrifuged as in step 4.
  - B. Zinc Sulfate flotation:
    1. Thoroughly mix a representative portion of the 10% formalin stool suspension or **fresh** unpreserved specimen in a 15 mL centrifuge tube and q.s. with tap water to approximately 10-12 mL. The amount of specimen to use will vary with its size and density.
    2. Centrifuge 1 minute at 1000-1200 xg.
    3. If sediment is about 1 milliliter in volume, decant the supernatant fluid. Otherwise adjust the density of the suspension by adding material from the 10% formalin suspension or diluting with more water. If adjustment of the sediment was necessary or if the stool is very oily, repeat the wash procedure.
    4. When using formalinized specimens, the specific gravity of the zinc sulfate solution must be adjusted to 1.2.<sup>1, 5</sup> Fill the tube about half full with zinc sulfate solution and resuspend the sediment by mixing thoroughly with applicator sticks.
    5. Add additional zinc sulfate solution to within 1 inch of the top.
    6. Centrifuge 1 minute at 1000-1200 xg.
    7. Carefully remove the tube from the centrifuge and, avoiding agitation, place it upright in a test tube rack or other suitable holder.
    8. Carefully fill the tube to the brim with zinc sulfate solution. Do not allow any overflow.
    9. A clean coverslip may now be placed on top of the tube. If the coverslip does not contact the meniscus of the liquid, carefully add more zinc sulfate solution until it does.
    10. Do not disturb the tube or coverslip for 10 minutes.
    11. With a quick, deft motion remove the coverslip straight upward so that a drop of liquid containing eggs and cysts adheres to the center of the coverslip.
    12. The coverslip may now be placed on a clean glass slide. If an iodine mount is desired, place a small drop of iodine on the slide prior to adding the coverslip. Sealing the edge of the slide with Vaspar (vaseline/paraffin 1:1 mixture) will prevent drying and distortion of the larger eggs.

#### QUALITY CONTROL

***This test should be performed per applicable local, state, or federal regulations or accrediting agencies.***

1. Visual inspection: Vials of Zn-PVA should contain approximately 15 mL of clear fluid, to insure 1:3 stool to preservative ratio.
2. If gelled, the fixative may be liquified by placing in a 50 C water bath until clear and fluid.
3. When a PVA fixed film of stock trophozoite or human buffy coat is stained, the organisms or cells should appear well fixed and defined.

**If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.**

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#### REFERENCES

1. Bartlett, Marilyn, et al. Comparative evaluation of a modified zinc sulfate flotation technique. J Clin Microbiol 1978;7:524-528.
2. Brook, MM and M. Goldman. Polyvinyl alcohol fixative as a preservative and adhesive for protozoa in dysenteric stools, and other liquid materials. J Lab and Clin Med 34: 1944;1554-1560.
3. Brooke, MM and Norman C. The effectiveness of the PVA fixative technique in revealing intestinal amebae in diagnostic cultures. Am J Trop Med Hyg 1955;4:479-482.
4. Brook, MM. PVA fixative technique in laboratory confirmation of amebiasis. Triangle 1960;4:326-335.
5. Brooke, MM. "Intestinal and Urogenital Protozoa" in Manual of Clinical Microbiology, ASM, Washington, D.C., 2<sup>nd</sup> ed., 1974;582-601.
6. Burrows, RB. *Microscopic Diagnosis of the Parasites of Man*, Yale University Press, New Haven, 1965.
7. Erdman, Dean 1981. Clinical comparison of ethyl acetate and diethyl ether in the formalin-ether sedimentation technique. J Clin Microbiol 14:483-485.
8. Garcia, LS. and L. Ash. *Diagnostic Parasitology*, C.V. Mosby Co., St. Louis, 1979;pp 9-24.
9. Harper, K. et al. Advantages of the PVA fixative two-bottle stool collection technique in the detection and Identification of intestinal parasites. Pub Health Lab 1957;15:96-108.
10. Horen, WP. Modification of Schaudin's fixative. J Clin Microbiol 1981;13:204-205.
11. Scholten, TH. and J. Yang. Evaluation of unpreserved and preserved stools for the detection and identification of intestinal parasites. Am J Clin Path 1974;62:563-567.
12. Swartzwalder, J. Clyde, GW. Hunter, WW. Frye. *Manual of Tropical Medicine*, W.B. Saunders Co., Philadelphia, 1966.
13. Simitch, T. and Z. Petrovich. Longevité de la forme vegetive de dysenterioe dans divers rñeaux. Arch Int Pasteur d'Algerie 1953;31:375-380.
14. Young, Kirk H. et al. Ethyl acetate as a substitute for diethyl ether in the formalin-ether sedimentation technique. J Clin Microbiol. 1979;10:852-853.

SN10715

REV. 11/25

|                                                                                                                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <br><b>Manufactured By</b>                | <p>Meridian Bioscience, Inc.<br/> 3471 River Hills Drive<br/> Cincinnati, OHIO - 45244 USA<br/> <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><u>Contacts:</u><br/> Main Telephone (+1) 513.271.3700<br/> Customer Service/Orders 800.543.1980<br/> Technical Support Center 800.343.3858<br/> Information Fax: 513.272.5432<br/> Ordering Fax: 513.271.0124<br/> E-mail: <a href="mailto:info@meridianbioscience.com">info@meridianbioscience.com</a></p>                                                                                                                                                                                                                                                                                                                                                     |
| <div data-bbox="469 577 654 632" style="border: 1px solid black; padding: 2px; display: inline-block;"> EC    REP </div>   | <p>Meridian Bioscience Europe, SRL<br/> Via Dell'Industria 7, 20035 Villa Cortese<br/> (Milano) ITALY<br/> <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><u>Contacts:</u><br/> Main Telephone (+39) 0331.433636<br/> E-mail: <a href="mailto:info@meridianbioscience.eu">info@meridianbioscience.eu</a><br/> Technical Support: <a href="mailto:MBE-TechService@meridianbioscience.eu">MBE-TechService@meridianbioscience.eu</a><br/> Customer Service/Orders:</p> <ul style="list-style-type: none"> <li>• For Italian Customers:<br/> <a href="mailto:ordini@meridianbioscience.com">ordini@meridianbioscience.com</a></li> <li>• For Distributors / International Customers:<br/> <a href="mailto:Export.CustomerService@meridianbioscience.eu">Export.CustomerService@meridianbioscience.eu</a></li> </ul> |
| <b>AUSTRALIAN<br/>SPONSOR</b>                                                                                              | <p>Emergo Australia<br/> Level 20, Tower II<br/> Darling Park<br/> 201 Sussex Street<br/> Sydney, NSW 2000<br/> Australia</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>UK Authorised<br/>Representative</b>                                                                                    | <p>Launch Diagnostics Limited<br/> Lakeview West<br/> Crossways Business Park<br/> Galleon Boulevard<br/> Dartford<br/> Kent DA2 6QE<br/> United Kingdom</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <div data-bbox="444 1163 738 1247" style="border: 1px solid black; padding: 2px; display: inline-block;"> CH    REP </div> | <p><b>MedEnvoy Switzerland</b><br/> Gotthardstrasse 28<br/> 6302 Zug<br/> Switzerland</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |

## INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.