# Para-Pak® CON-Trate® Stool Concentration Kit System

REF 960500

IVD

Rx Only

**INTENDED USE**  
The Meridian Para-Pak CON-Trate System is a complete system for concentrating and recovering helminth eggs, larvae and protozoan cysts from feces. Kit systems are designed for easy use by individuals not trained in microbiological procedures and afford an excellent means of minimizing the adverse effects of delay in specimen transportation.

**SUMMARY AND EXPLANATION OF THE TEST**  
Correct diagnosis of intestinal parasitic infection depends on proper collection, transport, detection, and identification of parasites in fecal specimens.<sup>1, 2, 4, 5, 7, 9</sup>

The detection and identification of small numbers of helminth eggs and protozoan cysts from large volumes of feces is of critical importance to the diagnosis of intestinal parasitic infection.

Historically, fecal specimens fixed in 10% formalin have been used in concentration procedures for detection and identification of parasitic elements.<sup>1, 2, 6, 9, 10, 12, 13</sup> Investigators have reported on the use of other fixatives such as SAF.<sup>12</sup> The Meridian Para-Pak CON-Trate system combines proven methodology with recent advances that optimize detection and identification, and minimizes specimen preparation in a safe, standardized, easily performed manner.<sup>8</sup>

**BIOLOGICAL PRINCIPLES**  
The Meridian Para-Pak CON-Trate System uses efficient, cost effective methods for recovering protozoan cysts, helminth eggs (including operculated eggs) and larvae, from preserved fecal specimens.

1. The addition of CON-Trate Reagent A and thorough mixing of the preserved specimen enhances the breakdown of fecal aggregates and mucus, thus freeing parasites.<sup>8</sup>
2. Filtering the stool-Reagent A suspension through the unique CON-Trate filtering device removes macroscopic fecal aggregates and debris.<sup>8</sup>

**REAGENTS/MATERIALS PROVIDED**  
*The maximum number of tests obtained from this test kit is listed on the outer box.*

**Catalog #960500**  
Reagent A (item #9601-03)  
CON-Trate tubes (item #3810A)  
Spin Con Funnel and filter (item #3814A)  
Caps for CON-Trate (item #3813)

**MATERIALS NOT PROVIDED**

1. Cotton-tipped applicator sticks
2. Microscope slides and coverslips
3. Physiological saline, 10% buffered formalin
4. Centrifuge
5. Microscope
6. Pipettes

**PRECAUTIONS**

1. All reagents are for in vitro diagnostic use only.
2. Observe good technique in handling and disposal of all biohazardous clinical and laboratory specimens and material.
3. Concentration of fecal specimens for parasitic examination is only an integral part of the overall scheme for identification of intestinal parasites. Additional tests and procedures may be found in appropriate references.
4. This product must not be used if:
  - a. The expiration date on the label has passed.
  - b. Proper storage conditions have not been observed.
5. **IMPORTANT:** See SDS for additional safety and hazard information.

Any serious incident that has occurred in relation to the device should be reported to Meridian Bioscience, Inc., 3471 River Hills Drive, Cincinnati, Ohio 45244 USA or Technical Support Center 800.343.3858 and competent authority of the EU Member State in which the clinician and/or patient is established.

## HAZARD AND PRECAUTIONARY STATEMENTS

Reagent A (MucoPenX), CON-Trate Reagent A  
and Macro-CON Surfactant

**Signal word**  
Warning  
**Hazard Statements**  
H302 - Harmful if swallowed  
Contains Polyethylene glycol 1,1,3,3-tetramethylbutylphenyl ether

**SHELF LIFE AND STORAGE**  
Shelf life of the Para-Pak CON-Trate system is indicated on the outer package label. Store at 15-30 C. Do not freeze.

## SPECIMEN COLLECTION AND PREPARATION

### A. SPECIMEN COLLECTION

1. The patient must be properly instructed and assisted in the collection of the fecal sample. Refer to Meridian Para-Pak collection/transport kit product inserts or appropriate references for recommended collection and transport methods.
2. Any specimen preserved in 10% formalin, SAF, or an unpreserved fecal sample may be used with the CON-Trate System.

### B. SPECIMEN QUALITY ASSURANCE

To assure a suitable clinical specimen using the CON-Trate System, observe the following:

1. The specimen/preservative mixture must be stored for a minimum of 30 minutes after collection for adequate fixation. **IMPORTANT:** Mix contents thoroughly.
2. The specimen must have been maintained at room temperature.
3. Adequate sample must be present. (2-3) gm of stool in 15 mL fixative is recommended.

### C. SPECIMEN PROCESSING

1. **Unpreserved Specimens** (For optimum results, it is recommended that specimens be preserved at the time of collection. Unpreserved specimens delayed in transport may have limited diagnostic value.<sup>1, 2, 7, 12</sup>)
  - a. Transfer 3-5 grams of unpreserved stool into 15 mL preservative of choice, for economy and convenience, use Meridian Para-Pak 10% Buffered Neutral Formalin (Catalog #900412) or SAF (Catalog #900212). Mix stool/preservative mixture thoroughly, break up any lumps or fecal masses (shaking for one minute is usually sufficient). The stool/preservative mixture should stand for a minimum of 30 minutes for adequate fixation. Follow procedure far use below for the preservative selected.
  2. **Immediate processing of unpreserved specimens**<sup>9</sup>
    - a. Transfer 5-6 grams of unpreserved stool into 10-15 mL of physiological saline. Mix stool/saline mixture thoroughly, break up any lumps or fecal masses.
    - b. Add 4 drops of CON-Trate Reagent A to the mixture. (Up to 8 drops of Reagent A may be added if the specimen is highly mucoid.)
    - c. Cap and mix the contents thoroughly by shaking.
    - d. Insert one of the CON-Trate filtering devices into the top of one of the disposable centrifuge tubes provided.
    - e. Pour fecal suspension through the filtering device into the centrifuge tube. Usually 3 mL in the centrifuge tube is sufficient unless the fecal suspension is thin.
    - f. Discard filtering device, add 10 mL physiological saline and centrifuge at 500 xg for 10 minutes (1800-2200 rpm). Decant the supernatant fluids, retaining the sediment. About 1 mL of sediment should be present. A portion of the sediment may be used for detection of *Cryptosporidium*. Consult an appropriate reference for proper preparation and examination. (See Hint 6)
    - g. Resuspend the sediment in 10 mL buffered formalin. Allow mixture to stand for at least 5 minutes before proceeding. (At this point, the mixture in the centrifuge tube may be capped and saved until a later time.)
    - h. Centrifuge the tube at 500 xg for 10 minutes (1800-2200 rpm). Examination of the tube after centrifugation should reveal four distinct layers from the top down:
      1. a "plug" of fecal debris
      2. a discolored aqueous layer
      3. a sediment layer, containing the parasites.The final sediment remaining should be approximately 0.25 mL.
    - i. Hold the tube in a vertical position. Free the plug of debris by ringing with a wooden applicator stick. Decant the upper layers, leaving the sediment. **DO NOT TURN THE TUBE UPRIGHT UNTIL THE SIDES OF THE TUBE HAVE BEEN CLEANED WITH COTTON TIPPED SWABS.**
    - j. Transfer a portion of the sediment to a clean glass microscope slide and prepare the mount of choice. Examine microscopically. Consult an appropriate reference for proper preparation and examination.<sup>5, 7, 9</sup>
  3. **Formalin preserved specimens**
    - a. Add 4 drops of CON-Trate Reagent A to the specimen in the fixative vial. (Up to 8 drops of Reagent A may be added if the specimen is highly mucoid)
    - b. Cap and mix the contents thoroughly shaking.
    - c. Insert one of the CON-Trate Filtering devices into the top of one of the disposable centrifuge tubes provided.
    - d. Pour fecal suspension through the filtering device into the centrifuge tube. Usually 3 mL in the centrifuge tube is sufficient unless the fecal suspension is thin.
    - e. Discard filtering device, add 10 mL saline and centrifuge at 500 xg for 10 minutes (1800-2200 rpm). Decant the supernatant fluids, retaining the sediment. A portion of the sediment may be used for detection of *Cryptosporidium*. Consult an appropriate reference for proper preparation and examination. (See Hint 6)
    - f. Resuspend the sediment in 9 mL of 10% formalin.
    - g. Centrifuge the tube at 500 xg for 10 minutes (1800-2200 rpm). Examination of the tube after centrifugation should reveal four distinct layers from the top down:
      1. a "plug" of fecal debris
      2. a discolored aqueous layer
      3. a sediment layer, containing the parasites.The final sediment remaining should be approximately 0.25 mL.
    - h. Hold the tube in a vertical position. Free the plug of debris by ringing with a wooden applicator stick. Decant the upper layers, leaving the sediment. **DO NOT TURN THE TUBE UPRIGHT UNTIL THE SIDES OF THE TUBE HAVE BEEN CLEANED WITH COTTON TIPPED SWABS.**
    - i. Transfer a portion of the sediment to a clean glass microscope slide and prepare the mount of choice. Examine microscopically. Consult an appropriate reference for proper preparation and examination.<sup>5, 7, 9</sup>
  4. **SAF preserved specimens**
    - a. Add 4 drops of CON-Trate Reagent A to the specimen in the fixative vial (Up to 8 drops of Reagent A may be added if the specimen is highly mucoid).
    - b. Cap and mix the contents thoroughly by shaking the vial several times.
    - c. Insert one of the CON-Trate filtering devices into the top of one of the disposable centrifuge tubes provided.
    - d. Pour fecal suspension through the filtering device into the centrifuge tube. Usually 3 mL in the centrifuge tube is sufficient unless the fecal suspension is thin.
    - e. Discard filtering device, add 10 mL saline and centrifuge at 500 xg for 10 minutes (1800-2200 rpm). Decant the supernatant fluids retaining the sediment. A portion of the sediment may be used for detection of *Cryptosporidium* species. Consult an appropriate reference for proper preparation and examination. (See Hint 6)
    - f. If permanent stained smears are to be prepared:
      1. Transfer some of the sediment to a drop of Mayer's Albumin on a slide.
      2. Mix well, then spread mixture over a clean, glass microscope slide producing an uneven film (thick and thin areas). This smear should be allowed to dry and can then be stained with a permanent stain such as Iron Hematoxylin or Wheatley's Trichrome (Catalog #400101).<sup>2</sup>
    - g. Resuspend the sediment in 9 mL of 10% formalin.
    - h. Centrifuge the tube at 500 xg for 10 minutes (1800-2200 rpm). Examination of the tube after centrifugation should reveal four distinct layers from the top down:
      1. a "plug" of fecal debris
      2. a discolored aqueous layer
      3. a sediment layer, containing the parasites. The final sediment remaining should be approximately 0.25 mL.
    - i. Hold the tube in a vertical position. Free the plug of debris by ringing with a wooden applicator stick. Decant the upper layers, leaving the sediment. **DO NOT TURN THE TUBE UPRIGHT UNTIL THE SIDES OF THE TUBE HAVE BEEN CLEANED WITH COTTON TIPPED SWABS.**
    - j. Transfer a portion of the sediment to a clean glass microscope slide and prepare the mount of choice. Examine microscopically. Consult an appropriate reference for proper preparation and examination.<sup>5, 7, 9</sup>

### SPECIFIC PERFORMANCE CHARACTERISTICS

The Meridian Para-Pak CON-Trate System, when used as directed under procedure for use, during clinical evaluation was found to provide results comparable/equivalent/superior<sup>8</sup> to the standard Ritchie-Formalin ether (Ethyl Acetate) sedimentation concentration procedure.

MucoPenX (Reagent A) has been formulated to break down the mucus present in stool specimens and will not interfere with staining procedures or cause distortion of parasites.

### HINTS

Experience will dictate appropriate techniques and volumes to assure an adequate sediment for microscopic examination. The following is a list of hints and suggestions which will accomplish this objective.

1. **Volume of preserved fecal suspension to add through the filter device:** After the first centrifugation, a sediment of 1.0 mL is optimum. In a dense fecal suspension, with a ratio of stool to preservative of 1:3-1:5, 3.0 mL of filtrate will provide the optimum sediment volume. In less dense fecal suspensions, such as one would see in a watery stool, larger volumes of filtrate are necessary. In some instances, as much as 10.0 - 12.0 mL of filtrate may be necessary. The hint is that the **decreasing** density of the preserved, fecal suspension dictates a proportional **increase** in the filtrate volume. If sufficient fecal material was present in the original fecal suspension and care was taken in adjusting the filtrate volume, the final sediment volume, after concentration, optimally, should be 0.25 mL. Attention to detail and experience will provide the anticipated results.
2. **Fecal suspension should not be forced through the filter device:** With very dense fecal suspensions, the flow rate through the filter device is slow. Do not force the material through the filter device by scraping with applicator sticks or washing the aqueous solutions. Such action often forces hard, grit-like material through the filter device. The presence of these materials in wet mount preparations results in a non-uniform flow or distribution of materials under the coverslip, making coverslipping difficult. The problem is minimized by adding the recommended amount of Reagent A and thoroughly mixing the specimen/fixative mixture. Vegetative matter may cover the screen in the filtration device; this may be removed by gently running an applicator stick through the material, or by tapping the side of the filtration device.
3. **Swabbing centrifuge tube:** Failure to swab the sides of the centrifuge tube after decanting the fecal debris plug and aqueous layer can result in a poor wet mount preparation.
4. **Sediment:** When the CON-Trate procedure is followed correctly, the sediment will appear dry and gritty. To facilitate reading, a drop of saline should be added to the sediment on the slide and the coverslip floated on top.
5. **Wet Preparation Theory:** Tradition has dictated that a direct microscopic examination be performed on stool specimens submitted to the laboratory. Tradition has perhaps obscured the rationale for this requisite. The direct microscopic examination of fresh, unpreserved, stool specimens or bowel materials facilitated the identification of protozoa by noting the characteristic motility of protozoan trophozoites. The common practice of submitting preserved stool specimens to the laboratory may eliminate the need for the direct microscopic examination. The detection rate of parasites may be increased by performing the first microscopic examination on the first centrifuged sediment, or by staining the sediment from an SAF preserved specimen with Wheatley's Trichrome or Iron Hematoxylin.
6. **Cryptosporidium Prep Using Para-Pak:**  
**For loose or watery stools:**  
After filtration of specimen, centrifuge sample at 500 xg for a full 10 minutes (1800-2200 rpm for most tabletop centrifuges). Decant the supernatant fluid. Approximately 0.5 to 1.0 mL of sediment should remain. Mix the sediment and prepare a smear for *Cryptosporidium*.  
**For semi-solid or solid stools:**  
Perform entire concentration procedure and prepare a slide from the sediment per lab protocol.

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**Catalogue n 960500**

Reagent A (article #9601-03)  
CON-Trate tubes (article #3810A)  
Spin-Con Funnel and filter (article #3814A)  
Caps for CON-Trate (article #3813)

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**REFERENCES**

1. ASMT. Recommended procedures for the examination of clinical specimens submitted for the diagnosis of parasitic infections. Am J Med Technol 1978;44:1101-1106.
2. American Society of Parasitologists. Procedures suggested for use in examination of clinical specimens for parasitic infection. J Parasitol 1977;63:959-960.
3. Blagg W et.al. A new concentration technique for the demonstration of protozoa and helminth eggs in feces. Am. Journ. Trop. Med. and Hyg. 4:23-28. 1955
4. Brook MM. Intestinal and urogenital protozoa. Manual of Clinical Microbiology. ASM Washington DC, 2nd ed. 582-601. 1974
5. Burrows RB. Microscopic diagnosis of the parasites of men. New Haven: Yale University Press;1905
6. Garcia LS, Shimizu R. Comparison of clinical results for the use of ethyl acetate and diethyl ether M the formalin-ether sedimentation technique performed on polyvinyl alcohol preserved specimens. J Clin Microbiol 1981;13:709-713.
7. Garcia LS, Voge M. Diagnostic clinical parasitology: I. proper specimen collection and processing. Am J Med Technol 1980;46:459-467.
8. Long EG, Tsift ff AT, Robinson A. Comparison of the FeKal CON-Trate® System with the formalin-ethyl acetate technique for detection of intestinal parasites. J Clin Microbiol 1985;22:210-211.
9. Melvin DM, Brooke MM. Laboratory procedures for the diagnosis of intestinal parasites. Atlanta, GA: US DHEW CDC 1980;80:8282
10. Ritchie LS. An ether sedimentation technique for routine stool examination. Bull. U.S. Army Med. Dept. 8:326.
11. Sapero, James J, Lawless DK C.P.A. Strome. 1951. An Improved Iodine - staining Technique for Routine Laboratory Diagnosis of Intestinal Protozoa. Science. 1948.114:550-551.
12. Scholten JL Whitby. Comparison of ethyl-acetate and DiEthyl Ether in the Formalin Ether Concentration Method Based on Proficiency Test Samples. Present ed at the 54th Canadian Association for Clinical Microbiology and Infectious Disease, Quebec, Canada. 1986.
13. Yang J and Scholten TH. A fixative for intestinal parasites permitting the use of concentration and permanent staining procedures. Am J Clin Pathol 1977;67-300-304.
14. Data On File

**SN10652**


**REV. 11/22**

 <p style="text-align: center;"><b>Manufactured By</b></p>	<p><b>Meridian Bioscience, Inc.</b>            3471 River Hills Drive            Cincinnati, OHIO - 45244 USA  <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><u>Contacts:</u>            Main Telephone (+1) 513.271.3700            Customer Service/Orders 800.543.1980            Technical Support Center 800.343.3858            Information Fax: 513.272.5432            Ordering Fax: 513.271.0124            E-mail: <a href="mailto:info@meridianbioscience.com">info@meridianbioscience.com</a></p>
<p style="text-align: center;"><b>AUSTRALIAN SPONSOR</b></p>	<p><b>Emergo Australia</b>            Level 20, Tower II            Darling Park            201 Sussex Street            Sydney, NSW 2000            Australia</p>

**INTERNATIONAL SYMBOL USAGE**

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.