Para-Pak[®] Macro-CON[®]

Stool Concentration System

REF 970120

IVD

Rx Only

INTENDED USE

Para-Pak Macro-CON is a system for the concentration of eggs, larvae, and protozoa from preserved fecal specimens. The filtration unit is designed to be used directly with the Para-Pak specimen collection vial. This creates a **completely closed system**, minimizing exposure of the user to potentially infectious agents. **Macro-CON** concentration utilizes the entire contents of the Para-Pak specimen collection vial, reducing variation due to specimen sampling. This is of particular advantage to the clinical parasitologist when a small number of organisms may be present in a large volume of stool. Kit systems are designed for easy use by individuals not trained in microbiologic procedures.

SUMMARY AND EXPLANATION OF THE TEST

Diagnosis of intestinal parasitic disease is confirmed by the recovery and identification of helminth eggs and larvae, or protozoan trophozoites and cysts. Concentration procedures facilitate this process and are of importance when small numbers of protozoa are present in large volumes of feces.^{1, 2, 4, 7, 8, 11, 13-15}

The Para-Pak Macro-CON System (patent # 4,675,110 June, 1987) exhibits several advantages over both the traditional gauze method and other available concentration devices. Fastening the concentrator directly to the patient specimen vial minimized user contamination by creating a closed system. The contents of the vial need never be directly manipulated by laboratory personnel. Concentration directly from the Para-Pak vial minimizes hands-on manipulation, facilitating rapid laboratory diagnosis.

Proper use of the Macro-CON filtration unit in conjunction with the Para-Pak specimen collection vial assures the parasitologist that diagnostic stages of intestinal parasites, if present in fecal material, will be concentrated.

BIOLOGICAL PRINCIPLES

The filtration unit consists of a 50 mL conical centrifuge tube capped with a flanged adapter incorporating a plastic mesh screen.

Preserved stool is first thoroughly mixed with a surfactant which helps reduce adhesive forces and break down fecal aggregates, thus freeing parasites.^{10, 14}

Filtering the specimen through the Macro-CON filter removes large clumps and debris. This is facilitated by the unique venting process which allows for air exchange between the tube and vial as the extracted material filters into the centrifuge tube. The simple design of Macro-CON insures trouble free performance and less clogging.

Emulsification of the sample with ethyl acetate extracts much of the fecal lipids.^{8, 15}

After centrifugation the lipid plug and excess ethyl acetate are decanted leaving a packed sediment enriched with parasite eggs, larvae, trophozoites and cysts.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. Conical centrifuge tubes (item# CS14375)
2. Macro-CON filtration units (item# CS4230)
3. Screw caps (item# CS2210A)
4. Bottle surfactant (item# 9601-02)

MATERIALS NOT PROVIDED

1. Ethyl acetate (suggested) or diethyl ether (optional)
2. 10% formalin, or physiological saline
3. Cotton tipped applicator sticks
4. Microscope slides and coverslips
5. Centrifuge
6. Microscope
7. Capillary pipette
8. Sheather's sugar solution (optional)

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Ethyl acetate and diethyl ether are flammable. Use in a well ventilated area. Avoid open flame. Avoid contact of the solution with skin and eyes. Should contact occur, flush with running water. Avoid breathing fumes.
3. PVA fixative is corrosive for metals. Avoid contact with metals.
4. Due to the infectious nature of unpreserved stools, use of gloves, care and handwashing should be employed when the specimen is obtained and manipulated.
5. Para-Pak Macro-CON is designed specifically for use with the Para-Pak line for parasitology products. Proper product performance, especially a leak free fit to the patient specimen can only be assured when using the Para-Pak product line.
6. The Para-Pak Macro-CON is designed for use with the provided surfactant. Proper product performance can only be assured when following the concentration procedure outlined in the package insert.
7. Any serious incident that has occurred in relation to the device should be reported to Meridian Bioscience, Inc., 3471 River Hills Drive, Cincinnati, Ohio 45244 USA or Technical Support Center 800-343-3858 and competent authority of the EU Member State in which the clinician and/or patient is established.
8. **IMPORTANT:** See SDS for additional safety and hazard information.

HAZARD and PRECAUTIONARY STATEMENTS

 Reagent A (MucoPenX), CON-Trate Reagent A and Macro-CON Surfactant	Signal word Warning Hazard Statements H302 - Harmful if swallowed Contains Polyethylene glycol 1,1,3,3-tetramethylbutylphenyl ether
---	---

SHELF LIFE AND STORAGE

The Macro-CON filtration components are stable for three years from date of manufacture, when stored at room temperature.

TEST PROCEDURE

Any specimen preserved in 10% formalin, sodium acetate acetic acid formalin (SAF), or polyvinyl alcohol (PVA) may be used. A procedure for the use of fresh stool is also described.^{5, 8, 13, 14}

The methodology is a modification of the Ritchie formalin-ether procedure.¹⁴


A. SPECIMEN COLLECTION

1. An appropriate specimen is of utmost importance in the diagnostic process.^{1-4, 8, 11} Para-Pak Macro-CON is designed specifically for use with the Para-Pak line of parasitology products (See Precautions #5). Collect the specimen as directed by the appropriate Para-Pak package insert supplied with the collection system.
2. The specimen must be well mixed in the preservative. The preservative-stool mixture must be allowed to stand for a minimum of 30 minutes at room temperature between the time of specimen collection and processing to ensure adequate fixation.
3. If a permanently stained smear is desired from Para-Pak SAF or LV-PVA preserved material, it is recommended that a portion of preserved stool be removed prior to concentration. The procedure for permanent staining can be found in the Para-Pak SAF or LV-PVA package insert.

B. SPECIMEN PREPARATION

1. Remove the cap from the patient specimen vial and add 10 drops of surfactant to the specimen vial. Replace the cap securely.
2. Shake vigorously for 60 seconds. Breakup of some formed or very mucoid specimens may be facilitated by vortexing.
3. Check that the 50 mL conical centrifuge tube is securely screwed onto the filtration unit. Remove the cap from the patient specimen vial and insert the open end of the Macro-CON filtration unit into the specimen vial, as shown in figure 1, with a gentle downward pressure until firmly seated.


Figure 1

Insert filtration unit tightly into top of Para-Pak vial.

4. Invert the Macro-CON system, and allow the specimen to filter through the mesh into the 50 mL conical centrifuge tube, as shown in Figure 2. If the specimen is thick, the centrifuge tube may be tapped sharply on a table top to facilitate the flow of material through the mesh.

Figure 2

Tap sharply to force solution into conical tube.

5. Grasping the serrated collar, loosen BUT DO NOT REMOVE the filtration unit from the conical centrifuge tube. Allow the Macro-CON system to sit undisturbed for 60 seconds. This will release pressure in the system and effect better drainage of the filter.
6. Re-tighten the filtration unit onto the conical centrifuge tube and tap sharply two or three times to force any remaining liquid into the conical centrifuge tube.
7. Tilt the Macro-CON system to approximately 30 degrees from horizontal and unscrew the filtration unit completely from the conical centrifuge tube.
8. Discard the filtration unit/Para-Pak specimen vial using the procedures established in your laboratory for disposal of fecal specimens.
9. Add 10% formalin to bring the level of the filtered specimen to the dotted line on the Para-Pak Macro-CON label.
10. Add 5 mL of ethyl acetate (diethyl ether, optional) to the conical centrifuge tube. Place a provided screw cap on the centrifuge tube securely. (Note: do not use the screw cap from the original specimen container of leakage and contamination may result.)
11. Shake vigorously for 60 seconds. If using diethyl ether, loosen the cap on the conical centrifuge tube to release pressure prior to centrifugation.
12. Retighten the cap on the conical centrifuge tube and centrifuge for 10 minutes at 500 xg (1800-2200 rpm for most table top centrifuges).
13. Pour off the supernatant fluid and debris layer and **KEEP THE TUBE INVERTED**.
14. **KEEPING THE TUBE INVERTED**, clean the inner walls with one or two cotton tipped applicator sticks to remove any remaining debris and ethyl acetate (this may be omitted if using diethyl ether). **Allowing excess ethyl acetate to run back into the pelleted sediment will result in a poor wet mount preparation due to the formation of solvent bubbles.** The packed sediment at the bottom of the tube contains the parasitic elements.
15. Add a few drops of 10% formalin or physiological saline and mix the pelleted sediment. If a permanently stained slide is desired from Para-Pak LV-PVA concentrated material, the pellet should be resuspended in a few drops of LV-PVA. Allow sediment to settle for a minute.
16. The specimen for examination should be drawn with a 5-9 inch capillary pipette from the top half of the resuspended material. This portion contains the greatest number of parasites and avoids larger, denser particles which may interfere with wet mount preparation. Place 1-2 drops from the capillary pipette on a microscope slide, place a coverslip on the suspension and examine immediately. When a wet mount is made with the resuspended pellet, newspaper should be legible through the prepared slide.
17. If iodine staining is desired, place a single drop of Lugol's iodine from a capillary pipette on the slide, and add a drop of resuspended sediment. Place a coverslip on the suspension and examine immediately.
18. The remaining sediment may be used for *Cryptosporidium* and *Isospora belli* examination. Thin smears may be made for acid fast staining. Sheather's sugar flotation may also be performed.^{6, 9, 12, 16}

NOTE: Rarely, a specimen may be extremely thick or mucoid necessitating a preliminary wash. The wash, in physiological saline or tap water, may be implemented between steps 8 and 9. The washed material should be centrifuged as in step 12.

C. PROCEDURE FOR FRESH STOOL SPECIMENS

This procedure permits the concentration technique to be performed on a Para-Pak Clean vial specimen after amoeba culture, rearing of hook worm larvae, bacterial culture, occult blood, or stool fat examination.

3-5 g of stool is suitable for the addition of 15 mL of 10% formalin.

1. Add 15 mL of 10% formalin to the 3-5 g stool specimen in a Para-Pak Clean vial.
2. Mix the contents thoroughly.
3. Allow to stand at room temperature for 30 minutes to ensure adequate fixation.
4. Follow the procedure in Section B, Specimen Processing, Steps 1-18.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

Macro-CON filtration units should be free of any cracks. The mesh screen should be uniform in appearance.

The surfactant should be free of any bacterial or fungal contamination.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

SPECIFIC PERFORMANCE CHARACTERISTICS

Clinical Data

Performance of the Para-Pak Macro-CON concentrator was evaluated by comparison to a common laboratory procedure in two separate locations. LV-PVA or 10% Formalin preserved specimens were concentrated using the Macro-CON system (MC) and a common laboratory procedure (CLP). After identifying the parasite or parasites present, each species was enumerated in each system.

The data are expressed as the number of times each organism was isolated. A range of counts, for each system, is also included.

In six cases only one system isolated an organism. In each case < 2 organisms per total coverslip were enumerated in the positive system.

ASSAY SPECIFICITY

IDENTIFICATION

	MC	RANGE	CLP	RANGE
<i>Ascaris lumbricoides</i>	24	16-390/CS	24	11-582/CS
<i>Giardia lamblia</i>	22	1/LPF-10+/HPF	22	.5/LPF-3/HPF
Hookworm	24	1-3.6 6/LPF	23	0-3.3/LPF
<i>Trichuris trichura</i>	20	0-100+/CS	23	1-100+/CS
<i>Entamoeba coli</i>	14	7/CS-25+/LPF	14	7/CS-20 LPF
<i>Endolimax nana</i>	9	.3-.9/HPF	9	2-5.2/HPF
<i>Entamoeba histolytica</i>	9	.6-13/HPF	9	.6-9/HPF
<i>Clinorchis sinensis</i>	6	4-50/CS	6	5-89/CS
<i>Strongyloides stercoralis</i>	5	0-15/CS	7	1-47/CS
<i>Hymenolepis nana</i>	7	15/CS-6.1/LPF	7	15/CS-19.8/LPF
<i>Taenia saginata</i>	4	1-248/CS	4	3-154/CS
<i>Chilomastix mesnili</i>	4	.9-11/HPF	4	1.2-20/HPF
<i>Cryptosporidium</i> sp.	1	3-HPF	1	3/HPF
<i>Entamoeba hartmanni</i>	6	.3-5/HPF	6	.3-4.5/HPF
<i>Opisthorchis viverrini</i>	2	31-147/CS	2	25-186/CS
<i>Diphyllobothrium latum</i>	3	46-800+/CS	3	100-800+/CS
<i>Isospora belli</i>	1	15/HPF	1	15/HPF
<i>Entamoeba polecki</i>	1	3.1/LPF	1	2.3/LPF

High Powered Field (HPF)

Low Powered Field (LPF)

Coverslip (CS)

In a Midwestern clinical laboratory, diarrheic stool specimens submitted for procedures other than ova and parasite examination, were concentrated using the methods previously described. Of 25 specimens examined, all were negative using MC and CLP.

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REFERENCES

1. ASMT. Recommended procedures for the examination of clinical specimens submitted for the diagnosis of parasitic infections. Am J Med Technol 1978;44:1101-1106.
2. American Society of Parasitologist. Procedures suggested for use in examination of clinical specimens to parasitic infection. J Parasitol 1977;63:959-960.
3. Brooke MM. Intestinal and urogenital protozoa, Manual of clinical microbiology 2nd ed. Washington DC: ASM; 1974 582-601.
4. Burrows RB. Microscopic diagnosis of the parasites of man. Yale University Pres, New haven, 1965 p. 319-381.
5. Garcia LS and Ash LR. Diagnostic parasitology. 2nd ed. The CV Mosby Co; 1979. p.1-18.
6. Garcia LS, Brewer TC, Bruckner DA and Schimizu RY. Clinical laboratory diagnosis of cryptosporidium from human fecal specimens. Clin Microbiol Newsl 1982;4:136-137.
7. Garcia LS and Shimizu R. Comparison of clinical results for the use of ethyl acetate and diethyl ether in the formalin-ether sedimentation technique performed on polyvinyl alcohol preserved specimens. J Clin Microbiol 1981;13:709-713.
8. Garcia LS and Voge M. Diagnostic clinical parasitology: I proper specimen collection and processing. Am J Med Technol 1980;46:459-467.
9. Ma P and Soave R. Three-step stool examination for cryptosporidiosis in 10 homosexual men with pretracted watery diarrhea. J Infect Dis 1983;147:824-828.
10. Maldonado JF, Acosta-Matienzo J and Velez-Herrera F. Comparative value of fecal examination procedures in the diagnosis of helminth infections. Exptl Parasitol 1953;2:403-416.
11. Melvin DM and Brooke MM. 1. Laboratory procedures for the diagnosis of intestinal parasites, US DHEW. 98079:8282, CDC Atlanta, GA, 23-65.
12. Ng E, Markell EK, Fleming RL, and Fried M. Demonstration of isospora belli by acid-fast stain in a patient with acquired immunodeficiency syndrome. J Clin Microbiol 1984;20:384-386.
13. Ritchie LS. An ether sedimentation technique for routine stool examination. Bull. US Army Med Dept 1948;8:326.
14. Rohrbacher GH. The recovery of nematode larvae by baermann apparatus as affected by a detergent. Proceedings of the helminth. Soc 1957;24:24-25.
15. Yang J and Scholten TH. A fixative for intestinal parasites permitting the use of concentration and permanent staining procedures. Am J Clin Pathol 1977;67:300-304.
16. Young KH, Bullock S, Melvin DM and Sprull CL. 1. Ethyl Acetate as a substitute for diethyl ether in the ether-formalin sedimentation technique. J Clin Microbiol 1979;10:852-853.
17. Zierdt WS.. Concentration and identification of cryptosporidium sp by use of a parasite concentrator. J Clin Microbiol 1984;20:860-861.

SN10600

REV. 04/24

 Manufactured By	Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com
EC REP	Meridian Bioscience Europe, SRL Via Dell'Industria 7, 20035 Villa Cortese (Milano) ITALY www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+39) 0331.433636 E-mail: info@meridianbioscience.eu Technical Support: MBE-TechService@meridianbioscience.eu Customer Service/Orders: - For Italian Customers: ordini@meridianbioscience.com - For Distributors / International Customers: Export.CustomerService@meridianbioscience.eu
UK Authorised Representative	Launch Diagnostics Ash House Ash Road Longfield DA3 8JD UK
CH REP	MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.