

immunocard
STAT![®]
Crypto/Giardia



Immunocard STAT![®] Crypto/Giardia

Rapid parasite detection when every minute matters

Immunocard STAT! Crypto/Giardia is a rapid immunoassay that detects *Cryptosporidium* and *Giardia* antigens in stool specimens. Reliable dual detection delivers clear results to support prompt diagnosis of parasitic infections. Early identification enables appropriate treatment decisions, helps reduce transmission risk, and supports efficient laboratory workflows.

meridian BIOSCIENCE[®]

Immunocard STAT! Crypto/Giardia

The Challenge

- *Cryptosporidium* and *Giardia* are among the leading causes of parasitic disease¹
- Microscopy is labor-intensive and operator-dependent
- Intermittent shedding can complicate detection
- Delayed identification may prolong symptoms and transmission

The Need

- Reliable dual-target detection in a single test
- Rapid turnaround to support clinical decisions
- Reduced subjectivity compared to traditional microscopy
- Efficient workflow for busy laboratory environments

Immunocard STAT! Crypto/Giardia Delivers

- Simultaneously detect *Cryptosporidium* and *Giardia* antigens in stool
- Rapid results support prompt diagnosis and treatment
- Simple procedure designed for routine laboratory use
- Clear visual readout enhances confidence in results

Specifications

CLIA Status

Moderate

Turnaround Time

10 minutes

Sample Type

Fresh, unpreserved, and preserved stool specimens

Sample Storage

Frozen (-20 C)

Refrigerated (2-8 C)

Room temperature (20-30 C)

Kit Storage

2-8 C

Performance

Cryptosporidium: Sensitivity 97.3%, Specificity 100%

Giardia: Sensitivity 100%, Specificity 100%

Catalog Number

750830

CPT Code

Cryptosporidium: 87328

Giardia: 87329



For more information on the Immunocard STAT! products, contact a specialist at meridianbioscience.com/contactus

References:

1. Centers for Disease Control and Prevention (CDC). Parasites – Causes. Available at: <https://www.cdc.gov/parasites/causes/>. Accessed March 17, 2026