

breathID[®]

IDkit Hp[®] Two



IDkit Hp[®] Two

A non-invasive urea breath test for the detection of *H. pylori* and confirmation of eradication.

Built for reliability and ease of use, the IDkit Hp Two simplifies *H. pylori* testing through fast, accurate sample collection and analysis with BreathID[®] systems. Its non-invasive testing method enhances patient comfort while delivering dependable results for both initial diagnosis and post-treatment confirmation of eradication. A one-way valve securely seals samples to avoid invalid results and maintain reliable diagnoses.

meridian BIOSCIENCE[®]

IDkit Hp Two

The Challenge

- *H. pylori* infection often presents with symptoms that mimic other gastrointestinal conditions, delaying accurate diagnosis
- Invasive procedures like endoscopy are time-consuming, uncomfortable, and unlikely to be used to confirm eradication
- If left undiagnosed and untreated, *H. pylori* infection can lead to chronic gastritis, peptic ulcer disease, or even gastric cancer¹

The Need

- A non-invasive test that reliably identifies active *H. pylori* infection and confirms successful eradication
- High-accuracy, guideline-consistent diagnostics to support efficient decision-making in GI clinics and reference labs
- Scalable, workflow-friendly solutions that meet the needs of both decentralized point-of-care testing and high-volume laboratory environments

BreathID Delivers

- Rapid, reliable results for active *H. pylori* infection and eradication confirmation
- Flexible testing options designed to meet the needs of any sized lab
- Fully automated analysis that minimizes operator steps and reduces turnaround time
- Non-invasive, patient-friendly testing backed by Meridian's trusted gastrointestinal diagnostic expertise

Specifications

CLIA Status

CLIA Non-regulated

Turnaround Time

2 minutes per test

Sample Storage

Room temp

Up to 14 day storage time

Kit Storage

Room temp

36 month shelf life

Performance

100% Sensitivity

97.9% Specificity

Catalog Number

IDkit Hp Two AC00063

CPT Code

83013 (analysis)

83014 (collection)



For more information on the BreathID products, contact a specialist at meridianbioscience.com/contactus

Reference:

1. Fennerty, M. B. *Helicobacter pylori*: why it still matters in 2005. *Cleveland Clinic Journal of Medicine*, 2005; 72(Suppl_2), S1-7. https://doi.org/10.3949/ccjm.72.suppl_2.s1.