

## Normal Mouse Serum, Heat Inactivated

Cat. No. 88131

STORE AT -30 to -10°C

<b>DESCRIPTION</b>	Normal Mouse Serum, Heat Inactivated
<b>TYPE</b>	Passive Blocker
<b>SHIPPING &amp; STORAGE INFORMATION</b>	Aliquot and store at -30 to -10°C. Avoid multiple freeze/thaw cycles.
<b>SAFETY PRECAUTIONS</b>	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
<b>NOTES</b>	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

### General Features/Characteristics

<b>SOURCE</b>	Internationally sourced, normal, healthy, non-infected and non-immunized 100% Swiss mice.
<b>PURIFICATION</b>	Pooled serum recovered from whole blood collection. Off-clot and frozen immediately. Product is heat inactivated at 56 ±1°C for 30 ± 2 minutes, 0.2 µm filtered.
<b>GLOBULIN</b>	≥ 0.6 g/dL
<b>TOTAL PROTEIN</b>	> 5.5 g dL
<b>BIOBURDEN</b>	Direct Inoculation Acceptable Range: ≤ 100 CFU/mL

### Applications

Specific methodologies have not been tested using this product.