

Human CA242 Calibrator Grade

Cat. No. A01267H

STORE AT -30 to -10°C

DESCRIPTION	Human CA242 Calibrator Grade Human Cancer Antigen CA242 (Gastrointestinal Cancer) Calibrator Grade
TYPE	Antigen
SHIPPING & STORAGE INFORMATION	Store at -30 to -10°C. Avoid multiple freeze/thaw cycles.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SOURCE	Cell Culture
FORMAT	Calibrator, Liquid
BUFFER	Phosphate Buffered Saline, pH 7.4 ± 0.1
PRESERVATIVES	0.1% Sodium Azide
PURIFICATION	Cell Culture Concentrate Product is 0.2 µm filtered
INACTIVATION	Not Applicable

Applications

Specific methodologies have not been tested using this product