

Sm/RNP Antigen

Cat. No. A07304B

STORE AT -95 to -65°C

DESCRIPTION	Sm/RNP Complex Antigen
TYPE	Antigen
SHIPPING & STORAGE INFORMATION	Store at -95 to -65°C. Avoid multiple freeze/thaw cycles.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SPECIFICITY	Other antigens (SS-B, SS-A, Scl-70, Jo-1) are not detected with ELISA assays
SOURCE	Rabbit and bovine thymus of US origin
FORMAT	Affinity Purified, Liquid
BUFFER	0.01 M Tris, 0.53 M Sodium Chloride
PRESERVATIVES	0.02% Sodium Azide
PURIFICATION	Immobilized human anti-RNP immune globulins
PURITY	SDS PAGE
INACTIVATION	Not Applicable

Applications

Recommended for use in standard solid phase Immunoassays. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.