

TRU Block™3

Cat. No. A66803H

STORE AT -30 to -10°C

DESCRIPTION	TRU Block™3 is a complex blocker that contains active ingredients and purified mouse IgG.
TYPE	Universal (Active & Passive) Blocker
SHIPPING & STORAGE INFORMATION	Store at -30 to -10°C. Avoid multiple freeze thaw cycles.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SOURCE	Mouse Serum Normal mouse serum plus proprietary ingredients.
FORMAT	Purified, Liquid
BUFFER	10 mM Sodium Phosphate, 150 mM Sodium Chloride. pH 7.0 - 7.4
PRESERVATIVES	0.05% Sodium Azide
PURIFICATION	Protein A purification
PURITY	SDS-PAGE (reduced & non-reduced) Acceptable Range: ≥ 95%
CONCENTRATION	OD 280nm Assay ($E^{0.1\%} = 1.4$) Acceptable Range: 24 – 26 mg/mL
ACTIVITY	ELISA (% Blocking)
BIOBURDEN	Direct Inoculation Acceptable Range: ≤ 5 CFU/mL

Applications

Active blocking of heterophilic antibody interference in immunoassays.