

## Normal Mouse Serum, Filtered

Cat. No. N14010M

STORE AT -30 to -10°C

DESCRIPTION	Normal Mouse Serum, Filtered
TYPE	Blocker
SHIPPING & STORAGE INFORMATION	Store at -30 to -10°C. Avoid multiple freeze / thaw cycles.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

### General Features/Characteristics

SOURCE	Mouse Serum Normal mouse serum internationally sourced, normal, healthy, non-infected and non-immunized 100% Swiss mouse. Processed in the USA.
FORMAT	Neat, Liquid
PRESERVATIVES	None
PURIFICATION	Pooled serum recovered from whole blood collection. Off-clot and frozen immediately. Product is 0.6 µm filtered
PROTEIN CONCENTRATION	Acceptable Range: 6.3 – 6.6 g/dL
IgG CONCENTRATION	RID Acceptable Range: ≥ 6 mg/mL
pH	Acceptable Range 7.5 - 8.0

### Applications

Specific methodologies have not been tested using this product. Each laboratory should determine an optimum working titer for use in its particular application.