

HBsAg (adr) Recombinant

Cat. No. R01652

STORE AT 2 to 8°C

DESCRIPTION	Hepatitis B Surface Antigen (HBsAg) (adr) Recombinant
MOLECULAR WEIGHT	23 - 27 kDa
TYPE	Recombinant Antigen
SHIPPING & STORAGE INFORMATION	Aliquot and store at 2 to 8°C. DO NOT FREEZE!
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SOURCE	CHO Cells
FORMAT	Purified, Liquid
BUFFER	20 mM Phosphate Buffer and 154 mM Sodium Chloride, pH 7.2 ± 0.2
PRESERVATIVES	None
PURIFICATION	RP-HPLC
CONCENTRATION	Lowry Method
PURITY	> 95%

Applications

Specific methodologies have not been tested using this product.