

HIV-1 p24 Antigen, Recombinant

Cat. No. VTI340

STORE AT -95 to -65°C

DESCRIPTION	HIV-1 p24 Antigen, Recombinant
TYPE	Recombinant Antigen
SHIPPING & STORAGE INFORMATION	Store at -95 to -65°C. Avoid multiple freeze thaw cycles.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SOURCE	Pichia pastoris Product manufactured in the United States using only materials processed in U.S. facilities. No materials of animal origin.
FORMAT	Purified, Liquid
CONCENTRATION	BCA Assay Acceptable Range: > 0.5 mg/mL
IDENTITY	SDS-PAGE followed by Coomassie staining. must demonstrate a major band migrating at 23 kDa. Immunospecificity verified by Western Blot.
BUFFER	0.02 M Sodium Phosphate, pH 6.8 - 7.2
PRESERVATIVES	None

Applications

Specific methodologies have not been tested using this product. Each laboratory should determine an optimum working titer for use in its particular application.